

Alcohol Testing Program - Instrument Processing Sheet

Agency: PANAMA CITY PD

Instrument Serial Number: 80-001204

Date In: 7/13/2026

DI Completion Date: 7/15/2026

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: WKP Date: 7/13/2026		Quality Checks By: WKP Date: 7/13/2026		Flow Adjustment By:	
☒ Annual ☐ Dropped Off ☐ Registration ☒ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: 177 ☒ Flow Verification (L/s) Flow Column #: ATP 105 32 mm 0.152 (.139-.169) 36 mm 0.167 (.156-.190) 53 mm 0.250 (.228-.278) 103 mm 0.519 (.447-.547) ☒ Barometric Pressure Check Gauge ID #: 34419 Gauge: 1017 Instrument: 1017 ☒ Stability Checks		Flow Column #: ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: ☐ Post Adjustment Verification (L/S) Flow Column #: 32 mm (.139-.169) 36 mm (.156-.190) 53 mm (.228-.278) 103 mm (.447-.547)	
		Simulator		Serial #	Lot#/Exp
		0.050		MP5088	202603A 3/3/2028
		0.080		MP5089	202603B 3/3/2028
		0.200		MP5090	202603D 3/4/2028
		0.080 DGS		N/A	23825080A1 10/5/2027
Optical Bench Adjustment By:		Department Inspection By: WKP			
Barometric Pressure Gauge: ID#:		Barometric Pressure ID#: 34420			
Simulator	Serial #	Lot #	Expiration	Gauge: 1016	Instrument: 1017
0.000		N/A	N/A	Mouth Alcohol Solution Lot #: 2025-D	Exp: 9/25/2027
0.040				Acetone Stock Solution Lot #: 2025-B	Exp: 9/22/2027
0.100				Simulator	Serial Number
0.200				0.000	MP5086
0.300				Interferent	MP5087
0.080 DGS	N/A			0.050	MP5088
☐ Post Optical Bench Adjustment Stability Checks				0.080	MP5089
Simulator	Serial #	Lot #	Expiration	0.200	MP5090
0.050				Attachments	
0.080				☒ Form 41	
0.200				☒ Stability Checks	
0.080 DGS	N/A			☒ Calibration Certificate	
Gauge ID #:				☐ Optical Bench Adjustment	
Gauge: Instrument:				☐ Post-Stability Checks	
Notes/Suggested Service: Admin Correction: Selected Form 40 attachment box. WKP 7/21/2026				☐ Flow Adjustment	
				☐ Form 40 WKP 7/21/26	
				☒ Other: Diagnostic Test	
				☒ Instrument Complies with Chapter 11D-8, FAC	
				☐ Instrument Does Not Comply with Chapter 11D-8, FAC	
				☒ Return to/Place into Evidentiary Use	
				☐ Remain Out of Evidentiary Use	
				☒ Conduct an Agency Inspection Before Evidentiary Use	
				Taylor Gutschow Digitally signed by Taylor Gutschow Date: 2026.07.20 10:31:38 -04'00'	LeAndra Higginbotham Digitally signed by LeAndra Higginbotham Date: 2026.07.21 11:22:46 -04'00'
				Tech Review	Admin Review

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: SOUTH OGDEN
Time of Inspection: 13:46

Date of Inspection: 07/13/2026

Serial Number: 80 001204
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:_____ Exp:_____	0.08g/210L Test (g/210L) Lot#:_____ Exp:_____	0.20g/210L Test (g/210L) Lot#:_____ Exp:_____	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:_____ Exp:_____

Number of Simulators Used: _____

Remarks:

BYPASS AI TO OPERATE INSTRUMENT. COMPLIANCE NOT DETERMINED.

Printed when AI bypass conducted, prior to agency name change. WKP 7/13/2026

WKP 7/15/2026

The above instrument complies (X) does not comply () with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Wen-Chi K Pierson

WEN-CHI K PIERSON

Signature and Printed Name

07/13/2026
Date

Stability Checks

0.05g/210L 0.047 to 0.053	0.08g/210L 0.077 to 0.083	0.20g/210L 0.194 to 0.206	DGS 0.08g/210L 0.077 to 0.083 ☒ ≤ 0.003 of Wet ☒																																																																																																																																																
PANAMA CITY PD Intoxilyzer - Alcotest Analyzer Model 8000 17/13/2026 Software: 8100.27 SN 80-101204 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.001</td><td>14:25</td></tr> <tr><td>Control Test</td><td>0.049</td><td>14:25</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:26</td></tr> <tr><td>Control Test</td><td>0.048</td><td>14:26</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:27</td></tr> <tr><td>Control Test</td><td>0.048</td><td>14:28</td></tr> <tr><td>Air Blank</td><td>0.001</td><td>14:28</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0483</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.1945</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.001	14:25	Control Test	0.049	14:25	Air Blank	0.000	14:26	Control Test	0.048	14:26	Air Blank	0.000	14:27	Control Test	0.048	14:28	Air Blank	0.001	14:28	Control Test Stats			Average	0.0483		Std Dev	0.0006		Rel Std Dev(%)	1.1945		PANAMA CITY PD Intoxilyzer - Alcotest Analyzer Model 8000 17/13/2026 Software: 8100.27 SN 80-101204 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>14:33</td></tr> <tr><td>Control Test</td><td>0.078</td><td>14:33</td></tr> <tr><td>Air Blank</td><td>0.001</td><td>14:34</td></tr> <tr><td>Control Test</td><td>0.077</td><td>14:35</td></tr> <tr><td>Air Blank</td><td>0.001</td><td>14:35</td></tr> <tr><td>Control Test</td><td>0.078</td><td>14:36</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:36</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0777</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7434</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	14:33	Control Test	0.078	14:33	Air Blank	0.001	14:34	Control Test	0.077	14:35	Air Blank	0.001	14:35	Control Test	0.078	14:36	Air Blank	0.000	14:36	Control Test Stats			Average	0.0777		Std Dev	0.0006		Rel Std Dev(%)	0.7434		PANAMA CITY PD Intoxilyzer - Alcotest Analyzer Model 8000 17/13/2026 Software: 8100.27 SN 80-101204 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>14:39</td></tr> <tr><td>Control Test</td><td>0.200</td><td>14:39</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:40</td></tr> <tr><td>Control Test</td><td>0.200</td><td>14:41</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:41</td></tr> <tr><td>Control Test</td><td>0.200</td><td>14:42</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:42</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.2000</td><td></td></tr> <tr><td>Std Dev</td><td>0.0001</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0001</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	14:39	Control Test	0.200	14:39	Air Blank	0.000	14:40	Control Test	0.200	14:41	Air Blank	0.000	14:41	Control Test	0.200	14:42	Air Blank	0.000	14:42	Control Test Stats			Average	0.2000		Std Dev	0.0001		Rel Std Dev(%)	0.0001		PANAMA CITY PD Intoxilyzer - Alcotest Analyzer Model 8000 17/13/2026 Software: 8100.27 SN 80-101204 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>13:54</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:54</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:54</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:55</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:55</td></tr> <tr><td>Control Test</td><td>0.080</td><td>13:55</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:56</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0793</td><td></td></tr> <tr><td>Std Dev</td><td>0.0005</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7277</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	13:54	Control Test	0.079	13:54	Air Blank	0.000	13:54	Control Test	0.079	13:55	Air Blank	0.000	13:55	Control Test	0.080	13:55	Air Blank	0.000	13:56	Control Test Stats			Average	0.0793		Std Dev	0.0005		Rel Std Dev(%)	0.7277	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.001	14:25																																																																																																																																																	
Control Test	0.049	14:25																																																																																																																																																	
Air Blank	0.000	14:26																																																																																																																																																	
Control Test	0.048	14:26																																																																																																																																																	
Air Blank	0.000	14:27																																																																																																																																																	
Control Test	0.048	14:28																																																																																																																																																	
Air Blank	0.001	14:28																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0483																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	1.1945																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	14:33																																																																																																																																																	
Control Test	0.078	14:33																																																																																																																																																	
Air Blank	0.001	14:34																																																																																																																																																	
Control Test	0.077	14:35																																																																																																																																																	
Air Blank	0.001	14:35																																																																																																																																																	
Control Test	0.078	14:36																																																																																																																																																	
Air Blank	0.000	14:36																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0777																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7434																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	14:39																																																																																																																																																	
Control Test	0.200	14:39																																																																																																																																																	
Air Blank	0.000	14:40																																																																																																																																																	
Control Test	0.200	14:41																																																																																																																																																	
Air Blank	0.000	14:41																																																																																																																																																	
Control Test	0.200	14:42																																																																																																																																																	
Air Blank	0.000	14:42																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.2000																																																																																																																																																		
Std Dev	0.0001																																																																																																																																																		
Rel Std Dev(%)	0.0001																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	13:54																																																																																																																																																	
Control Test	0.079	13:54																																																																																																																																																	
Air Blank	0.000	13:54																																																																																																																																																	
Control Test	0.079	13:55																																																																																																																																																	
Air Blank	0.000	13:55																																																																																																																																																	
Control Test	0.080	13:55																																																																																																																																																	
Air Blank	0.000	13:56																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0793																																																																																																																																																		
Std Dev	0.0005																																																																																																																																																		
Rel Std Dev(%)	0.7277																																																																																																																																																		
Operator's Signature <i>menpi</i>	Operator's Signature <i>menpi</i>	Operator's Signature <i>menpi</i>	Operator's Signature <i>menpi</i>																																																																																																																																																

PANAMA CITY PD
Intoxilyzer - Alcohol Analyzer
Model: 8000 SN: 60-101274
07/15/2026
Software: 8100.27

DIAGNOSTICS

Voltage/Current Test	OK
RAM Test	OK
EEPROM Checksum Test	OK
Real Time Clock Test	OK
DSP Test	OK
Analytical Stability Test	OK
Internal Printer Test	OK
Modem Test	OK
Temperature Regulation Test	OK

Inadvertently
conducted a diagnostics
test.

WMP 07/15/2026

menpi

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: PANAMA CITY PD

Time of Inspection: 10:58

Date of Inspection: 07/15/2026

Serial Number: 80-001204

Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23725080A1 Exp: 10/05/2027
0.000	0.048	0.077	0.197	0.079
0.000	0.048	0.077	0.198	0.079
0.000	0.049	0.076	0.198	0.079
0.000	0.048	0.076	0.198	0.079
0.000	0.048	0.077	0.198	0.079
0.000	0.048	0.077	0.198	0.079
0.000	0.048	0.077	0.199	0.079
0.000	0.049	0.077	0.198	0.080
0.000	0.049	0.077	0.198	0.079
0.000	0.049	0.077	0.199	0.079

Standard Deviations	0.0005	0.0004	0.0005	0.0003
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0004 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Wen-Chi K Pierson Signature and Printed Name

WEN-CHI K PIERSON

07/15/2026
Date



Calibration Certificate

Calibration Certificate

Serial Number:	<u>80-001204</u>	UNCERTAINTY* ±
Owning Agency:	<u>PANAMA CITY PD</u>	0.050 g/ 210 L
Calibration Date:	<u>07/15/2026</u>	0.080 g/ 210 L
Calibration Time:	<u>10:58</u>	0.200 g/ 210 L
		0.080 g/ 210 L Dry Gas Control
		0.004
		0.004
		0.008
		0.005

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

07/15/2026

Wen-Chi Pierson

Digitally signed by Wen-Chi Pierson
Date: 2026.07.15 12:53:54 -04'00'

Wen-Chi
Pierson

Pierson

Date _____

WEN-CHI K PIERSON,
Department Inspector

Service • Integrity • Respect • Quality