

INSTRUMENT PROCESSING SHEET

Agency **PORT ORANGE PD**

S/N 80-001245

Florida Department of
Law Enforcement

Date In 9/24/25

DI Completion Date N/A

☐ Ship ☐ P/U ☐ H/D ☒ CMI ☐ EE

Intake By <u>KTS</u> Date <u>9/24/25</u>	Quality Checks By <u>WKP</u> Date <u>09/25/2025</u>	Flow Calibration By _____ Date _____																			
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____ _____ _____ _____ _____ _____ _____ _____ _____	☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>181</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP 103</u> 32 mm <u>0.152</u> (.139 - .169) 36 mm <u>0.167</u> (.156 - .190) 53 mm <u>0.242</u> (.228 - .278) 103 mm <u>0.511</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28421</u> ☒ Stability Checks	Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																			
	<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td rowspan="2">0.050</td><td rowspan="2">MP5088</td><td>202406K</td></tr><tr><td>06/19/2026</td></tr><tr><td rowspan="2">0.080</td><td rowspan="2"></td><td></td></tr><tr><td></td></tr><tr><td rowspan="2">0.200</td><td rowspan="2"></td><td></td></tr><tr><td></td></tr><tr><td rowspan="2">0.080 DGS</td><td rowspan="2">N/A</td><td></td></tr><tr><td></td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050	MP5088	202406K	06/19/2026	0.080				0.200				0.080 DGS	N/A			Maintenance By _____ Date _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____ _____ _____ _____
Simulator	Serial #	Lot #/Exp																			
0.050	MP5088	202406K																			
		06/19/2026																			
0.080																					
0.200																					
0.080 DGS	N/A																				

Calibration Adjustment				By _____	
Barometric Pressure Gauge _____			ID # _____		
Simulator	Serial #	Lot #	Expiration		
0.000		N/A	N/A		
0.040					
0.100					
0.200					
0.300					
0.080 DGS	N/A				
☐ Post Calibration Adjustment Stability Checks					
Simulator	Serial #	Lot #	Expiration		
0.050					
0.080					
0.200					
0.080 DGS	N/A				
Notes/Suggested Service: Ambient fail during 0.05g/210L stability check. Attempted to rerun stabilities on 10/9/25 but instrument DSP failed during instrument start up diagnostics. WKP 10/9/2025. In addition to ambient fail, DVM Monitoring of channels equalled zero (DSP fail), will send to repair. SLH 10/9/2025					

Department Inspection		By _____	
Barometric Pressure ID# _____			
Gauge _____		Instrument _____	
Mouth Alcohol Solution Lot # _____			
Acetone Stock Solution Lot # _____			
Simulator	Serial Number		
0.000			
Interferent			
0.050			
0.080			
0.200			
Attachments			
☐ Form 41 ☒ Stability Checks ☐ Calibration Certificate ☐ Calibration Adjustment		☐ Post-Stability Checks ☐ Flow Calibration ☐ Form 40 ☒ Other <u>Form 51</u>	
☐ Instrument Complies with Chapter 11D-8, FAC ☒ Instrument Does Not Comply with Chapter 11D-8, FAC			
☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use			
☐ Conduct an Agency Inspection Before Evidentiary Use			
Taylor Gutschow		Shayla Platt Digitally signed by Shayla Platt Date: 2025.10.10 14:54:01 -04'00'	
Tech Review / Date _____		Admin Review / Date _____	

INTOXILYZER 8100
Instrument Initialization
10:06 09/25/2025

PORT ORANGE PD
Intoxilyzer - Alcohol Analyzer
Model 8100 SN 80-001245
09/25/2025
Software: 8100.27

Test	g/210L	Time
Air Blank	AMB*	14:12
Air Blank	AMB*	14:12

*Ambient Fail

PORT ORANGE PD
Intoxilyzer - Alcohol Analyzer
Model 8100 SN 80-001245
09/25/2025
Software: 8100.27

Test	g/210L	Time
Air Blank	AMB*	14:05
Air Blank	AMB*	14:06

*Ambient Fail

menpi
Operator's Signature

menpi
Operator's Signature

Instrument+ AMB*

Fail during 0.05g/210L

Stability. Directed a

fan towards instrument

to clear out ambient air
and repeated stability.

WKP 9/25/25

Repeated 0.05g/210L

stability check after

RCA. AMB* fail occurred

again. Stabilities paused

while awaiting further

guidance.

WKP 9/25/25

Return Material Authorization

Ship to: ☒ CMI, Inc.

☐ Enforcement Electronics

Shipment to repair facility authorized by: William Harrison on 10/9/2025

Items Returned: Instrument ☒ Supplies ☐ Other ☐ Describe: _____

Instrument Model: Intoxilyzer 8000 Serial Number: 80-001245

Bill To Address:

PORT ORANGE PD

4545 S Clyde Morris Blvd

Port Orange, FL 32129

Ship to Address:

FDLE Tallahassee

Reason for Return:

Ambient Failing during stability checks even after use of fan. DSP failed diagnostic after
stabilities started with zero readings at DVM Monitoring.

Please choose one of the following options:

☐ 1. I _____, authorize all repairs.

☐ 2. I _____, authorize repairs up to \$_____.

☒ 3. I require an estimate **BEFORE** any repairs will be authorized and/ or conducted.

Please contact: Name: William Harrison

Phone #: 386-334-1109 Email: wharrison@port-orange.org

ATP Contact Name: LeAndra Higginbotham ATP Email: LeAndraHigginbotham@fdle.state.fl.us