

Alcohol Testing Program - Instrument Processing Sheet

Agency: Hialeah Police Department **Instrument Serial Number:** 80-000876

Date In: 7/30/2026 **DI Completion Date:** 8/4/2026 ☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: <u>TDG</u> Date: <u>8/3/2026</u>		Quality Checks By: <u>TDG</u> Date: <u>8/4/2026</u>		Flow Adjustment By: _____																																																																									
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☒ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>115</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP106</u> 32 mm <u>0.140</u> (.139-.169) 36 mm <u>0.156</u> (.156-.190) 53 mm <u>0.238</u> (.228-.278) 103 mm <u>0.500</u> (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>33363</u> Gauge: <u>1016</u> Instrument: <u>1015</u> ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)																																																																									
						<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 20%;">Serial #</th> <th style="width: 25%;">Lot#/Exp</th> <th style="width: 40%;">Maintenance By: _____ Date: _____</th> </tr> <tr> <td>0.050</td> <td>MP6286</td> <td>202603A 3/3/2028</td> <td rowspan="5"> ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other: _____ </td> </tr> <tr> <td>0.080</td> <td>MP6287</td> <td>202603B 3/3/2028</td> </tr> <tr> <td>0.200</td> <td>MP6288</td> <td>202603D 3/4/2028</td> </tr> <tr> <td rowspan="2">0.080 DGS</td> <td rowspan="2">N/A</td> <td>28424080A3</td> </tr> <tr> <td>11/5/2026</td> </tr> </table>		Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____	0.050	MP6286	202603A 3/3/2028	☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other: _____	0.080	MP6287	202603B 3/3/2028	0.200	MP6288	202603D 3/4/2028	0.080 DGS	N/A	28424080A3	11/5/2026																																																				
						Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____																																																																				
						0.050	MP6286	202603A 3/3/2028	☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other: _____																																																																				
						0.080	MP6287	202603B 3/3/2028																																																																					
0.200	MP6288	202603D 3/4/2028																																																																											
0.080 DGS	N/A	28424080A3																																																																											
		11/5/2026																																																																											
Optical Bench Adjustment By: _____ Department Inspection By: <u>TDG</u> Barometric Pressure Gauge: _____ ID#: _____ Barometric Pressure ID#: <u>33363</u> Gauge: <u>1017</u> Instrument: <u>1016</u> Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>9/25/2027</u> Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>9/22/2027</u>																																																																													
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 15%;">Serial #</th> <th style="width: 15%;">Lot #</th> <th style="width: 15%;">Expiration</th> <th style="width: 40%;">Gauge: <u>1017</u> Instrument: <u>1016</u></th> </tr> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> <td>Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>9/25/2027</u></td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> <td>Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>9/22/2027</u></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> <td>Simulator Serial Number</td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> <td>0.000 MP6284</td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> <td>Interferent MP6285</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> <td>0.050 MP6286</td> </tr> <tr> <td colspan="4">☐ Post Optical Bench Adjustment Stability Checks</td> <td>0.080 MP6287</td> </tr> <tr> <td>Simulator</td> <td>Serial #</td> <td>Lot #</td> <td>Expiration</td> <td>0.200 MP6288</td> </tr> <tr> <td>0.050</td> <td></td> <td></td> <td></td> <td rowspan="5" style="background-color: #d3d3d3;"> Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment ☐ Post-Stability Checks ☐ Flow Adjustment ☐ Form 40 ☒ Other: Calculation Sheet </td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> <tr> <td colspan="4">Gauge ID #: _____ Gauge: _____ Instrument: _____</td> </tr> <tr> <td colspan="4" rowspan="3" style="background-color: #d3d3d3;"> Notes/Suggested Service: </td> <td colspan="2" style="background-color: #d3d3d3;"> ☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use </td> </tr> <tr> <td colspan="2" style="background-color: #d3d3d3;"> Digitally signed by Wen-Chi Pierson Date: 2026.08.06 12:47:30 -04'00' Digitally signed by Kaitlyn Spearin Date: 2026.08.07 06:29:11 -04'00' </td> </tr> <tr> <td colspan="2" style="background-color: #d3d3d3;"> Tech Review Admin Review </td> </tr> </table>		Simulator	Serial #	Lot #	Expiration	Gauge: <u>1017</u> Instrument: <u>1016</u>	0.000		N/A	N/A	Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>9/25/2027</u>	0.040				Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>9/22/2027</u>	0.100				Simulator Serial Number	0.200				0.000 MP6284	0.300				Interferent MP6285	0.080 DGS	N/A			0.050 MP6286	☐ Post Optical Bench Adjustment Stability Checks				0.080 MP6287	Simulator	Serial #	Lot #	Expiration	0.200 MP6288	0.050				Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment ☐ Post-Stability Checks ☐ Flow Adjustment ☐ Form 40 ☒ Other: Calculation Sheet	0.080				0.200				0.080 DGS	N/A			Gauge ID #: _____ Gauge: _____ Instrument: _____				Notes/Suggested Service:				☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use		Digitally signed by Wen-Chi Pierson Date: 2026.08.06 12:47:30 -04'00' Digitally signed by Kaitlyn Spearin Date: 2026.08.07 06:29:11 -04'00'		Tech Review Admin Review	
Simulator	Serial #	Lot #	Expiration	Gauge: <u>1017</u> Instrument: <u>1016</u>																																																																									
0.000		N/A	N/A	Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>9/25/2027</u>																																																																									
0.040				Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>9/22/2027</u>																																																																									
0.100				Simulator Serial Number																																																																									
0.200				0.000 MP6284																																																																									
0.300				Interferent MP6285																																																																									
0.080 DGS	N/A			0.050 MP6286																																																																									
☐ Post Optical Bench Adjustment Stability Checks				0.080 MP6287																																																																									
Simulator	Serial #	Lot #	Expiration	0.200 MP6288																																																																									
0.050				Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment ☐ Post-Stability Checks ☐ Flow Adjustment ☐ Form 40 ☒ Other: Calculation Sheet																																																																									
0.080																																																																													
0.200																																																																													
0.080 DGS	N/A																																																																												
Gauge ID #: _____ Gauge: _____ Instrument: _____																																																																													
Notes/Suggested Service:				☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																																																																									
				Digitally signed by Wen-Chi Pierson Date: 2026.08.06 12:47:30 -04'00' Digitally signed by Kaitlyn Spearin Date: 2026.08.07 06:29:11 -04'00'																																																																									
				Tech Review Admin Review																																																																									

Stability Checks

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L																																																																																																																																																
0.047 to 0.053 ☒	0.077 to 0.083 ☒	0.194 to 0.206 ☒	0.077 to 0.083 ☒ ≤ 0.003 of Wet ☒																																																																																																																																																
HIALEAH PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-000876 08/04/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>10:49</td></tr> <tr><td>Control Test</td><td>0.052</td><td>10:50</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>10:51</td></tr> <tr><td>Control Test</td><td>0.051</td><td>10:51</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>10:52</td></tr> <tr><td>Control Test</td><td>0.049</td><td>10:53</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>10:53</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0507</td><td></td></tr> <tr><td>Std Dev</td><td>0.0015</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>3.0149</td><td></td></tr> </tbody> </table> <i>MC</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	10:49	Control Test	0.052	10:50	Air Blank	0.000	10:51	Control Test	0.051	10:51	Air Blank	0.000	10:52	Control Test	0.049	10:53	Air Blank	0.000	10:53	Control Test Stats			Average	0.0507		Std Dev	0.0015		Rel Std Dev(%)	3.0149		HIALEAH PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-000876 08/04/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>10:56</td></tr> <tr><td>Control Test</td><td>0.080</td><td>10:57</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>10:57</td></tr> <tr><td>Control Test</td><td>0.080</td><td>10:58</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>10:59</td></tr> <tr><td>Control Test</td><td>0.079</td><td>10:59</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:00</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0797</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7247</td><td></td></tr> </tbody> </table> <i>MC</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	10:56	Control Test	0.080	10:57	Air Blank	0.000	10:57	Control Test	0.080	10:58	Air Blank	0.000	10:59	Control Test	0.079	10:59	Air Blank	0.000	11:00	Control Test Stats			Average	0.0797		Std Dev	0.0006		Rel Std Dev(%)	0.7247		HIALEAH PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-000876 08/04/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>11:03</td></tr> <tr><td>Control Test</td><td>0.203</td><td>11:03</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:04</td></tr> <tr><td>Control Test</td><td>0.199</td><td>11:05</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:05</td></tr> <tr><td>Control Test</td><td>0.200</td><td>11:06</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:07</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.2007</td><td></td></tr> <tr><td>Std Dev</td><td>0.0021</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.0374</td><td></td></tr> </tbody> </table> <i>MC</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	11:03	Control Test	0.203	11:03	Air Blank	0.000	11:04	Control Test	0.199	11:05	Air Blank	0.000	11:05	Control Test	0.200	11:06	Air Blank	0.000	11:07	Control Test Stats			Average	0.2007		Std Dev	0.0021		Rel Std Dev(%)	1.0374		HIALEAH PD Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-000876 08/04/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>11:09</td></tr> <tr><td>Control Test</td><td>0.080</td><td>11:09</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:10</td></tr> <tr><td>Control Test</td><td>0.079</td><td>11:10</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:11</td></tr> <tr><td>Control Test</td><td>0.080</td><td>11:11</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:12</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0797</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7247</td><td></td></tr> </tbody> </table> <i>MC</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	11:09	Control Test	0.080	11:09	Air Blank	0.000	11:10	Control Test	0.079	11:10	Air Blank	0.000	11:11	Control Test	0.080	11:11	Air Blank	0.000	11:12	Control Test Stats			Average	0.0797		Std Dev	0.0006		Rel Std Dev(%)	0.7247	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	10:49																																																																																																																																																	
Control Test	0.052	10:50																																																																																																																																																	
Air Blank	0.000	10:51																																																																																																																																																	
Control Test	0.051	10:51																																																																																																																																																	
Air Blank	0.000	10:52																																																																																																																																																	
Control Test	0.049	10:53																																																																																																																																																	
Air Blank	0.000	10:53																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0507																																																																																																																																																		
Std Dev	0.0015																																																																																																																																																		
Rel Std Dev(%)	3.0149																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	10:56																																																																																																																																																	
Control Test	0.080	10:57																																																																																																																																																	
Air Blank	0.000	10:57																																																																																																																																																	
Control Test	0.080	10:58																																																																																																																																																	
Air Blank	0.000	10:59																																																																																																																																																	
Control Test	0.079	10:59																																																																																																																																																	
Air Blank	0.000	11:00																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0797																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7247																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	11:03																																																																																																																																																	
Control Test	0.203	11:03																																																																																																																																																	
Air Blank	0.000	11:04																																																																																																																																																	
Control Test	0.199	11:05																																																																																																																																																	
Air Blank	0.000	11:05																																																																																																																																																	
Control Test	0.200	11:06																																																																																																																																																	
Air Blank	0.000	11:07																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.2007																																																																																																																																																		
Std Dev	0.0021																																																																																																																																																		
Rel Std Dev(%)	1.0374																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	11:09																																																																																																																																																	
Control Test	0.080	11:09																																																																																																																																																	
Air Blank	0.000	11:10																																																																																																																																																	
Control Test	0.079	11:10																																																																																																																																																	
Air Blank	0.000	11:11																																																																																																																																																	
Control Test	0.080	11:11																																																																																																																																																	
Air Blank	0.000	11:12																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0797																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7247																																																																																																																																																		

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: HIALEAH PD
Time of Inspection: 14:56

Date of Inspection: 08/04/2026

Serial Number: 80-000876
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:28424080A3 Exp: 11/05/2026
0.000	0.051	0.080	0.201	0.079
0.000	0.051	0.079	0.201	0.079
0.000	0.051	0.079	0.201	0.079
0.000	0.051	0.080	0.201	0.079
0.000	0.050	0.080	0.201	0.079
0.000	0.051	0.080	0.201	0.079
0.000	0.051	0.080	0.201	0.079
0.000	0.051	0.080	0.200	0.079
0.000	0.051	0.080	0.201	0.079
0.000	0.051	0.080	0.201	0.078

Standard Deviations	0.0003	0.0004	0.0003	0.0003
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Taylor D Gutschow TAYLOR D GUTSCHOW
Signature and Printed Name

08/04/2026
Date

Florida Department of Law Enforcement
Alcohol Testing Program - Average Calculations for Calibration Certificate

Calculations By: TDG

Serial Number: 80-000876

Date: 8/5/2026

Solution:

0.05 g/210 L Test: 0.509 / 10 = 0.051 g/210L
(SUM) (AVERAGE)

0.08 g/210L Test: 0.798 / 10 = 0.080 g/210L
(SUM) (AVERAGE)

0.20 g/210L Test: 2.009 / 10 = 0.201 g/210L
(SUM) (AVERAGE)

0.08 g/210L Dry Gas Std Test: 0.789 / 10 = 0.079 g/210L
(SUM) (AVERAGE)

Comments:

Average = (Sum of observations) / (Total number of observations)

All results are reported in g/210L



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000876, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000876</u>	Average of $y \pm$ Uncertainty*	
Owning Agency:	<u>HIALEAH PD</u>	0.051 g/ 210 L	0.004
Calibration Date:	<u>08/04/2026</u>	0.080 g/210 L	0.004
Calibration Time:	<u>14:56</u>	0.201 g/ 210 L	0.008
		0.079 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

FDLE/ATP Form 69 August 2026

Issuing Authority: Alcohol Testing Program

08/04/2026

Date

Taylor
Gutschow

TAYLOR D GUTSCHOW,
Department Inspector

Digitally signed by Taylor
Gutschow
Date: 2026.08.05 13:34:57
-04'00'

Service • Integrity • Respect • Quality