

Alcohol Testing Program - Instrument Processing Sheet

Agency: Florida Highway Patrol

Instrument Serial Number: 80-003407

Date In: 5/7/2026

DI Completion Date: 5/8/2026

☐ Ship ☒ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: <u>IAS</u> Date: <u>5/7/2026</u>		Quality Checks By: <u>WKP</u> Date: <u>5/8/2026</u>		Flow Adjustment By: <u>WKP</u>																																															
☒ Annual ☒ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes: Arrived with back two missing feet. Added two new feet. IAS 5/7/2026		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>248</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP 105</u> 32 mm <u>0.125</u> (.139-.169) 36 mm <u>0.144</u> (.156-.190) 53 mm <u>0.218</u> (.228-.278) 103 mm <u>0.500</u> (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>34419</u> Gauge: <u>1012</u> Instrument: <u>1008</u> ☒ Stability Checks		Flow Column #: <u>ATP 103</u> ☒ 5L/min – 17mm ☒ 15L/min – 53mm ☒ 30L/min – 103mm ☒ R-Value: <u>246</u> ☒ Post Adjustment Verification (L/S) Flow Column #: <u>ATP 105</u> 32 mm <u>0.156</u> (.139-.169) 36 mm <u>0.171</u> (.156-.190) 53 mm <u>0.238</u> (.228-.278) 103 mm <u>0.507</u> (.447-.547)																																															
						<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 20%;">Serial #</th> <th style="width: 25%;">Lot#/Exp</th> <th style="width: 40%;">Maintenance By: _____ Date: _____</th> </tr> <tr> <td rowspan="2">0.050</td> <td rowspan="2">MP5088</td> <td>202603A</td> <td>☐ Battery Replacement</td> </tr> <tr> <td>3/3/2028</td> <td>☐ Dry Gas Regulator Replacement</td> </tr> <tr> <td rowspan="2">0.080</td> <td rowspan="2">MP5089</td> <td>202603B</td> <td>☐ Tank Sensor Tare</td> </tr> <tr> <td>3/3/2028</td> <td>☐ Breath Tube Replacement</td> </tr> <tr> <td rowspan="2">0.200</td> <td rowspan="2">MP5090</td> <td>202603D</td> <td>☐ Other:</td> </tr> <tr> <td>3/4/2028</td> <td></td> </tr> <tr> <td rowspan="2">0.080 DGS</td> <td rowspan="2">N/A</td> <td>23725080A1</td> <td></td> </tr> <tr> <td>10/5/2027</td> <td></td> </tr> </table>		Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____	0.050	MP5088	202603A	☐ Battery Replacement	3/3/2028	☐ Dry Gas Regulator Replacement	0.080	MP5089	202603B	☐ Tank Sensor Tare	3/3/2028	☐ Breath Tube Replacement	0.200	MP5090	202603D	☐ Other:	3/4/2028		0.080 DGS	N/A	23725080A1		10/5/2027																	
						Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____																																										
						0.050	MP5088	202603A	☐ Battery Replacement																																										
								3/3/2028	☐ Dry Gas Regulator Replacement																																										
0.080	MP5089	202603B	☐ Tank Sensor Tare																																																
		3/3/2028	☐ Breath Tube Replacement																																																
0.200	MP5090	202603D	☐ Other:																																																
		3/4/2028																																																	
0.080 DGS	N/A	23725080A1																																																	
		10/5/2027																																																	
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 20%;">Serial #</th> <th style="width: 25%;">Lot #</th> <th style="width: 40%;">Expiration</th> </tr> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 50%;">Gauge: <u>1013</u> Instrument: <u>1009</u></th> <th style="width: 50%;">Barometric Pressure ID#: <u>34419</u></th> </tr> <tr> <td>Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>9/25/2027</u></td> <td></td> </tr> <tr> <td>Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>9/22/2027</u></td> <td></td> </tr> <tr> <td colspan="2" style="background-color: #f2f2f2;"> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 50%;">Simulator</th> <th style="width: 50%;">Serial Number</th> </tr> <tr> <td>0.000</td> <td>MP5086</td> </tr> <tr> <td>Interferent</td> <td>MP5087</td> </tr> <tr> <td>0.050</td> <td>MP5088</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> </tr> </table> </td> </tr> </table>		Gauge: <u>1013</u> Instrument: <u>1009</u>	Barometric Pressure ID#: <u>34419</u>	Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>9/25/2027</u>		Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>9/22/2027</u>		<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 50%;">Simulator</th> <th style="width: 50%;">Serial Number</th> </tr> <tr> <td>0.000</td> <td>MP5086</td> </tr> <tr> <td>Interferent</td> <td>MP5087</td> </tr> <tr> <td>0.050</td> <td>MP5088</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> </tr> </table>		Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090
Simulator	Serial #	Lot #	Expiration																																																
0.000		N/A	N/A																																																
0.040																																																			
0.100																																																			
0.200																																																			
0.300																																																			
0.080 DGS	N/A																																																		
Gauge: <u>1013</u> Instrument: <u>1009</u>	Barometric Pressure ID#: <u>34419</u>																																																		
Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>9/25/2027</u>																																																			
Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>9/22/2027</u>																																																			
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 50%;">Simulator</th> <th style="width: 50%;">Serial Number</th> </tr> <tr> <td>0.000</td> <td>MP5086</td> </tr> <tr> <td>Interferent</td> <td>MP5087</td> </tr> <tr> <td>0.050</td> <td>MP5088</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> </tr> </table>		Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090																																						
Simulator	Serial Number																																																		
0.000	MP5086																																																		
Interferent	MP5087																																																		
0.050	MP5088																																																		
0.080	MP5089																																																		
0.200	MP5090																																																		
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 20%;">Serial #</th> <th style="width: 25%;">Lot #</th> <th style="width: 40%;">Expiration</th> </tr> <tr> <td>0.050</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </table>		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th colspan="2" style="background-color: #f2f2f2;">Attachments</th> </tr> <tr> <td style="width: 50%; vertical-align: top;"> ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment </td> <td style="width: 50%; vertical-align: top;"> ☐ Post-Stability Checks ☒ Flow Adjustment ☐ Form 40 ☐ Other: </td> </tr> </table>		Attachments		☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment	☐ Post-Stability Checks ☒ Flow Adjustment ☐ Form 40 ☐ Other:																								
Simulator	Serial #	Lot #	Expiration																																																
0.050																																																			
0.080																																																			
0.200																																																			
0.080 DGS	N/A																																																		
Attachments																																																			
☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment	☐ Post-Stability Checks ☒ Flow Adjustment ☐ Form 40 ☐ Other:																																																		
Gauge ID #: _____ Gauge: _____ Instrument: _____		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																																																	
Notes/Suggested Service:		<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; vertical-align: top;"> LeAndra Higginbotham Tech Review </td> <td style="width: 50%; vertical-align: top;"> Digitally signed by LeAndra Higginbotham Date: 2026.05.15 10:00:22 -04'00' Kaitlyn Spearin Digitally signed by Kaitlyn Spearin Date: 2026.05.18 12:09:45 -04'00' Admin Review </td> </tr> </table>		LeAndra Higginbotham Tech Review	Digitally signed by LeAndra Higginbotham Date: 2026.05.15 10:00:22 -04'00' Kaitlyn Spearin Digitally signed by Kaitlyn Spearin Date: 2026.05.18 12:09:45 -04'00' Admin Review																																														
LeAndra Higginbotham Tech Review	Digitally signed by LeAndra Higginbotham Date: 2026.05.15 10:00:22 -04'00' Kaitlyn Spearin Digitally signed by Kaitlyn Spearin Date: 2026.05.18 12:09:45 -04'00' Admin Review																																																		

FL HIGHWAY PATROL
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-003407
05/08/2026
Software: 8100.27

Flow Rate Calibration*****

1: Rate (Liters/min) = 5

$\text{SQRT}(\text{Diff}) = 6.707$

2: Rate (Liters/min) = 15

$\text{SQRT}(\text{Diff}) = 11.355$

3: Rate (Liters/min) = 30

$\text{SQRT}(\text{Diff}) = 21.117$

Dependent Data Scale Factor = 100000 L/min

Independent Data Scale Factor = 256

Rounded Slope = 666

Rounded Intercept = -559253

Correlation = 0.99621

Flow adjustment

WKP 5/8/2026

0.05g/210L 0.047 to 0.053		0.08g/210L 0.077 to 0.083		0.20g/210L 0.194 to 0.206		DGS 0.08g/210L 0.077 to 0.083 ≤0.003 of Wet																																																																																																																																														
FL HIGHWAY PATROL Intoxilyzer - Alcohol Analyzer Model: 8000 05/08/2026 Software: 8100.27 SN 81-003407 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>09:09</td></tr><tr><td>Control Test</td><td>0.049</td><td>09:10</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:11</td></tr><tr><td>Control Test</td><td>0.050</td><td>09:11</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:12</td></tr><tr><td>Control Test</td><td>0.049</td><td>09:12</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:13</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0493</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>1.1703</td><td></td></tr></tbody></table>	Test	g/210L	Time	Air Blank	0.000	09:09	Control Test	0.049	09:10	Air Blank	0.000	09:11	Control Test	0.050	09:11	Air Blank	0.000	09:12	Control Test	0.049	09:12	Air Blank	0.000	09:13	Control Test Stats			Average	0.0493		Std Dev	0.0006		Rel Std Dev(%)	1.1703		FL HIGHWAY PATROL Intoxilyzer - Alcohol Analyzer Model: 8000 05/08/2026 Software: 8100.27 SN 81-003407 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>09:14</td></tr><tr><td>Control Test</td><td>0.078</td><td>09:15</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:15</td></tr><tr><td>Control Test</td><td>0.079</td><td>09:16</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:16</td></tr><tr><td>Control Test</td><td>0.078</td><td>09:17</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:18</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0783</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7370</td><td></td></tr></tbody></table>	Test	g/210L	Time	Air Blank	0.000	09:14	Control Test	0.078	09:15	Air Blank	0.000	09:15	Control Test	0.079	09:16	Air Blank	0.000	09:16	Control Test	0.078	09:17	Air Blank	0.000	09:18	Control Test Stats			Average	0.0783		Std Dev	0.0006		Rel Std Dev(%)	0.7370		FL HIGHWAY PATROL Intoxilyzer - Alcohol Analyzer Model: 8000 05/08/2026 Software: 8100.27 SN 80-003407 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>09:19</td></tr><tr><td>Control Test</td><td>0.197</td><td>09:19</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:20</td></tr><tr><td>Control Test</td><td>0.196</td><td>09:20</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:21</td></tr><tr><td>Control Test</td><td>0.196</td><td>09:22</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:22</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.1963</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.2941</td><td></td></tr></tbody></table>	Test	g/210L	Time	Air Blank	0.000	09:19	Control Test	0.197	09:19	Air Blank	0.000	09:20	Control Test	0.196	09:20	Air Blank	0.000	09:21	Control Test	0.196	09:22	Air Blank	0.000	09:22	Control Test Stats			Average	0.1963		Std Dev	0.0006		Rel Std Dev(%)	0.2941		FL HIGHWAY PATROL Intoxilyzer - Alcohol Analyzer Model: 8000 05/08/2026 Software: 8100.27 SN 80-003407 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>09:24</td></tr><tr><td>Control Test</td><td>0.080</td><td>09:24</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:24</td></tr><tr><td>Control Test</td><td>0.081</td><td>09:25</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:25</td></tr><tr><td>Control Test</td><td>0.081</td><td>09:26</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:26</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0807</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7157</td><td></td></tr></tbody></table>	Test	g/210L	Time	Air Blank	0.000	09:24	Control Test	0.080	09:24	Air Blank	0.000	09:24	Control Test	0.081	09:25	Air Blank	0.000	09:25	Control Test	0.081	09:26	Air Blank	0.000	09:26	Control Test Stats			Average	0.0807		Std Dev	0.0006		Rel Std Dev(%)	0.7157		Operator's Signature Operator's Signature Operator's Signature Operator's Signature
Test	g/210L	Time																																																																																																																																																		
Air Blank	0.000	09:09																																																																																																																																																		
Control Test	0.049	09:10																																																																																																																																																		
Air Blank	0.000	09:11																																																																																																																																																		
Control Test	0.050	09:11																																																																																																																																																		
Air Blank	0.000	09:12																																																																																																																																																		
Control Test	0.049	09:12																																																																																																																																																		
Air Blank	0.000	09:13																																																																																																																																																		
Control Test Stats																																																																																																																																																				
Average	0.0493																																																																																																																																																			
Std Dev	0.0006																																																																																																																																																			
Rel Std Dev(%)	1.1703																																																																																																																																																			
Test	g/210L	Time																																																																																																																																																		
Air Blank	0.000	09:14																																																																																																																																																		
Control Test	0.078	09:15																																																																																																																																																		
Air Blank	0.000	09:15																																																																																																																																																		
Control Test	0.079	09:16																																																																																																																																																		
Air Blank	0.000	09:16																																																																																																																																																		
Control Test	0.078	09:17																																																																																																																																																		
Air Blank	0.000	09:18																																																																																																																																																		
Control Test Stats																																																																																																																																																				
Average	0.0783																																																																																																																																																			
Std Dev	0.0006																																																																																																																																																			
Rel Std Dev(%)	0.7370																																																																																																																																																			
Test	g/210L	Time																																																																																																																																																		
Air Blank	0.000	09:19																																																																																																																																																		
Control Test	0.197	09:19																																																																																																																																																		
Air Blank	0.000	09:20																																																																																																																																																		
Control Test	0.196	09:20																																																																																																																																																		
Air Blank	0.000	09:21																																																																																																																																																		
Control Test	0.196	09:22																																																																																																																																																		
Air Blank	0.000	09:22																																																																																																																																																		
Control Test Stats																																																																																																																																																				
Average	0.1963																																																																																																																																																			
Std Dev	0.0006																																																																																																																																																			
Rel Std Dev(%)	0.2941																																																																																																																																																			
Test	g/210L	Time																																																																																																																																																		
Air Blank	0.000	09:24																																																																																																																																																		
Control Test	0.080	09:24																																																																																																																																																		
Air Blank	0.000	09:24																																																																																																																																																		
Control Test	0.081	09:25																																																																																																																																																		
Air Blank	0.000	09:25																																																																																																																																																		
Control Test	0.081	09:26																																																																																																																																																		
Air Blank	0.000	09:26																																																																																																																																																		
Control Test Stats																																																																																																																																																				
Average	0.0807																																																																																																																																																			
Std Dev	0.0006																																																																																																																																																			
Rel Std Dev(%)	0.7157																																																																																																																																																			

Florida Department of Law Enforcement

Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: FL HIGHWAY PATROL
Time of Inspection: 11:12

Date of Inspection: 05/08/2026

Serial Number: 80-003407
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23725080A1 Exp: 10/05/2027
0.000	0.049	0.078	0.196	0.080
0.000	0.049	0.077	0.195	0.080
0.000	0.050	0.078	0.198	0.080
0.000	0.050	0.078	0.197	0.081
0.000	0.049	0.078	0.197	0.080
0.000	0.050	0.078	0.197	0.080
0.000	0.050	0.078	0.197	0.080
0.000	0.050	0.079	0.197	0.080
0.000	0.049	0.078	0.198	0.080
0.000	0.050	0.078	0.197	0.080

Standard Deviations	0.0005	0.0004	0.0008	0.0003
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0005 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

 _____
Signature and Printed Name

WEN-CHI K PIERSON

05/08/2026
Date



Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

Calibration Certificate

This is to certify the calibration of Intoxilyzer 8000 serial number 80-003407, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-003407</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>FL HIGHWAY PATROL</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>05/08/2026</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>11:12</u>	0.200 g/ 210 L 0.008
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

**Wen-Chi
Pierson**
Digitally signed by
Wen-Chi Pierson
Date: 2026.05.08
11:45:17 -04'00'

05/08/2026

Date

**WEN-CHI K PIERSON,
Department Inspector**

FDLE/ATP Form 69 January 2026

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality