

Return Material Authorization

Ship to:

☐ CMI, Inc.

☒ Enforcement Electronics

Shipment to repair facility authorized by: Ailisa Starnes on 1/31/2022

Items Returned: Instrument ☒ Supplies ☐ Other ☐ Describe: _____

Instrument Model: Intoxilyzer 8000 Serial Number: 80-001076

Bill To Address:

Zephyrhills Police Department

Attn: Sgt. Ailisa Starnes

Ship to Address:

Florida Department of Law Enforcement

Fort Myers Regional Operations Center

Attn: Alcohol Testing Program

4700 Terminal Drive, Suite 1

Fort Myers, FL 33907

Reason for Return:

Exterior of the analog phone jack broke off, and the interior retracted. Cannot perform uploads in the field. Uploaded all records via direct-connect. Verified no records are on the instrument.

Please choose one of the following options:

☐ 1. I _____, authorize all repairs.

☐ 2. I _____, authorize repairs up to \$_____.

☒ 3. I require an estimate **BEFORE** any repairs will be authorized and/ or conducted.

Please contact: Name: Ailisa Starnes

Phone #: 813-679-4826 Email: AStarnes@police.zephyrhills.fl.us

ATP Contact Name: Taylor Gutschow ATP Email: TaylorGutschow@fdle.state.fl.us



INSTRUMENT PROCESSING SHEET

Agency Zephyrhills PDS/N 80-001076Florida Department of
Law EnforcementDate In 2/1/2022

DI Completion Date _____

☐ Ship☐ P/U☐ H/D☐ CMI☒ EE

Intake	By TDG	Quality Checks	By	Date	Flow Calibration	By	Date																			
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☒ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: <u>Outside of analog phone jack is broken off (AI packaged broken piece in a baggie).</u> _____ _____ _____ _____ _____ _____ _____		☐ Breath Tube Screen ☐ Replace External O-Rings ☐ Instrument Set Up Verified ☐ R-Value _____ ☐ Flow Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547) ☐ Barometric Pressure Check Gauge ID # _____ ☐ Stability Checks			Flow Column # _____ ☐ 5L/min - 17mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																					
		<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050			0.080			0.200			0.080 DGS	N/A				<table border="1"><thead><tr><th>Maintenance</th><th>By</th></tr></thead><tbody><tr><td>☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____ _____</td><td></td></tr></tbody></table>	Maintenance	By	☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____ _____			
Simulator	Serial #	Lot #/Exp																								
0.050																										
0.080																										
0.200																										
0.080 DGS	N/A																									
Maintenance	By																									
☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____ _____																										

Calibration Adjustment	By	Department Inspection	By																																																													
Barometric Pressure Gauge _____ ID # _____ <table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table> ☐ Post Calibration Adjustment Stability Checks <table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>	Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A				Barometric Pressure ID# _____ Gauge _____ Instrument _____ Mouth Alcohol Solution Lot # _____ Acetone Stock Solution Lot # _____ <table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td></td></tr><tr><td>Interferent</td><td></td></tr><tr><td>0.050</td><td></td></tr><tr><td>0.080</td><td></td></tr><tr><td>0.200</td><td></td></tr></tbody></table>	Simulator	Serial Number	0.000		Interferent		0.050		0.080		0.200			
Simulator	Serial #	Lot #	Expiration																																																													
0.000		N/A	N/A																																																													
0.040																																																																
0.100																																																																
0.200																																																																
0.300																																																																
0.080 DGS	N/A																																																															
Simulator	Serial #	Lot #	Expiration																																																													
0.050																																																																
0.080																																																																
0.200																																																																
0.080 DGS	N/A																																																															
Simulator	Serial Number																																																															
0.000																																																																
Interferent																																																																
0.050																																																																
0.080																																																																
0.200																																																																
		<table border="1"><thead><tr><th>Attachments</th></tr></thead><tbody><tr><td>☐ Form 41 ☐ Stability Checks ☐ Calibration Certificate ☐ Calibration Adjustment</td><td>☐ Post-Stability Checks ☐ Flow Calibration ☐ Form 40 ☒ Other <u>Form 51</u></td></tr></tbody></table>	Attachments	☐ Form 41 ☐ Stability Checks ☐ Calibration Certificate ☐ Calibration Adjustment	☐ Post-Stability Checks ☐ Flow Calibration ☐ Form 40 ☒ Other <u>Form 51</u>																																																											
Attachments																																																																
☐ Form 41 ☐ Stability Checks ☐ Calibration Certificate ☐ Calibration Adjustment	☐ Post-Stability Checks ☐ Flow Calibration ☐ Form 40 ☒ Other <u>Form 51</u>																																																															

Notes/Suggested Service: <u>AI is unable to upload with the exterior of the analog phone jack broken off and the rest of it retracted into the instrument. Direct-connected to the instrument and uploaded records. Will send to EE for repair (TDG)</u> _____ _____ _____ _____	☐ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use
	Tech Review / Date _____ Admin Review / Date _____

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: ZEPHYRHILLS PD
Time of Inspection: 11:32

Date of Inspection: 11/02/2022

Serial Number: 80-001076
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Number of Simulators Used: _____

Remarks:

BYPASSED TO OPERATE INSTRUMENT. AI NOT CONDUCTED.

Not determined

MG
11/02/2022

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Taylor D Gutschow

TAYLOR D GUTSCHOW
Signature and Printed Name

11/02/2022
Date

Type of Test	Serial Number	Agency	Date	Performed By
Stabilities	80-001076	Zephyrhills PD	11/02/2022	TDG <i>MG</i>

0.05g/210L		0.08g/210L		0.20g/210L		DGS 0.08g/210L																																																																																																																																																	
0.047 to 0.053		0.077 to 0.083		0.194 to 0.206		0.077 to 0.083																																																																																																																																																	
✓		✓		✓		✓																																																																																																																																																	
ZEPHYRHILLS PD Intoxilyzer - Alcohol Analyzer Model 8000 11/02/2022 Software: 8100.27 SN 80-001076 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>11:58</td></tr><tr><td>Control Test</td><td>0.050</td><td>11:59</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:00</td></tr><tr><td>Control Test</td><td>0.050</td><td>12:00</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:01</td></tr><tr><td>Control Test</td><td>0.050</td><td>12:02</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:02</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0500</td><td></td></tr><tr><td>Std Dev</td><td>0.0000</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr></tbody></table> <i>MG</i> Operator's Signature		Test	g/210L	Time	Air Blank	0.000	11:58	Control Test	0.050	11:59	Air Blank	0.000	12:00	Control Test	0.050	12:00	Air Blank	0.000	12:01	Control Test	0.050	12:02	Air Blank	0.000	12:02	Control Test Stats			Average	0.0500		Std Dev	0.0000		Rel Std Dev(%)	0.0000		ZEPHYRHILLS PD Intoxilyzer - Alcohol Analyzer Model 8000 11/02/2022 Software: 8100.27 SN 80-001076 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>12:06</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:07</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:07</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:08</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:08</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:09</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:10</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0790</td><td></td></tr><tr><td>Std Dev</td><td>0.0000</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr></tbody></table> <i>MG</i> Operator's Signature		Test	g/210L	Time	Air Blank	0.000	12:06	Control Test	0.079	12:07	Air Blank	0.000	12:07	Control Test	0.079	12:08	Air Blank	0.000	12:08	Control Test	0.079	12:09	Air Blank	0.000	12:10	Control Test Stats			Average	0.0790		Std Dev	0.0000		Rel Std Dev(%)	0.0000		ZEPHYRHILLS PD Intoxilyzer - Alcohol Analyzer Model 8000 11/02/2022 Software: 8100.27 SN 80-001076 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>12:14</td></tr><tr><td>Control Test</td><td>0.200</td><td>12:15</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:15</td></tr><tr><td>Control Test</td><td>0.201</td><td>12:16</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:16</td></tr><tr><td>Control Test</td><td>0.200</td><td>12:17</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:18</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.2003</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.2882</td><td></td></tr></tbody></table> <i>MG</i> Operator's Signature		Test	g/210L	Time	Air Blank	0.000	12:14	Control Test	0.200	12:15	Air Blank	0.000	12:15	Control Test	0.201	12:16	Air Blank	0.000	12:16	Control Test	0.200	12:17	Air Blank	0.000	12:18	Control Test Stats			Average	0.2003		Std Dev	0.0006		Rel Std Dev(%)	0.2882		<i>DGS</i> ZEPHYRHILLS PD Intoxilyzer - Alcohol Analyzer Model 8000 11/02/2022 Software: 8100.27 SN 80-001076 <table><thead><tr><th>Test</th><th>g/210L</th><th>Time</th></tr></thead><tbody><tr><td>Air Blank</td><td>0.000</td><td>12:22</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:22</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:22</td></tr><tr><td>Control Test</td><td>0.078</td><td>12:23</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:23</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:24</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:24</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0787</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7339</td><td></td></tr></tbody></table> <i>MG</i> Operator's Signature		Test	g/210L	Time	Air Blank	0.000	12:22	Control Test	0.079	12:22	Air Blank	0.000	12:22	Control Test	0.078	12:23	Air Blank	0.000	12:23	Control Test	0.079	12:24	Air Blank	0.000	12:24	Control Test Stats			Average	0.0787		Std Dev	0.0006		Rel Std Dev(%)	0.7339	
Test	g/210L	Time																																																																																																																																																					
Air Blank	0.000	11:58																																																																																																																																																					
Control Test	0.050	11:59																																																																																																																																																					
Air Blank	0.000	12:00																																																																																																																																																					
Control Test	0.050	12:00																																																																																																																																																					
Air Blank	0.000	12:01																																																																																																																																																					
Control Test	0.050	12:02																																																																																																																																																					
Air Blank	0.000	12:02																																																																																																																																																					
Control Test Stats																																																																																																																																																							
Average	0.0500																																																																																																																																																						
Std Dev	0.0000																																																																																																																																																						
Rel Std Dev(%)	0.0000																																																																																																																																																						
Test	g/210L	Time																																																																																																																																																					
Air Blank	0.000	12:06																																																																																																																																																					
Control Test	0.079	12:07																																																																																																																																																					
Air Blank	0.000	12:07																																																																																																																																																					
Control Test	0.079	12:08																																																																																																																																																					
Air Blank	0.000	12:08																																																																																																																																																					
Control Test	0.079	12:09																																																																																																																																																					
Air Blank	0.000	12:10																																																																																																																																																					
Control Test Stats																																																																																																																																																							
Average	0.0790																																																																																																																																																						
Std Dev	0.0000																																																																																																																																																						
Rel Std Dev(%)	0.0000																																																																																																																																																						
Test	g/210L	Time																																																																																																																																																					
Air Blank	0.000	12:14																																																																																																																																																					
Control Test	0.200	12:15																																																																																																																																																					
Air Blank	0.000	12:15																																																																																																																																																					
Control Test	0.201	12:16																																																																																																																																																					
Air Blank	0.000	12:16																																																																																																																																																					
Control Test	0.200	12:17																																																																																																																																																					
Air Blank	0.000	12:18																																																																																																																																																					
Control Test Stats																																																																																																																																																							
Average	0.2003																																																																																																																																																						
Std Dev	0.0006																																																																																																																																																						
Rel Std Dev(%)	0.2882																																																																																																																																																						
Test	g/210L	Time																																																																																																																																																					
Air Blank	0.000	12:22																																																																																																																																																					
Control Test	0.079	12:22																																																																																																																																																					
Air Blank	0.000	12:22																																																																																																																																																					
Control Test	0.078	12:23																																																																																																																																																					
Air Blank	0.000	12:23																																																																																																																																																					
Control Test	0.079	12:24																																																																																																																																																					
Air Blank	0.000	12:24																																																																																																																																																					
Control Test Stats																																																																																																																																																							
Average	0.0787																																																																																																																																																						
Std Dev	0.0006																																																																																																																																																						
Rel Std Dev(%)	0.7339																																																																																																																																																						

Comments:

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: ZEPHYRHILLS PD
Time of Inspection: 15:06

Date of Inspection: 11/02/2022

Serial Number: 80-001076
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202201C Exp: 01/11/2024	0.08g/210L Test (g/210L) Lot#:202201D Exp: 01/18/2024	0.20g/210L Test (g/210L) Lot#:202201E Exp: 01/18/2024	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:00521080A2 Exp: 02/05/2023
0.000	0.050	0.079	0.201	0.079
0.000	0.049	0.079	0.200	0.079
0.000	0.049	0.079	0.200	0.079
0.000	0.050	0.079	0.201	0.079
0.000	0.050	0.079	0.201	0.079
0.000	0.050	0.079	0.201	0.079
0.000	0.049	0.078	0.201	0.079
0.000	0.049	0.079	0.200	0.078
0.000	0.050	0.079	0.201	0.079
0.000	0.050	0.079	0.201	0.079

Standard Deviations	0.0005	0.0003	0.0004	0.0003
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.


I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



Signature and Printed Name

TAYLOR D GUTSCHOW

11/02/2022
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001076, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001076</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>ZEPHYRHILLS PD</u>	0.050 g/210 L 0.004
Calibration Date:	<u>11/02/2022</u>	0.080 g/210 L 0.004
Calibration Time:	<u>15:06</u>	0.200 g/210 L 0.007
		0.080 g/210 L Dry Gas Control 0.005

All results are reported in g/210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

11/02/2022

Date


TAYLOR D GUTSCHOW,

Department Inspector

FDLE/ATP Form 69 December 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality