

Alcohol Testing Program - Instrument Processing Sheet

Agency: Haines City PD

Instrument Serial Number: 80-000796

Date In: 9/16/2026

DI Completion Date: 9/16/2026

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: <u>IAS</u> Date: <u>9/16/2026</u>		Quality Checks By: <u>IAS</u> Date: <u>9/16/2026</u>		Flow Adjustment By: _____																																									
☒ Annual ☐ Dropped Off ☐ Registration ☒ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>194</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP103</u> 32 mm <u>.140</u> (.139-.169) 36 mm <u>.160</u> (.156-.190) 53 mm <u>.230</u> (.228-.278) 103 mm <u>.496</u> (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>28427</u> Gauge: <u>1018</u> Instrument: <u>1021</u> ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)																																									
		<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot#/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202603A 3/3/2028</td></tr><tr><td>0.080</td><td>MP5089</td><td>202603B 3/3/2028</td></tr><tr><td>0.200</td><td>MP5090</td><td>202603D 3/4/2028</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>23825080A1 10/5/2027</td></tr></tbody></table>		Simulator	Serial #	Lot#/Exp	0.050	MP5088	202603A 3/3/2028	0.080	MP5089	202603B 3/3/2028	0.200	MP5090	202603D 3/4/2028	0.080 DGS	N/A	23825080A1 10/5/2027	<table border="1"><thead><tr><th>Maintenance</th><th>By:</th><th>Date:</th></tr></thead><tbody><tr><td>☐ Battery Replacement</td><td></td><td></td></tr><tr><td>☐ Dry Gas Regulator Replacement</td><td></td><td></td></tr><tr><td>☐ Tank Sensor Tare</td><td></td><td></td></tr><tr><td>☐ Breath Tube Replacement</td><td></td><td></td></tr><tr><td>☐ Other:</td><td></td><td></td></tr></tbody></table>		Maintenance	By:	Date:	☐ Battery Replacement			☐ Dry Gas Regulator Replacement			☐ Tank Sensor Tare			☐ Breath Tube Replacement			☐ Other:									
Simulator	Serial #	Lot#/Exp																																											
0.050	MP5088	202603A 3/3/2028																																											
0.080	MP5089	202603B 3/3/2028																																											
0.200	MP5090	202603D 3/4/2028																																											
0.080 DGS	N/A	23825080A1 10/5/2027																																											
Maintenance	By:	Date:																																											
☐ Battery Replacement																																													
☐ Dry Gas Regulator Replacement																																													
☐ Tank Sensor Tare																																													
☐ Breath Tube Replacement																																													
☐ Other:																																													
Optical Bench Adjustment By: _____		Department Inspection By: <u>IAS</u>																																											
Barometric Pressure Gauge: _____ ID#: _____		Barometric Pressure ID#: <u>28662</u>																																											
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Gauge: <u>1019</u> Instrument: <u>1021</u> Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>9/25/2027</u> Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>9/22/2027</u>		<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5086</td></tr><tr><td>Interferent</td><td>MP5087</td></tr><tr><td>0.050</td><td>MP5088</td></tr><tr><td>0.080</td><td>MP5089</td></tr><tr><td>0.200</td><td>MP5090</td></tr></tbody></table>		Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090
Simulator	Serial #	Lot #	Expiration																																										
0.000		N/A	N/A																																										
0.040																																													
0.100																																													
0.200																																													
0.300																																													
0.080 DGS	N/A																																												
Simulator	Serial Number																																												
0.000	MP5086																																												
Interferent	MP5087																																												
0.050	MP5088																																												
0.080	MP5089																																												
0.200	MP5090																																												
☐ Post Optical Bench Adjustment Stability Checks		Attachments																																											
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment ☒ Average Calculations		☐ Post-Stability Checks ☐ Flow Adjustment ☐ Form 40 ☐ Other:																					
Simulator	Serial #	Lot #	Expiration																																										
0.050																																													
0.080																																													
0.200																																													
0.080 DGS	N/A																																												
Gauge ID #: _____ Gauge: _____ Instrument: _____																																													
Notes/Suggested Service:		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																																											
		<table border="0"><tr><td>S. LeAndra Higginbotham</td><td>Digitally signed by S. LeAndra Higginbotham Date: 2026.09.17 09:40:32 -04'00'</td><td>Wen-Chi Pierson</td><td>Digitally signed by Wen-Chi Pierson Date: 2026.09.21 13:21:54 -04'00'</td></tr></table>				S. LeAndra Higginbotham	Digitally signed by S. LeAndra Higginbotham Date: 2026.09.17 09:40:32 -04'00'	Wen-Chi Pierson	Digitally signed by Wen-Chi Pierson Date: 2026.09.21 13:21:54 -04'00'																																				
S. LeAndra Higginbotham	Digitally signed by S. LeAndra Higginbotham Date: 2026.09.17 09:40:32 -04'00'	Wen-Chi Pierson	Digitally signed by Wen-Chi Pierson Date: 2026.09.21 13:21:54 -04'00'																																										
		<table border="0"><tr><td>Tech Review</td><td>Admin Review</td></tr></table>				Tech Review	Admin Review																																						
Tech Review	Admin Review																																												

Stability Checks

0.050 g/210L	0.080 g/210L	0.200 g/210L	DGS 0.080 g/210L
0.047 to 0.053 g/210L	0.077 to 0.083 g/210L	0.194 to 0.206 g/210L	0.077 to 0.083 g/210L ≤0.003 g/210L of Wet
HAINES CITY PD Intoxilyzer - Alcohol Analyzer Model 8000 19/16/2026 Software: 8.00.27 SN 80-000755 Test g/210L Time Air Blank0.00014:30 Control Test0.05214:31 Air Blank0.00014:31 Control Test0.05014:32 Air Blank0.00014:33 Control Test0.05014:33 Air Blank0.00014:34 Control Test Stats Average0.0507 Std Dev0.0012 Rel Std Dev(%)2.2790	HAINES CITY PD Intoxilyzer - Alcohol Analyzer Model 8000 09/16/2026 Software: 8.00.27 SN 80-000755 Test g/210L Time Air Blank0.00014:35 Control Test0.07914:36 Air Blank0.00014:36 Control Test0.07914:37 Air Blank0.00014:37 Control Test0.07914:38 Air Blank0.00014:38 Control Test Stats Average0.0791 Std Dev0.0001 Rel Std Dev(%)0.0000	HAINES CITY PD Intoxilyzer - Alcohol Analyzer Model 8000 09/16/2026 Software: 8.00.27 SN 80-000755 Test g/210L Time Air Blank0.00014:39 Control Test0.19714:40 Air Blank0.00014:41 Control Test0.19714:41 Air Blank0.00014:41 Control Test0.19614:42 Air Blank0.00014:43 Control Test Stats Average0.1967 Std Dev0.0006 Rel Std Dev(%)0.2936	HAINES CITY PD Intoxilyzer - Alcohol Analyzer Model 8000 09/16/2026 Software: 8.00.27 SN 80-000755 Test g/210L Time Air Blank0.00014:43 Control Test0.08014:44 Air Blank0.00014:44 Control Test0.07914:44 Air Blank0.00014:45 Control Test0.08014:45 Air Blank0.00014:46 Control Test Stats Average0.0797 Std Dev0.0006 Rel Std Dev(%)0.7247

Isabella Szymanski
Operator's Signature

Isabella Szymanski
Operator's Signature

Isabella Szymanski
Operator's Signature

Isabella Szymanski
Operator's Signature

Florida Department of Law Enforcement

Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: HAINES CITY PD

Time of Inspection: 16:28

Date of Inspection: 09/16/2026

Serial Number: 80-000796

Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23825080A1 Exp: 10/05/2027
/ 0.000	0.050	0.078	0.196	0.080
/ 0.000	0.050	0.078	0.196	0.080
/ 0.000	0.050	0.078	0.195	0.080
/ 0.000	0.050	0.078	0.195	0.080
/ 0.000	0.050	0.078	0.196	0.079
/ 0.000	0.050	0.077	0.195	0.080
/ 0.000	0.050	0.077	0.195	0.080
/ 0.000	0.050	0.078	0.195	0.080
/ 0.000	0.050	0.078	0.195	0.079
/ 0.000	0.050	0.078	0.195	0.079

Standard Deviations	0.0000	0.0004	0.0004	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

00: Ambient Fail. *Cleared area around instrument, turned on fans and reran. IAS 9/16/26*

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Isabella Smykowski ISABELLA A SMYKOWSKI
Signature and Printed Name

09/16/2026
Date

Florida Department of Law Enforcement
Alcohol Testing Program - 8000 Average Calculations for Calibration Certificate

Dept. Inspector: **IAS**

Serial Number: **80-000796**

Date: **9/16/2026**

Values are from Department Inspection Report - Intoxilyzer 8000 (Form 41)

Solutions:	Alcohol Free Test	0.05g/210L Test	0.08g/210L Test	0.20g/210L Test	0.08 g/210L Dry Gas Std Test
	0.000	0.050	0.078	0.196	0.080
	0.000	0.050	0.078	0.196	0.080
	0.000	0.050	0.078	0.195	0.080
	0.000	0.050	0.078	0.195	0.080
	0.000	0.050	0.078	0.196	0.079
	0.000	0.050	0.077	0.195	0.080
	0.000	0.050	0.077	0.195	0.080
	0.000	0.050	0.078	0.195	0.080
	0.000	0.050	0.078	0.195	0.079
	0.000	0.050	0.078	0.195	0.079

SUM:

0.000	0.500	0.778	1.953	0.797
-------	-------	-------	-------	-------

AVERAGE:

0.000	0.050	0.078	0.195	0.080
--------------	--------------	--------------	--------------	--------------

Average = (Sum of observations) / (Total number of observations)

All results are reported in g/210L

Comments:

--



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000796, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000796</u>	Average of $y \pm$ Uncertainty*
Owning Agency:	<u>HAINES CITY PD</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>09/16/2026</u>	0.078 g/ 210 L 0.004
Calibration Time:	<u>16:28</u>	0.195 g/ 210 L 0.008
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Isabella
Smykowski
Digitally signed by Isabella
Smykowski
Date: 2026.09.16 17:00:50 -04'00'

09/16/2026

Date

ISABELLA A SMYKOWSKI,
Department Inspector

FDLE/ATP Form 69 August 2026

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Alcohol Testing Program - Instrument Processing Sheet

Agency: HAINES CITY PD

Instrument Serial Number: 80-000796

Date In: 4/30/2026

DI Completion Date: N/A

☒ Ship

☐ P/U

☐ H/D

☐ CMI

☒ EE

Intake By: <u>WKP</u> Date: <u>4/30/2026</u>		Quality Checks By: <u>IAS</u> Date: <u>5/1/2026</u>		Flow Adjustment By: _____																
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>203</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP102</u> 32 mm <u>0.140</u> (.139-.169) 36 mm <u>0.156</u> (.156-.190) 53 mm <u>0.230</u> (.228-.278) 103 mm <u>0.496</u> (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>34419</u> Gauge: <u>1010</u> Instrument: <u>1008</u> ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)																
		<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 20%;">Serial #</th> <th style="width: 25%;">Lot#/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP5088</td> <td>202603A 3/3/2028</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> <td>202603B 3/3/2028</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> <td>202603D 3/4/2028</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>23825080A1 10/5/2027</td> </tr> </tbody> </table>		Simulator	Serial #	Lot#/Exp	0.050	MP5088	202603A 3/3/2028	0.080	MP5089	202603B 3/3/2028	0.200	MP5090	202603D 3/4/2028	0.080 DGS	N/A	23825080A1 10/5/2027	Maintenance By: _____ Date: _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other:	
Simulator	Serial #	Lot#/Exp																		
0.050	MP5088	202603A 3/3/2028																		
0.080	MP5089	202603B 3/3/2028																		
0.200	MP5090	202603D 3/4/2028																		
0.080 DGS	N/A	23825080A1 10/5/2027																		
Optical Bench Adjustment By: _____		Department Inspection By: <u>KTS</u>																		
Barometric Pressure Gauge: _____ ID#: _____		Barometric Pressure ID#: _____																		
Simulator	Serial #	Lot #	Expiration	Gauge: _____ Instrument: _____																
0.000		N/A	N/A	Mouth Alcohol Solution Lot #: _____ Exp: _____																
0.040				Acetone Stock Solution Lot #: _____ Exp: _____																
0.100				Simulator	Serial Number															
0.200				0.000																
0.300				Interferent																
0.080 DGS	N/A			0.050																
☐ Post Optical Bench Adjustment Stability Checks				0.080																
Simulator	Serial #	Lot #	Expiration	0.200																
0.050				Attachments																
0.080				☒ Form 41 ☒ Stability Checks ☐ Calibration Certificate ☐ Optical Bench Adjustment ☐ Post-Stability Checks ☐ Flow Adjustment ☐ Form 40 ☒ Other: Return Material Authorization.																
0.200																				
0.080 DGS	N/A																			
Gauge ID #: _____ Gauge: _____ Instrument: _____																				
Notes/Suggested Service: Restarted instrument after pre-inspection diagnostic fail. Instrument will not go into ready mode- still DSP failing. KTS 5/4/26 Tech Review: added stabilities and check marked stabilities on attachments. 5/18/26 Added Tallahassee to repair form. KTS 5/19/26				☐ Instrument Complies with Chapter 11D-8, FAC ☒ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use																
				LeAndra Higginbotham Digitally signed by LeAndra Higginbotham Date: 2026.05.20 10:23:14 -04'00' Wen-Chi Pierson Digitally signed by Wen-Chi Pierson Date: 2026.05.20 11:44:38 -04'00'																
Tech Review				Admin Review																

Stability Checks

0.050 g/210L 0.047 to 0.053 g/210L	0.080 g/210L 0.077 to 0.083 g/210L	0.200 g/210L 0.194 to 0.206 g/210L	DGS 0.080 g/210L 0.077 to 0.083 g/210L ±0.003 g/210L of Wet																																																																																																																																													
HAINES CITY PD Intoxilyzer - Alcohol Analyzer Model: 8000 05/01/2026 Software: 8100.27 SN 80-000756 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:33</td></tr> <tr><td>Control Test</td><td>0.050</td><td>12:39</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:40</td></tr> <tr><td>Control Test</td><td>0.051</td><td>12:41</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:42</td></tr> <tr><td>Control Test</td><td>0.050</td><td>12:42</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0503</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.1471</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	12:33	Control Test	0.050	12:39	Air Blank	0.000	12:40	Control Test	0.051	12:41	Air Blank	0.000	12:42	Control Test	0.050	12:42	Control Test Stats			Average	0.0503		Std Dev	0.0006		Rel Std Dev(%)	1.1471		HAINES CITY PD Intoxilyzer - Alcohol Analyzer Model: 8000 05/01/2026 Software: 8100.27 SN 80-000756 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:59</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:00</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:00</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:01</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:02</td></tr> <tr><td>Control Test</td><td>0.081</td><td>13:02</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:03</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0797</td><td></td></tr> <tr><td>Std Dev</td><td>0.0012</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.4494</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	12:59	Control Test	0.079	13:00	Air Blank	0.000	13:00	Control Test	0.079	13:01	Air Blank	0.000	13:02	Control Test	0.081	13:02	Air Blank	0.000	13:03	Control Test Stats			Average	0.0797		Std Dev	0.0012		Rel Std Dev(%)	1.4494		HAINES CITY PD Intoxilyzer - Alcohol Analyzer Model: 8000 05/01/2026 Software: 8100.27 SN 80-000756 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:47</td></tr> <tr><td>Control Test</td><td>0.198</td><td>12:49</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:49</td></tr> <tr><td>Control Test</td><td>0.199</td><td>12:49</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:49</td></tr> <tr><td>Control Test</td><td>0.199</td><td>12:50</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:50</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.1987</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2906</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	12:47	Control Test	0.198	12:49	Air Blank	0.000	12:49	Control Test	0.199	12:49	Air Blank	0.000	12:49	Control Test	0.199	12:50	Air Blank	0.000	12:50	Control Test Stats			Average	0.1987		Std Dev	0.0006		Rel Std Dev(%)	0.2906		HAINES CITY PD Intoxilyzer - Alcohol Analyzer Model: 8000 05/01/2026 Software: 8100.27 SN 80-000756 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>13:52</td></tr> <tr><td>Control Test</td><td>0.080</td><td>13:52</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:52</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:53</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:53</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:53</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:54</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0793</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7277</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	13:52	Control Test	0.080	13:52	Air Blank	0.000	13:52	Control Test	0.079	13:53	Air Blank	0.000	13:53	Control Test	0.079	13:53	Air Blank	0.000	13:54	Control Test Stats			Average	0.0793		Std Dev	0.0006		Rel Std Dev(%)	0.7277	
Test	g/210L	Time																																																																																																																																														
Air Blank	0.000	12:33																																																																																																																																														
Control Test	0.050	12:39																																																																																																																																														
Air Blank	0.000	12:40																																																																																																																																														
Control Test	0.051	12:41																																																																																																																																														
Air Blank	0.000	12:42																																																																																																																																														
Control Test	0.050	12:42																																																																																																																																														
Control Test Stats																																																																																																																																																
Average	0.0503																																																																																																																																															
Std Dev	0.0006																																																																																																																																															
Rel Std Dev(%)	1.1471																																																																																																																																															
Test	g/210L	Time																																																																																																																																														
Air Blank	0.000	12:59																																																																																																																																														
Control Test	0.079	13:00																																																																																																																																														
Air Blank	0.000	13:00																																																																																																																																														
Control Test	0.079	13:01																																																																																																																																														
Air Blank	0.000	13:02																																																																																																																																														
Control Test	0.081	13:02																																																																																																																																														
Air Blank	0.000	13:03																																																																																																																																														
Control Test Stats																																																																																																																																																
Average	0.0797																																																																																																																																															
Std Dev	0.0012																																																																																																																																															
Rel Std Dev(%)	1.4494																																																																																																																																															
Test	g/210L	Time																																																																																																																																														
Air Blank	0.000	12:47																																																																																																																																														
Control Test	0.198	12:49																																																																																																																																														
Air Blank	0.000	12:49																																																																																																																																														
Control Test	0.199	12:49																																																																																																																																														
Air Blank	0.000	12:49																																																																																																																																														
Control Test	0.199	12:50																																																																																																																																														
Air Blank	0.000	12:50																																																																																																																																														
Control Test Stats																																																																																																																																																
Average	0.1987																																																																																																																																															
Std Dev	0.0006																																																																																																																																															
Rel Std Dev(%)	0.2906																																																																																																																																															
Test	g/210L	Time																																																																																																																																														
Air Blank	0.000	13:52																																																																																																																																														
Control Test	0.080	13:52																																																																																																																																														
Air Blank	0.000	13:52																																																																																																																																														
Control Test	0.079	13:53																																																																																																																																														
Air Blank	0.000	13:53																																																																																																																																														
Control Test	0.079	13:53																																																																																																																																														
Air Blank	0.000	13:54																																																																																																																																														
Control Test Stats																																																																																																																																																
Average	0.0793																																																																																																																																															
Std Dev	0.0006																																																																																																																																															
Rel Std Dev(%)	0.7277																																																																																																																																															
<i>Isabella Smythowski</i> Operator's Signature	<i>Isabella Smythowski</i> Operator's Signature	<i>Isabella Smythowski</i> Operator's Signature	<i>Isabella Smythowski</i> Operator's Signature																																																																																																																																													

Cleared area and
reran analysis. IAS 5/11/26

HOLMES CITY MO
Intoxilyzer - Alcohol analyzer
Model 8000 SN 80-007396
12/01/2026
Software 8000.27

Test	9/20/26	Time
01- Blank	99%	12:52
02- Blank	99%	12:53
03- Blank	99%	12:53
04- Blank	99%	12:53
05- Blank	99%	12:53
06- Blank	99%	12:53
07- Blank	99%	12:53
08- Blank	99%	12:53
09- Blank	99%	12:53
10- Blank	99%	12:53
11- Blank	99%	12:53
12- Blank	99%	12:53
13- Blank	99%	12:53
14- Blank	99%	12:53
15- Blank	99%	12:53
16- Blank	99%	12:53
17- Blank	99%	12:53
18- Blank	99%	12:53
19- Blank	99%	12:53
20- Blank	99%	12:53
21- Blank	99%	12:53
22- Blank	99%	12:53
23- Blank	99%	12:53
24- Blank	99%	12:53
25- Blank	99%	12:53
26- Blank	99%	12:53
27- Blank	99%	12:53
28- Blank	99%	12:53
29- Blank	99%	12:53
30- Blank	99%	12:53
31- Blank	99%	12:53
32- Blank	99%	12:53
33- Blank	99%	12:53
34- Blank	99%	12:53
35- Blank	99%	12:53
36- Blank	99%	12:53
37- Blank	99%	12:53
38- Blank	99%	12:53
39- Blank	99%	12:53
40- Blank	99%	12:53
41- Blank	99%	12:53
42- Blank	99%	12:53
43- Blank	99%	12:53
44- Blank	99%	12:53
45- Blank	99%	12:53
46- Blank	99%	12:53
47- Blank	99%	12:53
48- Blank	99%	12:53
49- Blank	99%	12:53
50- Blank	99%	12:53
51- Blank	99%	12:53
52- Blank	99%	12:53
53- Blank	99%	12:53
54- Blank	99%	12:53
55- Blank	99%	12:53
56- Blank	99%	12:53
57- Blank	99%	12:53
58- Blank	99%	12:53
59- Blank	99%	12:53
60- Blank	99%	12:53
61- Blank	99%	12:53
62- Blank	99%	12:53
63- Blank	99%	12:53
64- Blank	99%	12:53
65- Blank	99%	12:53
66- Blank	99%	12:53
67- Blank	99%	12:53
68- Blank	99%	12:53
69- Blank	99%	12:53
70- Blank	99%	12:53
71- Blank	99%	12:53
72- Blank	99%	12:53
73- Blank	99%	12:53
74- Blank	99%	12:53
75- Blank	99%	12:53
76- Blank	99%	12:53
77- Blank	99%	12:53
78- Blank	99%	12:53
79- Blank	99%	12:53
80- Blank	99%	12:53
81- Blank	99%	12:53
82- Blank	99%	12:53
83- Blank	99%	12:53
84- Blank	99%	12:53
85- Blank	99%	12:53
86- Blank	99%	12:53
87- Blank	99%	12:53
88- Blank	99%	12:53
89- Blank	99%	12:53
90- Blank	99%	12:53
91- Blank	99%	12:53
92- Blank	99%	12:53
93- Blank	99%	12:53
94- Blank	99%	12:53
95- Blank	99%	12:53
96- Blank	99%	12:53
97- Blank	99%	12:53
98- Blank	99%	12:53
99- Blank	99%	12:53
100- Blank	99%	12:53

Abella Hughes
Operator's Signature

Stability Check 0.08 g/210 L: 80-000796

Florida Department of Law Enforcement

Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: HAINES CITY PD
Time of Inspection: 08:00

Date of Inspection: 05/04/2026

Serial Number: 80-000796
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK		No	Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000		No	Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No	Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Standard Deviations				
---------------------	--	--	--	--

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: _____ Number of Simulators Used: _____


Remarks:

DSP FAILNon-compliance: .

DSP failed pre-inspection diagnostic check. KTS 5/4/26

The above instrument complies () does not comply (X) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



KATIE T SPEARIN

Signature and Printed Name

05/04/2026

Date

Return Material Authorization

Ship to:☐ CMI, Inc.☒ Enforcement ElectronicsShipment to repair facility authorized by: Sean Bruner on 5/7/2026Items Returned: ☒ Instrument ☐ Supplies ☐ Other Describe: _____Instrument Model: Intoxilyzer 8000 Serial Number: 80-000796

Bill To Address:

Sean Burner

Haines City PD

35400 Highway 27, Haines City, FL 33850

Ship to Address:

Katie Spearin

FL Dept of Law Enforcement

813 B Lake Bradford Rd, FL 32304 Tallahassee

Reason for Return:

DSP Fail- DVM monitor measuring zeros for filters.

I require an estimate **BEFORE** any repairs will be authorized and/or conductedPlease contact: Sean Bruner Phone #: 727-641-2144Email: Sean.Bruner@hainescity.com

ATP Contact Name: Katie Spearin

ATP Email: KaitlynSpearin@fdle.state.fl.us