

Alcohol Testing Program - Instrument Processing Sheet

Agency: CITRUS COUNTY SO

Instrument Serial Number: 80-000819

Date In: 7/22/2026

DI Completion Date: N/A

☐ Ship

☐ P/U

☐ H/D

☒ CMI

☐ EE

Intake By: WKP Date: 7/22/2026		Quality Checks By: _____ Date: _____		Flow Adjustment By: _____																																							
☒ Annual ☒ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: Instrument dropped off by agency with DSP Failure during the instrument start-up diagnostics. WKP 7/22/2026		☐ Breath Tube Screen ☐ Replace External O-Rings ☐ Instrument Set Up Verified ☐ R-Value: _____ ☐ Flow Verification (L/s) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547) ☐ Barometric Pressure Check Gauge ID #: _____ Gauge: _____ Instrument: _____ ☐ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)																																							
						<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <th style="width:15%;">Simulator</th> <th style="width:20%;">Serial #</th> <th style="width:15%;">Lot#/Exp</th> <th style="width:50%;">Maintenance By: _____ Date: _____</th> </tr> <tr> <td>0.050</td> <td></td> <td></td> <td rowspan="5"> ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other: _____ </td> </tr> <tr><td>0.080</td><td></td><td></td></tr> <tr><td>0.200</td><td></td><td></td></tr> <tr><td>0.080 DGS</td><td>N/A</td><td></td></tr> <tr><td> </td><td></td><td></td></tr> </table>		Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____	0.050			☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other: _____	0.080			0.200			0.080 DGS	N/A																				
						Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____																																		
						0.050			☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other: _____																																		
						0.080																																					
0.200																																											
0.080 DGS	N/A																																										
Optical Bench Adjustment By: _____		Department Inspection By: _____																																									
Barometric Pressure Gauge: _____ ID#: _____		Barometric Pressure ID#: _____																																									
<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <th style="width:15%;">Simulator</th> <th style="width:20%;">Serial #</th> <th style="width:15%;">Lot #</th> <th style="width:50%;">Expiration</th> </tr> <tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr> <tr><td>0.040</td><td></td><td></td><td></td></tr> <tr><td>0.100</td><td></td><td></td><td></td></tr> <tr><td>0.200</td><td></td><td></td><td></td></tr> <tr><td>0.300</td><td></td><td></td><td></td></tr> <tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr> </table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Gauge: _____ Instrument: _____ Mouth Alcohol Solution Lot #: _____ Exp: _____ Acetone Stock Solution Lot #: _____ Exp: _____ <table border="1" style="width:100%; border-collapse: collapse;"> <tr> <th style="width:50%;">Simulator</th> <th style="width:50%;">Serial Number</th> </tr> <tr><td>0.000</td><td></td></tr> <tr><td>Interferent</td><td></td></tr> <tr><td>0.050</td><td></td></tr> <tr><td>0.080</td><td></td></tr> <tr><td>0.200</td><td></td></tr> </table>		Simulator	Serial Number	0.000		Interferent		0.050		0.080		0.200	
Simulator	Serial #	Lot #	Expiration																																								
0.000		N/A	N/A																																								
0.040																																											
0.100																																											
0.200																																											
0.300																																											
0.080 DGS	N/A																																										
Simulator	Serial Number																																										
0.000																																											
Interferent																																											
0.050																																											
0.080																																											
0.200																																											
☐ Post Optical Bench Adjustment Stability Checks		Attachments																																									
<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <th style="width:15%;">Simulator</th> <th style="width:20%;">Serial #</th> <th style="width:15%;">Lot #</th> <th style="width:50%;">Expiration</th> </tr> <tr><td>0.050</td><td></td><td></td><td></td></tr> <tr><td>0.080</td><td></td><td></td><td></td></tr> <tr><td>0.200</td><td></td><td></td><td></td></tr> <tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr> </table>		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			☐ Form 41 ☐ Post-Stability Checks ☐ Stability Checks ☐ Flow Adjustment ☐ Calibration Certificate ☐ Form 40 ☐ Optical Bench Adjustment ☒ Other: Return Material Authorization																					
Simulator	Serial #	Lot #	Expiration																																								
0.050																																											
0.080																																											
0.200																																											
0.080 DGS	N/A																																										
Gauge ID #: _____ Gauge: _____ Instrument: _____																																											
Notes/Suggested Service: Powered on instrument and the start up diagnostic check had a DSP Failure. 3 and 9 micron channels are reading zeros. Sending instrument to repair. WKP 7/22/2026		☐ Instrument Complies with Chapter 11D-8, FAC ☒ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use																																									
		LeAndra Higginbotham Digitally signed by LeAndra Higginbotham Date: 2026.08.03 17:13:03 -04'00' Kaitlyn Spearin Digitally signed by Kaitlyn Spearin Date: 2026.08.10 12:52:46 -04'00'																																									
		Tech Review Admin Review																																									

Return Material Authorization

Ship to:☒ CMI, Inc.☐ Enforcement ElectronicsShipment to repair facility authorized by: Laura Anstead on 7/22/2026Items Returned: ☒ Instrument ☐ Supplies ☐ Other Describe: Click or tap here to enter text.Instrument Model: Intoxilyzer 8000 Serial Number: 80-000819

Bill To Address:

Laura Anstead
2604 W. Woodland Ridge Dr.
Lecanto, FL, 34461

Ship to Address:

FDLE Offsite Mail Facility
Alcohol Testing Program
813-B Lake Bradford Road
Tallahassee, FL 32304

Reason for Return:

Instrument is DSP failing during instrument start up diagnostic test. The instrument 3 and 9 micron filters are reading '0'.

I require an estimate **BEFORE** any repairs will be authorized and/or conducted

Please contact: Laura Anstead Phone #: 352-424-5941Email: Lanstead@sheriffcitrus.org

ATP Contact Name: Wen-Chi Pierson

ATP Email: Wen-ChiPierson@fdle.state.fl.us

Alcohol Testing Program - Instrument Processing Sheet

Agency: CITRUS COUNTY SO

Instrument Serial Number: 80-000819

Date In: 2/4/2026

DI Completion Date: 2/9/2026

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: WKP Date: 2/4/2026		Quality Checks By: WKP Date: 2/6/2026		Flow Adjustment By: _____																													
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: 213 ☒ Flow Verification (L/s) Flow Column #: ATP 105 32 mm 0.160 (.139-.169) 36 mm 0.175 (.156-.190) 53 mm 0.242 (.228-.278) 103 mm 0.500 (.447-.547) ☒ Barometric Pressure Check Gauge ID #: 28421 Gauge: 1012 Instrument: 1010 ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm (.139-.169) 36 mm (.156-.190) 53 mm (.228-.278) 103 mm (.447-.547)																													
		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot#/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP5088</td> <td>202406K 6/19/2026</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> <td>202406L 6/19/2026</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> <td>202406N 6/20/2026</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>AG510701 4/17/2027</td> </tr> </tbody> </table>		Simulator	Serial #	Lot#/Exp	0.050	MP5088	202406K 6/19/2026	0.080	MP5089	202406L 6/19/2026	0.200	MP5090	202406N 6/20/2026	0.080 DGS	N/A	AG510701 4/17/2027	Maintenance By: _____ Date: _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Tank Sensor Tare ☐ Breath Tube Replacement ☐ Other:														
Simulator	Serial #	Lot#/Exp																															
0.050	MP5088	202406K 6/19/2026																															
0.080	MP5089	202406L 6/19/2026																															
0.200	MP5090	202406N 6/20/2026																															
0.080 DGS	N/A	AG510701 4/17/2027																															
Optical Bench Adjustment By: _____		Department Inspection By: WKP																															
Barometric Pressure Gauge: _____ ID#: _____		Barometric Pressure ID#: 34418																															
<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Gauge: 1023 Instrument: 1022 Mouth Alcohol Solution Lot #: 2025-D Exp: 9/25/2027 Acetone Stock Solution Lot #: 2025-B Exp: 9/22/2027			
Simulator	Serial #	Lot #	Expiration																														
0.000		N/A	N/A																														
0.040																																	
0.100																																	
0.200																																	
0.300																																	
0.080 DGS	N/A																																
☐ Post Optical Bench Adjustment Stability Checks		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>				Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A										
Simulator	Serial #	Lot #	Expiration																														
0.050																																	
0.080																																	
0.200																																	
0.080 DGS	N/A																																
Gauge ID #: _____		Gauge: _____ Instrument: _____																															
Notes/Suggested Service:		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																															
		Digitally signed by Taylor Gutschow Date: 2026.02.10 09:45:55 -05'00'																															
		Digitally signed by Kaitlyn Spearin Date: 2026.02.17 07:47:01 -05'00'																															
		Tech Review Admin Review																															

Stability Checks

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083
✓ CITRUS COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 02/15/2025 Software: 8100.27 Test g/210L Time Air Blank 0.000 14:03 Control Test 0.048 14:04 Air Blank 0.000 14:04 Control Test 0.047 14:05 Air Blank 0.000 14:06 Control Test 0.047 14:06 Air Blank 0.000 14:07 Control Test Stats Average 0.0473 Std Dev 0.0006 Rel Std Dev(%) 1.2198	✓ CITRUS COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 02/15/2025 Software: 8100.27 Test g/210L Time Air Blank 0.000 14:15 Control Test 0.077 14:15 Air Blank 0.000 14:16 Control Test 0.077 14:16 Air Blank 0.000 14:17 Control Test 0.077 14:18 Air Blank 0.000 14:19 Control Test Stats Average 0.0770 Std Dev 0.0000 Rel Std Dev(%) 0.0000	✓ CITRUS COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 02/15/2025 Software: 8100.27 Test g/210L Time Air Blank 0.000 14:24 Control Test 0.198 14:24 Air Blank 0.000 14:25 Control Test 0.197 14:25 Air Blank 0.000 14:26 Control Test 0.198 14:26 Air Blank 0.000 14:27 Control Test Stats Average 0.1977 Std Dev 0.0006 Rel Std Dev(%) 0.2921	✓ CITRUS COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 02/15/2025 Software: 8100.27 Test g/210L Time Air Blank 0.000 13:57 Control Test 0.080 13:58 Air Blank 0.000 13:58 Control Test 0.080 13:58 Air Blank 0.000 13:59 Control Test 0.080 13:59 Air Blank 0.000 14:00 Control Test Stats Average 0.0800 Std Dev 0.0000 Rel Std Dev(%) 0.0000
<i>ment</i> Operator's Signature	<i>ment</i> Operator's Signature	<i>ment</i> Operator's Signature	<i>ment</i> Operator's Signature

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: CITRUS COUNTY SO
Time of Inspection: 12:32

Date of Inspection: 02/09/2026

Serial Number: 80-000819
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202406K Exp: 06/19/2026	0.08g/210L Test (g/210L) Lot#:202406L Exp: 06/19/2026	0.20g/210L Test (g/210L) Lot#:202406N Exp: 06/20/2026	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG510701 Exp: 04/17/2027
0.000	0.048	0.078	0.200	0.080
0.000	0.048	0.078	0.200	0.080
0.000	0.048	0.078	0.200	0.080
0.000	0.048	0.079	0.200	0.080
0.000	0.048	0.078	0.200	0.080
0.000	0.048	0.079	0.200	0.080
0.000	0.048	0.079	0.200	0.080
0.000	0.049	0.079	0.200	0.080
0.000	0.048	0.079	0.200	0.080
0.000	0.049	0.079	0.200	0.080

Standard Deviations	0.0004	0.0005	0.0000	0.0000
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0002 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



WEN-CHI K PIERSON
Signature and Printed Name

02/09/2026
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000819, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000819</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>CITRUS COUNTY SO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>02/09/2026</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>12:32</u>	0.200 g/ 210 L 0.008
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Wen-Chi

Chi Pierson

Date: 2026.02.10
09:26:23 -05'00'

Pierson

02/09/2026

Date

WEN-CHI K PIERSON,

Department Inspector

FDLE/ATP Form 69 January 2026

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality