

INSTRUMENT PROCESSING SHEET

Agency Lake County SO S/N 80-000823
 Date In 2/29/16 Date Out 3/8/16 Ship P/U H/D CMI EE

Intake Performed By TP

Registration
 Annual
 Return from CMI
 Return from Enforcement Electronics
 Other _____

Visual Inspection:
OK Case OK Handle
OK Dry Gas Holder OK Feet
OK Keyboard/Plug OK Back/Plugs
OK Screws tight OK Breath Hose

Other Equipment:
 Power cord
 Printer Cable
 Other Static Bag

Notes: _____

Quality Checks Performed By DUB

Breath Tube Screen
 Replace O-Rings
 Instrument Set Up Verified
 R-Value 196
 Flow Verification (L/s)
 Flow Column # ATP103
 32mm 0.148 (.139 - .169)
 36mm 0.164 (.156 - .190)
 53mm 0.234 (.228 - .278)
 103mm 0.496 (.447 - .547)

Barometric Pressure Check
 Gauge ID # 28427

Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	SD1018	201507A 7/14/17
0.08	SD1011	201502G 2/24/17
0.20	G4444	201505A 5/12/17
0.08 DGS	N/A	AGL05301 2/22/18

Flow Calibration Performed By DUB

Flow Calibration N/A
 Flow Calibration Complete
 Flow Column # _____
 5L/min - 17mm
 15L/min - 53mm
 30L/min - 103mm
 R-Value _____
 Post Calibration Verification (L/s)
 Flow Column # _____
 32mm _____ (.139 - .169)
 36mm _____ (.156 - .190)
 53mm _____ (.228 - .278)
 103mm _____ (.447 - .547)

Maintenance Performed By _____

Battery Replacement
 Dry Gas Regulator Replacement
 Breath Tube Replacement
 Other _____

Suggested Service

RECEIVED
 MAR 09 2016
 FDLE
 Alcohol Testing Program

Optical Bench Calibration Performed By _____

Optical Bench Calibration N/A
 Optical Bench Calibration Complete

Barometric Pressure Gauge ID # _____

Simulator	Serial Number	Lot Number	Expiration
0.000		N/A	N/A
0.040			
0.100			
0.200			
0.400			
0.080 DGS	N/A		

Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05			
0.08			
0.20			
0.08 DGS	N/A		

Department Inspection Performed By DUB

Barometric Pressure 1022 Gauge
 ID# 28427 1021 Instrument

Mouth Alcohol Solution Lot # 2015-A
 Acetone Stock Solution Lot # 2016-B

Simulator	Serial Number
0.00	SD1022
Interferent	SD1021
0.05	SD1018
0.08	SD1011
0.20	G4444

Attachments

Form 41
 Pre-Stability Tests
 Flow Calibration
 Optical Bench Cal
 Post-Stability Tests
 Other _____

Notes: QC-13K

Instrument Complies with Chapter 11D-8, FAC
 Instrument Does Not Comply with Chapter 11D-8, FAC
 Return to/Place into Evidentiary Use
 Remain Out of Evidentiary Use
 Conduct an Agency Inspection Before Evidentiary Use

Patrick Murphy

Quality Control Review

3/9/16

Date

Stability Checks 80-000823 Lake County SO. 3/8/16 RMB

RMB

RMB

LAKE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000823
03/08/2016
Software: 8100.27

LAKE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000823
03/08/2016
Software: 8100.27

LAKE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000823
03/08/2016
Software: 8100.27

LAKE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-000823
03/08/2016
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:42
Control Test	0.051	12:42
Air Blank	0.000	12:43
Control Test	0.050	12:43
Air Blank	0.000	12:44
Control Test	0.050	12:45
Air Blank	0.000	12:45
Control Test Stats		
Average	0.0503	
Std Dev	0.0006	
Rel Std Dev(%)	1.1471	

Test	g/210L	Time
Air Blank	0.000	12:50
Control Test	0.080	12:50
Air Blank	0.000	12:51
Control Test	0.080	12:51
Air Blank	0.000	12:52
Control Test	0.080	12:53
Air Blank	0.000	12:53
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Test	g/210L	Time
Air Blank	0.000	12:33
Control Test	0.198	12:34
Air Blank	0.000	12:34
Control Test	0.200	12:35
Air Blank	0.000	12:36
Control Test	0.200	12:36
Air Blank	0.000	12:37
Control Test Stats		
Average	0.1933	
Std Dev	0.0012	
Rel Std Dev(%)	0.5793	

Test	g/210L	Time
Air Blank	0.000	12:46
Control Test	0.081	12:46
Air Blank	0.000	12:47
Control Test	0.081	12:47
Air Blank	0.000	12:48
Control Test	0.081	12:48
Air Blank	0.000	12:48
Control Test Stats		
Average	0.0810	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

BK

RMB

Operator's Signature

RMB

Operator's Signature

RMB

Operator's Signature

RMB

Operator's Signature