

## INSTRUMENT PROCESSING SHEET

Agency New Smyrna Beach PD S/N 80-002231

Date In 7/1/16 Date Out 7/6/16

Ship  P/U  H/D  CMI  EE

RECEIVED  
JUL 06 2016

Alcohol Testing Program  
FDLE

**Intake** Performed By AP

Registration  
 Annual  
 Return from CMI  
 Return from Enforcement Electronics  
 Other \_\_\_\_\_

Visual Inspection:  
 Case  Handle  
 Dry Gas Holder  Feet  
 Keyboard/Plug  Back/Plugs  
 Screws tight  Breath Hose

Other Equipment:  
 Power cord  
 Printer Cable  
 Other Static Bag

Notes: \_\_\_\_\_

**Quality Checks** Performed By AP

Breath Tube Screen  
 Replace O-Rings  
 Instrument Set Up Verified  
 R-Value 225  
 Flow Verification (L/s)  
 Flow Column # A-7P102  
 32mm 150 (.139 - .169)  
 36mm 177 (.156 - .190)  
 53mm 246 (.228 - .278)  
 103mm 511 (.447 - .547)

Barometric Pressure Check  
 Gauge ID # 28427

Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	SD1018	201507A 7/14/17
0.08	SD1011	201601P 1/26/18
0.20	SD1025	201505A 5/12/17
0.08 DGS	N/A	A660530 2/12/18

**Flow Calibration** Performed By \_\_\_\_\_

Flow Calibration N/A  
 Flow Calibration Complete  
 Flow Column # \_\_\_\_\_  
 5L/min - 17mm  
 15L/min - 53mm  
 30L/min - 103mm

R-Value \_\_\_\_\_  
 Post Calibration Verification (L/s)  
 Flow Column # \_\_\_\_\_  
 32mm \_\_\_\_\_ (.139 - .169)  
 36mm \_\_\_\_\_ (.156 - .190)  
 53mm \_\_\_\_\_ (.228 - .278)  
 103mm \_\_\_\_\_ (.447 - .547)

**Maintenance** Performed By \_\_\_\_\_

Battery Replacement  
 Dry Gas Regulator Replacement  
 Breath Tube Replacement  
 Other \_\_\_\_\_

**Suggested Service**

**Optical Bench Calibration** Performed By \_\_\_\_\_

Optical Bench Calibration N/A  
 Optical Bench Calibration Complete

**Barometric Pressure Gauge** ID # \_\_\_\_\_

Simulator	Serial Number	Lot Number	Expiration
0.000		N/A	N/A
0.040			
0.100			
0.200			
0.400			
0.080 DGS	N/A		

Post Calibration Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05			
0.08			
0.20			
0.08 DGS	N/A		

**Department Inspection** Performed By AP

Barometric Pressure  
 ID# 28427 Gauge 1019  
 Instrument 1014

Mouth Alcohol Solution Lot # 2016-A  
 Acetone Stock Solution Lot # 2016-B

Simulator	Serial Number
0.00	SD1019
Interferent	SD1021
0.05	SD1018
0.08	SD1011
0.20	SD1025

**Attachments**

Form 41  
 Pre-Stability Tests  
 Flow Calibration  
 Optical Bench Cal  
 Post-Stability Tests  
 Other Form 40 & ATL Paperwork

Notes: AI received AMBIENT FAIL Messages. Took instrument out of service and sent to ATP. Ran Stabilities and received no Ambient fail messages. Department inspection conducted in compliance.  
- Advise AI to check room conditions.  
Breath Rebalanced

Instrument Complies with Chapter 11D-8, FAC  
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place Into Evidentiary Use  
 Remain Out of Evidentiary Use  
 Conduct an Agency Inspection Before Evidentiary Use

Quality Control Review

Date

GC-PWS

7/6/16

80-002231  
 Stambli, Mycles, Kess  
 7/6/16

INTOXILYZER 8000  
 Instrument Initialization  
 07:14 07/06/2016

NEW SMYRNA BEACH PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-002231  
 07/06/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	07:55
Control Test	0.078	07:55
Air Blank	0.000	07:55
Control Test	0.079	07:56
Air Blank	0.000	07:56
Control Test	0.079	07:57
Air Blank	0.000	07:57
Control Test Stats		
Average	0.0787	
Std Dev	0.0006	
Rel Std Dev(%)	0.7339	

DGS

Operator's Signature

NEW SMYRNA BEACH PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-002231  
 07/06/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:07
Control Test	0.078	08:08
Air Blank	0.000	08:09
Control Test	0.079	08:09
Air Blank	0.000	08:10
Control Test	0.078	08:11
Air Blank	0.000	08:11
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Operator's Signature

NEW SMYRNA BEACH PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-002231  
 07/06/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:12
Control Test	0.048	08:13
Air Blank	0.000	08:13
Control Test	0.048	08:14
Air Blank	0.000	08:14
Control Test	0.048	08:15
Air Blank	0.000	08:15
Control Test Stats		
Average	0.0480	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Operator's Signature

NEW SMYRNA BEACH PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-002231  
 07/06/2016  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:17
Control Test	0.198	08:17
Air Blank	0.000	08:18
Control Test	0.198	08:18
Air Blank	0.000	08:19
Control Test	0.199	08:20
Air Blank	0.000	08:20
Control Test Stats		
Average	0.1983	
Std Dev	0.0006	
Rel Std Dev(%)	0.2911	

MS  
 MSK

Operator's Signature