



Alcohol Testing Program

**INSTRUMENT PROCESSING SHEET**

Agency Winter Park PD S/N 80-000994  
Date In 7/27/17 Date Out 9/26/17  Ship  P/U  H/D  CMI  EE

**Intake** Performed By [Signature]

Registration  
 Annual  
 Return from CMI  
 Return from Enforcement  
 Electronics  
 Other \_\_\_\_\_

Visual Inspection:  
OK Case OK Handle  
OK Dry Gas Holder OK Feet  
OK Keyboard/Plug OK Back/Plugs  
OK Screws tight OK Breath Hose

Other Equipment:  
 Power cord  
 Printer Cable  
 Other Static Bag

Notes: \_\_\_\_\_

**Quality Checks** Performed By [Signature]

Lab Temp °C 21.9  
 Breath Tube Screen  
 Replace External O-Rings  
 Instrument Set Up Verified  
 R-Value 195  
 Flow Verification (L/s) ATP103  
 Flow Column # ATP102  
 32mm 143 (.139 - .169)  
 36mm 160 (.156 - .190)  
 53mm 235 (.228 - .278)  
 103mm 500 (.447 - .547)

Barometric Pressure Check  
 Gauge ID # 26932

Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	<u>SD3962</u>	<u>201707D</u> <u>7/25/19</u>
0.08	<u>SD1013</u>	<u>201707E</u> <u>7/25/19</u>
0.20	<u>SD2935</u> <u>DR3856</u>	<u>201707C</u> <u>7/24/19</u>
0.08 DGS	N/A	<u>AG617605</u> <u>7/14/18</u>

**Flow Calibration** Performed By \_\_\_\_\_

Flow Calibration N/A  
 Flow Calibration Complete  
 Flow Column # \_\_\_\_\_  
 5L/min - 17mm  
 15L/min - 53mm  
 30L/min - 103mm

R-Value \_\_\_\_\_  
 Post Calibration Verification (L/s)  
 Flow Column # \_\_\_\_\_  
 32mm \_\_\_\_\_ (.139 - .169)  
 36mm \_\_\_\_\_ (.156 - .190)  
 53mm \_\_\_\_\_ (.228 - .278)  
 103mm \_\_\_\_\_ (.447 - .547)

**Maintenance** Performed By [Signature]

Battery Replacement  
 Dry Gas Regulator Replacement  
 Breath Tube Replacement  
 Other \_\_\_\_\_

**Quality Checks Cont.** Performed By [Signature]

**Simulator Temperatures °C**  
 External Digital Therm. ID #: 300505  
 34°C +/- .2 Serial #: SD3962  
 34°C +/- .2 Serial #: SD1013  
 34°C +/- .2 Serial #: DR3856

**Calibration Adjustment** Performed By \_\_\_\_\_

Calibration Adjustment N/A  
 Calibration Adjustment Complete  
 Barometric Pressure Gauge \_\_\_\_\_ ID # \_\_\_\_\_

Simulator	Serial Number	Lot Number	Expiration
0.000		N/A	N/A
0.040			
0.100			
0.200			
0.300			
0.080 DGS	N/A		

Post Calibration Adjustment Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05			
0.08			
0.20			
0.08 DGS	N/A		

**Department Inspection** Performed By [Signature]

Barometric Pressure 1012 Gauge  
 ID# 26932 1012 Instrument

Mouth Alcohol Solution Lot # 2016-C  
 Acetone Stock Solution Lot # 2017-A

Simulator	Serial Number
0.00	<u>SD1019</u>
Interferent	<u>SD1021</u>
0.05	<u>SD3962</u>
0.08	<u>SD1013</u>
0.20	<u>SD2935</u> <u>DR3856</u>

**Attachments**

Form 41  
 Pre-Stability Tests  
 Flow Calibration

Calibration Adjustment  
 Post-Stability Tests  
 Other Form 40

Notes/Suggested Service:  
QA/QC 9/26/17 [Signature]

[Signature] 9/26/17  
 Date

Instrument Complies with Chapter 11D-8, FAC  
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use  
 Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

Quality Control Review

Date

80-000994  
 Stability checks  
 9/26/17

INTOXILYZER 8000  
 Instrument Initialization  
 08:17 09/26/2017

WINTER PARK P.D.  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-000994  
 09/26/2017  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:51
Control Test	0.080	08:52
Air Blank	0.000	08:52
Control Test	0.079	08:52
Air Blank	0.000	08:53
Control Test	0.080	08:53
Air Blank	0.000	08:54
Control Test Stats		
Average	0.0797	
Std Dev	0.0006	
Rel Std Dev(%)	0.7247	

DGS

Operator's Signature

WINTER PARK P.D.  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-000994  
 09/26/2017  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	08:55
Control Test	0.048	08:55
Air Blank	0.000	08:56
Control Test	0.048	08:57
Air Blank	0.000	08:57
Control Test	0.048	08:58
Air Blank	0.000	08:58
Control Test Stats		
Average	0.0480	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Operator's Signature

WINTER PARK P.D.  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-000994  
 09/26/2017  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:01
Control Test	0.197	09:01
Air Blank	0.000	09:02
Control Test	0.197	09:03
Air Blank	0.000	09:03
Control Test	0.197	09:04
Air Blank	0.000	09:04
Control Test Stats		
Average	0.1970	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Operator's Signature

WINTER PARK P.D.  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-000994  
 09/26/2017  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:06
Control Test	0.079	09:06
Air Blank	0.000	09:07
Control Test	0.079	09:08
Air Blank	0.000	09:08
Control Test	0.079	09:09
Air Blank	0.000	09:09
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Operator's Signature

9/26/17  
 [Signature]