



Alcohol Testing Program

### INSTRUMENT PROCESSING SHEET

Agency Indian Harbour Beach, PD S/N 80-001263

Date In 9/25/17

Date Out 9/27/17

Ship  P/U  H/D  CMI  EE

**Intake** Performed By TS

Registration  
 Annual  
 Return from CMI  
 Return from Enforcement Electronics  
 Other \_\_\_\_\_

Visual Inspection:  
 Case  Handle  
 Dry Gas Holder  Feet  
 Keyboard/Plug  Back/Plugs  
 Screws tight  Breath Hose

Other Equipment:  
 Power cord  
 Printer Cable  
 Other Static Bag

Notes: \_\_\_\_\_

**Quality Checks** Performed By SP

Lab Temp °C 21.9  
 Breath Tube Screen  
 Replace External O-Rings  
 Instrument Set Up Verified  
 R-Value 213  
 Flow Verification (L/s)  
 Flow Column # ATP103  
 32mm .144 (.139 - .169)  
 36mm .164 (.156 - .190)  
 53mm .234 (.228 - .278)  
 103mm .496 (.447 - .547)

Barometric Pressure Check  
 Gauge ID # 28427

Stability Checks

Simulator	Serial #	Lot #/Exp
0.05	<u>G2078</u>	<u>201707D</u> <u>7-25-19</u>
0.08	<u>G8149</u>	<u>201707E</u> <u>7-25-19</u>
0.20	<u>G2408</u>	<u>201707C</u> <u>7-24-19</u>
0.08 DGS	N/A	<u>AG626604</u> <u>9-22-18</u>

**Flow Calibration** Performed By \_\_\_\_\_

Flow Calibration N/A  
 Flow Calibration Complete  
 Flow Column # \_\_\_\_\_  
 5L/min - 17mm  
 15L/min - 53mm  
 30L/min - 103mm  
 R-Value \_\_\_\_\_  
 Post Calibration Verification (L/s)  
 Flow Column # \_\_\_\_\_  
 32mm \_\_\_\_\_ (.139 - .169)  
 36mm \_\_\_\_\_ (.156 - .190)  
 53mm \_\_\_\_\_ (.228 - .278)  
 103mm \_\_\_\_\_ (.447 - .547)

**Maintenance** Performed By \_\_\_\_\_

Battery Replacement  
 Dry Gas Regulator Replacement  
 Breath Tube Replacement  
 Other \_\_\_\_\_

**Quality Checks Cont.** Performed By SP  
**Simulator Temperatures °C**  
 External Digital Therm. ID#: 300504  
 34°C +/-2 Serial #: G2078  
 34°C +/-2 Serial #: G8149  
 34°C +/-2 Serial #: G2408

**Calibration Adjustment** Performed By \_\_\_\_\_

Calibration Adjustment N/A  
 Calibration Adjustment Complete

Barometric Pressure Gauge \_\_\_\_\_ ID # \_\_\_\_\_

Simulator	Serial Number	Lot Number	Expiration
0.000		N/A	N/A
0.040			
0.100			
0.200			
0.300			
0.080 DGS	N/A		

Post Calibration Adjustment Stability Checks

Simulator	Serial Number	Lot Number	Expiration
0.05			
0.08			
0.20			
0.08 DGS	N/A		

**Department Inspection** Performed By SP

Barometric Pressure 1010 Gauge  
 ID# 28427 1010 Instrument

Mouth Alcohol Solution Lot # 2016-C  
 Acetone Stock Solution Lot # 2017-A

Simulator	Serial Number
0.00	<u>G2879</u>
Interferent	<u>G8144</u>
0.05	<u>G2078</u>
0.08	<u>G8149</u>
0.20	<u>G2408</u>

**Attachments**

Form 41  
 Pre-Stability Tests  
 Flow Calibration  
 Calibration Adjustment  
 Post-Stability Tests  
 Other \_\_\_\_\_

Notes/Suggested Service:  
QA/QC by SP 9/27/17

[Signature]  
 Quality Control Review

9/27/17  
 Date

Instrument Complies with Chapter 11D-8, FAC  
 Instrument Does Not Comply with Chapter 11D-8, FAC

Return to/Place into Evidentiary Use  
 Remain Out of Evidentiary Use

Conduct an Agency Inspection Before Evidentiary Use

RECEIVED  
FDLE  
Alcohol Testing Program  
SEP 27 2017

# STABILITY CHECKS #80-001263 - INDIAN HARBOUR PD.

INDIAN HARBOUR PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001263  
 09/27/2017  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:32
Control Test	0.049	09:33
Air Blank	0.000	09:34
Control Test	0.049	09:34
Air Blank	0.000	09:35
Control Test	0.049	09:35
Air Blank	0.000	09:36
Control Test	0.000	09:36
Control Test Stats		
Average	0.0490	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP

Operator's Signature

SP  
 1/12/16  
 DR

INDIAN HARBOUR PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001263  
 09/27/2017  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:13
Control Test	0.080	09:14
Air Blank	0.000	09:15
Control Test	0.080	09:15
Air Blank	0.000	09:16
Control Test	0.080	09:16
Air Blank	0.000	09:17
Control Test	0.000	09:17
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

SP

Operator's Signature

INDIAN HARBOUR PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001263  
 09/27/2017  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:18
Control Test	0.200	09:19
Air Blank	0.000	09:19
Control Test	0.200	09:20
Air Blank	0.000	09:21
Control Test	0.199	09:21
Air Blank	0.000	09:22
Control Test	0.000	09:22
Control Test Stats		
Average	0.1997	
Std Dev	0.0006	
Rel Std Dev(%)	0.2892	

SP

Operator's Signature

INDIAN HARBOUR PD  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-001263  
 09/27/2017  
 Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	09:24
Control Test	0.080	09:24
Air Blank	0.000	09:25
Control Test	0.080	09:25
Air Blank	0.000	09:26
Control Test	0.080	09:26
Air Blank	0.000	09:27
Control Test	0.000	09:27
Control Test Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

DGS

SP

Operator's Signature

SP