



INSTRUMENT PROCESSING SHEET

Agency Charlotte County SO

S/N 80-007083

Florida Department of
Law Enforcement

Date In 7/26/2021 DI Completion Date 8/2/2021

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By MH	Quality Checks	By MH	Date 07/28/2021	Flow Calibration	By	Date																																																										
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: AI Note: <u>DGS Outside Tolerance</u>		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>245</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP 106</u> 32 mm <u>0.156</u> (.139 - .169) 36 mm <u>0.171</u> (.156 - .190) 53 mm <u>0.242</u> (.228 - .278) 103 mm <u>0.496</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28199</u> ☒ Stability Checks			Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																																												
		<table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP5092</td> <td>202010A 10/05/2022</td> </tr> <tr> <td>0.080</td> <td>MP5093</td> <td>202010B 10/05/2022</td> </tr> <tr> <td>0.200</td> <td>MP5097</td> <td>202010D 10/06/2022</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>AG026705 09/23/2022</td> </tr> </tbody> </table>		Simulator	Serial #	Lot #/Exp	0.050	MP5092	202010A 10/05/2022	0.080	MP5093	202010B 10/05/2022	0.200	MP5097	202010D 10/06/2022	0.080 DGS	N/A	AG026705 09/23/2022	Maintenance By _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ DI Temp. Checks By MH _____ ☒ Lab Temp °C <u>21.56</u> External Digital Therm. ID#: <u>300503</u> ☒ 34°C +/- .2 Serial #: <u>MP5092</u> ☒ 34°C +/- .2 Serial #: <u>MP5093</u> ☒ 34°C +/- .2 Serial #: <u>MP5097</u>																																														
Simulator	Serial #	Lot #/Exp																																																															
0.050	MP5092	202010A 10/05/2022																																																															
0.080	MP5093	202010B 10/05/2022																																																															
0.200	MP5097	202010D 10/06/2022																																																															
0.080 DGS	N/A	AG026705 09/23/2022																																																															
Calibration Adjustment By MH _____ Barometric Pressure Gauge <u>1017</u> ID # <u>28663</u> <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>MP5095</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td>MP5096</td> <td>20060</td> <td>02/10/2022</td> </tr> <tr> <td>0.100</td> <td>MP5098</td> <td>20420</td> <td>09/22/2022</td> </tr> <tr> <td>0.200</td> <td>MP5100</td> <td>20160</td> <td>03/18/2022</td> </tr> <tr> <td>0.300</td> <td>MP5101</td> <td>20030</td> <td>01/21/2022</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>00521080A2</td> <td>02/05/2023</td> </tr> </tbody> </table> ☒ Post Calibration Adjustment Stability Checks <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP5092</td> <td>202010A</td> <td>10/05/2022</td> </tr> <tr> <td>0.080</td> <td>MP5093</td> <td>202010B</td> <td>10/05/2022</td> </tr> <tr> <td>0.200</td> <td>MP5097</td> <td>202010D</td> <td>10/06/2022</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>AG026705</td> <td>09/23/2022</td> </tr> </tbody> </table>				Simulator	Serial #	Lot #	Expiration	0.000	MP5095	N/A	N/A	0.040	MP5096	20060	02/10/2022	0.100	MP5098	20420	09/22/2022	0.200	MP5100	20160	03/18/2022	0.300	MP5101	20030	01/21/2022	0.080 DGS	N/A	00521080A2	02/05/2023	Simulator	Serial #	Lot #	Expiration	0.050	MP5092	202010A	10/05/2022	0.080	MP5093	202010B	10/05/2022	0.200	MP5097	202010D	10/06/2022	0.080 DGS	N/A	AG026705	09/23/2022	Department Inspection By MH _____ Barometric Pressure ID# <u>68639</u> Gauge <u>1017</u> Instrument <u>1017</u> Mouth Alcohol Solution Lot # <u>2021-B</u> Acetone Stock Solution Lot # <u>2020-A</u> <table border="1"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>SD1014</td> </tr> <tr> <td>Interferent</td> <td>SD1015</td> </tr> <tr> <td>0.050</td> <td>MP5092</td> </tr> <tr> <td>0.080</td> <td>MP5093</td> </tr> <tr> <td>0.200</td> <td>MP5097</td> </tr> </tbody> </table>		Simulator	Serial Number	0.000	SD1014	Interferent	SD1015	0.050	MP5092	0.080	MP5093	0.200	MP5097
Simulator	Serial #	Lot #	Expiration																																																														
0.000	MP5095	N/A	N/A																																																														
0.040	MP5096	20060	02/10/2022																																																														
0.100	MP5098	20420	09/22/2022																																																														
0.200	MP5100	20160	03/18/2022																																																														
0.300	MP5101	20030	01/21/2022																																																														
0.080 DGS	N/A	00521080A2	02/05/2023																																																														
Simulator	Serial #	Lot #	Expiration																																																														
0.050	MP5092	202010A	10/05/2022																																																														
0.080	MP5093	202010B	10/05/2022																																																														
0.200	MP5097	202010D	10/06/2022																																																														
0.080 DGS	N/A	AG026705	09/23/2022																																																														
Simulator	Serial Number																																																																
0.000	SD1014																																																																
Interferent	SD1015																																																																
0.050	MP5092																																																																
0.080	MP5093																																																																
0.200	MP5097																																																																
Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☒ Calibration Adjustment ☐ Post-Stability Checks ☐ Flow Calibration ☒ Form 40 ☒ Other <u>AI Note</u>				☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place Into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use 2021.08.																																																													
Notes/Suggested Service: _____ _____ _____ _____ _____ _____ _____ _____				Taylor Gutschow Digitally signed by Taylor Gutschow Date: 2021.08.03 11:48:04 -0400 04 09:59:57 Tech Review / Date Admin Review / Date																																																													

Florida Department of Law Enforcement

Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: CHARLOTTE COUNTY SO
Time of Inspection: 11:21

Date of Inspection: 08/02/2021

Serial Number: 80-007083
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202010A Exp: 10/05/2022	0.08g/210L Test (g/210L) Lot#:202010B Exp: 10/05/2022	0.20g/210L Test (g/210L) Lot#:202010D Exp: 10/06/2022	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG026705 Exp: 09/23/2022
0.000	0.050	0.079	0.197	0.080
0.000	0.051	0.079	0.198	0.080
0.000	0.050	0.079	0.199	0.079
0.000	0.050	0.079	0.198	0.079
0.000	0.051	0.080	0.199	0.079
0.000	0.050	0.079	0.199	0.079
0.000	0.050	0.079	0.198	0.079
0.000	0.051	0.080	0.199	0.079
0.000	0.051	0.080	0.199	0.079
0.000	0.051	0.080	0.198	0.079


Standard Deviations	0.0005	0.0005	0.0006	0.0004
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0005 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



Signature and Printed Name

MICHAEL D HAUGHEY

08/02/2021
Date

TYPE OF TEST	SERIAL NUMBER	AGENCY	DATE	PERFORMED BY
Stabilities	80-007083	Charlotte County SO	07/28/21	MK

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L																																																																																																																																																
0.047 to 0.053 ☒	0.077 to 0.083 ☒	0.194 to 0.206 ☒	0.077 to 0.083 ☒																																																																																																																																																
CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model: 8000 07/28/2021 Software: 8100.27 <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>09:17</td></tr><tr><td>Control Test</td><td>0.046</td><td>09:18</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:19</td></tr><tr><td>Control Test</td><td>0.046</td><td>09:19</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:20</td></tr><tr><td>Control Test</td><td>0.047</td><td>09:20</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:21</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0463</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>1.2461</td><td></td></tr></table> MK Operator's Signature	Test	g/210L	Time	Air Blank	0.000	09:17	Control Test	0.046	09:18	Air Blank	0.000	09:19	Control Test	0.046	09:19	Air Blank	0.000	09:20	Control Test	0.047	09:20	Air Blank	0.000	09:21	Control Test Stats			Average	0.0463		Std Dev	0.0006		Rel Std Dev(%)	1.2461		CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model: 8000 07/28/2021 Software: 8100.27 <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>09:22</td></tr><tr><td>Control Test</td><td>0.077</td><td>09:23</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:23</td></tr><tr><td>Control Test</td><td>0.076</td><td>09:24</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:25</td></tr><tr><td>Control Test</td><td>0.076</td><td>09:25</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:26</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0763</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7564</td><td></td></tr></table> MK Operator's Signature	Test	g/210L	Time	Air Blank	0.000	09:22	Control Test	0.077	09:23	Air Blank	0.000	09:23	Control Test	0.076	09:24	Air Blank	0.000	09:25	Control Test	0.076	09:25	Air Blank	0.000	09:26	Control Test Stats			Average	0.0763		Std Dev	0.0006		Rel Std Dev(%)	0.7564		CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model: 8000 07/28/2021 Software: 8100.27 <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>09:30</td></tr><tr><td>Control Test</td><td>0.193</td><td>09:31</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:32</td></tr><tr><td>Control Test</td><td>0.194</td><td>09:32</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:33</td></tr><tr><td>Control Test</td><td>0.195</td><td>09:34</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:34</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.1940</td><td></td></tr><tr><td>Std Dev</td><td>0.0910</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.5155</td><td></td></tr></table> MK Operator's Signature	Test	g/210L	Time	Air Blank	0.000	09:30	Control Test	0.193	09:31	Air Blank	0.000	09:32	Control Test	0.194	09:32	Air Blank	0.000	09:33	Control Test	0.195	09:34	Air Blank	0.000	09:34	Control Test Stats			Average	0.1940		Std Dev	0.0910		Rel Std Dev(%)	0.5155		CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model: 8000 07/28/2021 Software: 8100.27 <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>09:35</td></tr><tr><td>Control Test</td><td>0.077</td><td>09:35</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:36</td></tr><tr><td>Control Test</td><td>0.077</td><td>09:36</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:37</td></tr><tr><td>Control Test</td><td>0.077</td><td>09:37</td></tr><tr><td>Air Blank</td><td>0.000</td><td>09:38</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0770</td><td></td></tr><tr><td>Std Dev</td><td>0.0000</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr></table> MK Operator's Signature	Test	g/210L	Time	Air Blank	0.000	09:35	Control Test	0.077	09:35	Air Blank	0.000	09:36	Control Test	0.077	09:36	Air Blank	0.000	09:37	Control Test	0.077	09:37	Air Blank	0.000	09:38	Control Test Stats			Average	0.0770		Std Dev	0.0000		Rel Std Dev(%)	0.0000	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:17																																																																																																																																																	
Control Test	0.046	09:18																																																																																																																																																	
Air Blank	0.000	09:19																																																																																																																																																	
Control Test	0.046	09:19																																																																																																																																																	
Air Blank	0.000	09:20																																																																																																																																																	
Control Test	0.047	09:20																																																																																																																																																	
Air Blank	0.000	09:21																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0463																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	1.2461																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:22																																																																																																																																																	
Control Test	0.077	09:23																																																																																																																																																	
Air Blank	0.000	09:23																																																																																																																																																	
Control Test	0.076	09:24																																																																																																																																																	
Air Blank	0.000	09:25																																																																																																																																																	
Control Test	0.076	09:25																																																																																																																																																	
Air Blank	0.000	09:26																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0763																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7564																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:30																																																																																																																																																	
Control Test	0.193	09:31																																																																																																																																																	
Air Blank	0.000	09:32																																																																																																																																																	
Control Test	0.194	09:32																																																																																																																																																	
Air Blank	0.000	09:33																																																																																																																																																	
Control Test	0.195	09:34																																																																																																																																																	
Air Blank	0.000	09:34																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.1940																																																																																																																																																		
Std Dev	0.0910																																																																																																																																																		
Rel Std Dev(%)	0.5155																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:35																																																																																																																																																	
Control Test	0.077	09:35																																																																																																																																																	
Air Blank	0.000	09:36																																																																																																																																																	
Control Test	0.077	09:36																																																																																																																																																	
Air Blank	0.000	09:37																																																																																																																																																	
Control Test	0.077	09:37																																																																																																																																																	
Air Blank	0.000	09:38																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0770																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-007083, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-007083</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>CHARLOTTE COUNTY SO</u>	0.050 g/ 210 L	0.005
Calibration Date:	<u>08/02/2021</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>11:21</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Thermometer temperatures are checked with NIST traceable Eutechnics 4400 digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the uses of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

08/02/2021

Date


MICHAEL D HAUGHEY,
Department Inspector

FDLE/ATP Form 69 January 2021
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

CHARLOTTE COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000
08/12/2021

Auto Calibration

Max Power Res Value = 75

Auto Range Res Value = 72

CHANCELL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 1.4830 (-0.0080)
Sample #2 = 1.4540 (-0.0130)
Sample #3 = 1.4760 (-0.0070)
Sample #4 = 1.4810 (-0.0140)
Avg % Abs = 1.4703 (-0.0113)
STD DEV = 0.0144 (-0.0038)
REL STD DEV = 0.977 (33.405)

Sol Value = 0.100 g/210L ***
Fit value = 0.4762 mg/l ***
Samples Taken = 4, Discarded = 1
Sum lo = 12780, Sum lo = 12940

CHANCELL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1520 (-0.0030)
Sample #2 = 0.1510 (-0.0030)
Sample #3 = 0.1290 (-0.0370)
Sample #4 = 0.0860 (-0.0540)
Avg % Abs = 0.1220 (-0.0293)
STD DEV = 0.0331 (-0.0293)
REL STD DEV = 27.099 (59.761)

CHANCELL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.1060 (-0.0020)
Sample #2 = 0.1120 (-0.0000)
Sample #3 = 0.1200 (-0.0180)
Sample #4 = 0.0910 (-0.0270)
Avg % Abs = 0.1077 (-0.0150)
STD DEV = 0.0150 (-0.0137)
REL STD DEV = 13.911 (51.652)

Sol Value = 0.040 g/210L ***
Fit value = 0.1905 mg/l ***
Samples Taken = 4, Discarded = 1
Sum lo = 12790, Sum lo = 12945

CHANCELL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 0.8990 (-0.0080)
Sample #2 = 0.8440 (-0.0390)
Sample #3 = 0.8950 (-0.0270)
Sample #4 = 0.9010 (-0.0370)
Avg % Abs = 0.8800 (-0.0343)
STD DEV = 0.0313 (-0.0064)
REL STD DEV = 3.559 (18.726)

CHANCELL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 6.6230 (-0.0140)
Sample #2 = 6.6110 (-0.0080)
Sample #3 = 6.6260 (-0.0090)
Sample #4 = 6.5920 (-0.0250)
Avg % Abs = 6.6097 (-0.0140)
STD DEV = 0.0170 (-0.0095)
REL STD DEV = 0.258 (68.139)

Sol Value = 0.300 g/210L ***
Fit value = 1.4286 mg/l ***
Samples Taken = 4, Discarded = 1
Sum lo = 12767, Sum lo = 12932

CHANCELL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 5.5220 (-0.0060)
Sample #2 = 5.5020 (-0.0270)
Sample #3 = 5.5080 (-0.0370)
Sample #4 = 5.5840 (-0.0410)
Avg % Abs = 5.5147 (-0.0350)
STD DEV = 0.0031 (-0.0072)
REL STD DEV = 0.055 (20.603)

CHANCELL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 5.5670 (-0.0100)
Sample #2 = 9.5500 (-0.0150)
Sample #3 = 9.5380 (-0.0270)
Sample #4 = 9.5270 (-0.0380)
Avg % Abs = 9.5383 (-0.0267)
STD DEV = 0.0115 (-0.0115)
REL STD DEV = 0.121 (43.139)

Sol Value = 0.200 g/210L ***
Fit value = 0.5524 mg/l ***
Samples Taken = 4, Discarded = 1
Sum lo = 12773, Sum lo = 12935

CHANCELL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3.8130 (-0.0140)
Sample #2 = 3.8220 (-0.0060)
Sample #3 = 3.8110 (-0.0170)
Sample #4 = 3.7500 (-0.0620)
Avg % Abs = 3.7943 (-0.0283)
STD DEV = 0.0388 (-0.0297)
REL STD DEV = 1.022 (104.719)

AUTO CAL DATA >>>>
CHANCELL 1 >>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0
Std Dev = 0.03 Rel Std Dev = 27.10
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 0.880
Std Dev = 0.03 Rel Std Dev = 3.56
Sol Val = 0.4762 mg/l or 0.100 g/210L
% Abs = 1.975
Std Dev = 0.04 Rel Std Dev = 1.80
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 3.794
Std Dev = 0.04 Rel Std Dev = 1.02
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 5.505
Std Dev = 0.00 Rel Std Dev = 0.05
Zero Order Coef = -297.49
First Order Coef = 2485.72
Second Order Coef = 29.33
Standard Deviation = 25.33891

CHANCELL 2 >>>>
Sol Val = 0.0000 mg/l or 0.000 g/210L
% Abs = 0.108
Std Dev = 0.01 Rel Std Dev = 13.91
Sol Val = 0.1905 mg/l or 0.040 g/210L
% Abs = 3.432
Std Dev = 0.01 Rel Std Dev = 0.36
Sol Val = 0.9524 mg/l or 0.200 g/210L
% Abs = 6.610
Std Dev = 0.02 Rel Std Dev = 0.26
Sol Val = 1.4286 mg/l or 0.300 g/210L
% Abs = 9.538
Std Dev = 0.01 Rel Std Dev = 0.12
Zero Order Coef = -137.59
First Order Coef = 1366.79
Second Order Coef = 15.10
Standard Deviation = 23.699081

CHANCELL 1 >>>>
Sol Value = 0.080 g/210L ***
Fit value = 0.3810 mg/l ***
Samples Taken = 4, Discarded = 1
Sum lo = 12795, Sum lo = 12940

CHANCELL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3301.00
Sample #2 = 3310.00
Sample #3 = 3295.00
Sample #4 = 3263.00
Average Result = 3289.3333
STD DEV = 24.0069
REL STD DEV = 0.730

Dry Gas H2O Adjust Results >>>>
Barometric Pressure = 1017
3 um H2O Adjust (mg/l x 10,000) = 1027
9 um H2O Adjust (mg/l x 10,000) = 520
AUTO CAL PASS

Solution Stats Quadratic Fit Chan 2
Act Fit Residual
g/210L g/210L g/210L
0.000 0.000 -0.0002
0.040 0.040 0.0000
0.100 0.099 0.0006
0.200 0.201 -0.0007
0.300 0.300 0.0003

Sol Value = 0.080 g/210L ***
Fit value = 0.3810 mg/l ***
Samples Taken = 4, Discarded = 1
Sum lo = 12795, Sum lo = 12940

CHANCELL 1 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 2795.00
Sample #2 = 2847.00
Sample #3 = 2758.00
Sample #4 = 2741.00
Average Result = 2782.0000
STD DEV = 56.3298
REL STD DEV = 2.046

CHANCELL 2 >>>>
Sample % Abs (% Abs Ref)
Sample #1 = 3301.00
Sample #2 = 3310.00
Sample #3 = 3295.00
Sample #4 = 3263.00
Average Result = 3289.3333
STD DEV = 24.0069
REL STD DEV = 0.730

Dry Gas H2O Adjust Results >>>>
Barometric Pressure = 1017
3 um H2O Adjust (mg/l x 10,000) = 1027
9 um H2O Adjust (mg/l x 10,000) = 520
AUTO CAL PASS

Optical Calibration

SN: 80-007083
Agency: Charlotte County SO
Date: 08/02/21
By: MUK

TYPE OF TEST	SERIAL NUMBER	AGENCY	DATE	PERFORMED BY
Stabilities - PDSI	80-007683	Charlotte County SO	08/02/21	MMX

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L
0.047 to 0.053 ☒	0.077 to 0.083 ☒	0.194 to 0.206 ☒	0.077 to 0.083 ☒
CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-007083 08/02/2021 Software: 8100.27 Test g/210L Time Air Blank 0.000 09:36 Control Test 0.049 09:36 Air Blank 0.000 09:37 Control Test 0.049 09:38 Air Blank 0.000 09:38 Control Test 0.050 09:39 Air Blank 0.000 09:39 Control Test Stats Average 0.0493 Std Dev 0.0006 Rel Std Dev(%) 1.1703	CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-007083 08/02/2021 Software: 8100.27 Test g/210L Time Air Blank 0.000 09:40 Control Test 0.080 09:41 Air Blank 0.000 09:42 Control Test 0.079 09:42 Air Blank 0.000 09:43 Control Test 0.080 09:43 Air Blank 0.000 09:44 Control Test Stats Average 0.0797 Std Dev 0.0006 Rel Std Dev(%) 0.7247	CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-007083 08/02/2021 Software: 8100.27 Test g/210L Time Air Blank 0.000 09:45 Control Test 0.197 09:46 Air Blank 0.000 09:47 Control Test 0.197 09:47 Air Blank 0.000 09:48 Control Test 0.198 09:49 Air Blank 0.000 09:49 Control Test Stats Average 0.1973 Std Dev 0.0006 Rel Std Dev(%) 0.2926	CHARLOTTE COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-007083 08/02/2021 Software: 8100.27 Test g/210L Time Air Blank 0.000 09:51 Control Test 0.079 09:51 Air Blank 0.000 09:51 Control Test 0.079 09:52 Air Blank 0.000 09:52 Control Test 0.079 09:53 Air Blank 0.000 09:53 Control Test Stats Average 0.0790 Std Dev 0.0000 Rel Std Dev(%) 0.0000
MMX Operator's Signature	MMX Operator's Signature	MMX Operator's Signature	MMX Operator's Signature

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: CHARLOTTE COUNTY SO

Time of Inspection: 08:51

Date of Inspection: 08/02/2021

Serial Number: 80-007083

Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:

Number of Simulators Used: _____

Remarks:

COMPLAINT NOT DETERMINED. AI BYPASS

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



MICHAEL D HAUGHEY

Signature and Printed Name

08/02/2021
Date

Sir or Ma'am,

Please be advised that this instrument is being sent to you due to the dry gas reading being outside of tolerance.

Regards,

Cpl. Paul Dempster #0776

Charlotte County Sheriff's Office

MMX
7/26/21

Return Material Authorization

Ship to:

☒ CMI, Inc.

☐ Enforcement Electronics

Shipment to repair facility authorized by: Michael Kern on 5/27/2020

Items Returned: Instrument ☒ Supplies ☐ Other ☐ Describe: _____

Instrument Model: Intoxilyzer 8000 Serial Number: 80-007083

Bill To Address:

Charlotte County SO

Ship to Address:

Please return to FDLE Ft. Myers

Reason for Return:

Flow Sensor Replacement

Please choose one of the following options:

☐ 1. I _____, authorize all repairs.

☐ 2. I _____, authorize repairs up to \$_____.

☐ 3. I require an estimate **BEFORE** any repairs will be authorized and/ or conducted.

Please contact: Name: Michael Kern

Phone #: 941-833-6328

Email: mkern@ccsofl.net

ATP Contact Name: _____ ATP Email: _____



INSTRUMENT PROCESSING SHEET

Agency Charlotte County Sheriff's OfficeS/N 80-007083Florida Department of
Law EnforcementDate In 05/22/2020 DI Completion Date _____☐ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	Performed By <u>RAW</u>	Quality Checks	Performed By <u>XX</u>	Flow Calibration	Performed By _____																												
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>a7</u> ☐ Flow Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28427</u> ☐ Stability Checks		Flow Column # _____ ☐ 5L/min - 17mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																													
Final Release Date		<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050			0.080			0.200			0.080 DGS	N/A			Maintenance Performed By _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____														
Simulator	Serial #	Lot #/Exp																															
0.050																																	
0.080																																	
0.200																																	
0.080 DGS	N/A																																
				Temperature Checks Performed By _____ ☐ Lab Temp °C _____ External Digital Therm. ID#: _____ ☐ 34°C +- .2 Serial #: _____ ☐ 34°C +- .2 Serial #: _____ ☐ 34°C +- .2 Serial #: _____																													
Calibration Adjustment	Performed By _____	Department Inspection Performed By _____																															
Barometric Pressure Gauge ID # _____		Barometric Pressure ID# _____																															
<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th><th>Lot Number</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>	Simulator	Serial Number	Lot Number	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A				Gauge _____ Instrument _____			
Simulator	Serial Number	Lot Number	Expiration																														
0.000		N/A	N/A																														
0.040																																	
0.100																																	
0.200																																	
0.300																																	
0.080 DGS	N/A																																
		Mouth Alcohol Solution Lot # _____																															
		Acetone Stock Solution Lot # _____																															
☐ Post Calibration Adjustment Stability Checks		<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td></td></tr><tr><td>Interferent</td><td></td></tr><tr><td>0.050</td><td></td></tr><tr><td>0.080</td><td></td></tr><tr><td>0.200</td><td></td></tr></tbody></table>		Simulator	Serial Number	0.000		Interferent		0.050		0.080		0.200																			
Simulator	Serial Number																																
0.000																																	
Interferent																																	
0.050																																	
0.080																																	
0.200																																	
<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th><th>Lot Number</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>	Simulator	Serial Number	Lot Number	Expiration	0.050				0.080				0.200				0.080 DGS	N/A				Attachments											
Simulator	Serial Number	Lot Number	Expiration																														
0.050																																	
0.080																																	
0.200																																	
0.080 DGS	N/A																																
		☐ Form 41		☐ Post-Stability Checks																													
		☐ Stability Checks		☐ Flow Calibration																													
		☐ Calibration Certificate		☐ Form 40																													
		☐ Calibration Adjustment		☐ Other _____																													
Notes/Suggested Service: _____ _____ _____ _____ _____ _____		☐ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use																															
		Tech Review / Date _____		Admin Review / Date _____																													

