



INSTRUMENT PROCESSING SHEET

Agency Clewiston Police Department


S/N 80-000955

Florida Department of
Law Enforcement

Date In 08/01/2022 DI Completion Date 08/02/2022

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By <u>DERR</u>	Quality Checks	By <u>DERR</u>	Date <u>8/2/2022</u>	Flow Calibration	By _____	Date _____																																																												
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____ _____ _____ _____ _____		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>148</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP101</u> 32 mm <u>0.156</u> (.139 - .169) 36 mm <u>0.164</u> (.156 - .190) 53 mm <u>0.238</u> (.228 - .278) 103 mm <u>0.511</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28199</u> ☒ Stability Checks			Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																																														
		<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP5092</td> <td>202201C 01/11/2024</td> </tr> <tr> <td>0.080</td> <td>MP5093</td> <td>202201D 01/18/2024</td> </tr> <tr> <td>0.200</td> <td>MP5094</td> <td>202201E 01/18/2024</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>00521080A2 02/05/2023</td> </tr> </tbody> </table>			Simulator	Serial #	Lot #/Exp	0.050	MP5092	202201C 01/11/2024	0.080	MP5093	202201D 01/18/2024	0.200	MP5094	202201E 01/18/2024	0.080 DGS	N/A	00521080A2 02/05/2023	Maintenance By _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____																																															
Simulator	Serial #	Lot #/Exp																																																																	
0.050	MP5092	202201C 01/11/2024																																																																	
0.080	MP5093	202201D 01/18/2024																																																																	
0.200	MP5094	202201E 01/18/2024																																																																	
0.080 DGS	N/A	00521080A2 02/05/2023																																																																	
Calibration Adjustment By _____				Department Inspection By <u>DERR</u>																																																															
Barometric Pressure Gauge _____ ID # _____ <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> ☐ Post Calibration Adjustment Stability Checks <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table>				Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			Barometric Pressure ID# <u>26932</u> Gauge <u>1019</u> Instrument <u>1019</u> Mouth Alcohol Solution Lot # <u>2021-D</u> Acetone Stock Solution Lot # <u>2021-C</u> <table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>MP5095</td> </tr> <tr> <td>Interferent</td> <td>MP5097</td> </tr> <tr> <td>0.050</td> <td>MP5092</td> </tr> <tr> <td>0.080</td> <td>MP5093</td> </tr> <tr> <td>0.200</td> <td>MP5094</td> </tr> </tbody> </table>				Simulator	Serial Number	0.000	MP5095	Interferent	MP5097	0.050	MP5092	0.080	MP5093	0.200	MP5094
Simulator	Serial #	Lot #	Expiration																																																																
0.000		N/A	N/A																																																																
0.040																																																																			
0.100																																																																			
0.200																																																																			
0.300																																																																			
0.080 DGS	N/A																																																																		
Simulator	Serial #	Lot #	Expiration																																																																
0.050																																																																			
0.080																																																																			
0.200																																																																			
0.080 DGS	N/A																																																																		
Simulator	Serial Number																																																																		
0.000	MP5095																																																																		
Interferent	MP5097																																																																		
0.050	MP5092																																																																		
0.080	MP5093																																																																		
0.200	MP5094																																																																		
Notes/Suggested Service: <u>Bypassed AI to access the instrument from dissabled mode. DERR</u> <u>Tech review noted wrong serial number in the stabilities backdrop. DERR 8/2/22</u> _____ _____ _____				Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Calibration Adjustment ☐ Post-Stability Checks ☐ Flow Calibration ☒ Form 40 ☐ Other _____ ☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use Taylor Gutschow Digitally signed by Taylor Gutschow Date: 2022.08.02 11:32:16 -0400 04 13:34:30 Tech Review / Date Admin Review Date																																																															



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000955, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000955</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>CLEWISTON P.D.</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>08/02/2022</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>10:13</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

FDLE/ATP Form 69 December 2021
Issuing Authority: Alcohol Testing Program

08/02/2022

Date

DAVID E REYES-RIVERA,

Department Inspector

Service • Integrity • Respect • Quality

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: CLEWISTON P.D.
Time of Inspection: 10:13

Date of Inspection: 08/02/2022

Serial Number: 80-000955
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202201C Exp: 01/11/2024	0.08g/210L Test (g/210L) Lot#:202201D Exp: 01/18/2024	0.20g/210L Test (g/210L) Lot#:202201E Exp: 01/18/2024	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:00521080A2 Exp: 02/05/2023
0.000	0.049	0.079	0.201	0.081
0.000	0.049	0.078	0.201	0.081
0.000	0.049	0.078	0.200	0.081
0.000	0.049	0.078	0.200	0.081
0.000	0.050	0.078	0.200	0.081
0.000	0.049	0.079	0.200	0.081
0.000	0.050	0.079	0.200	0.081
0.000	0.050	0.078	0.201	0.081
0.000	0.050	0.078	0.200	0.081
0.000	0.050	0.078	0.201	0.081


Standard Deviations	0.0005	0.0004	0.0005	0.0000
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.


DAVID E REYES-RIVERA
 Signature and Printed Name

08/02/2022
 Date

Type of Test	Serial Number	Agency	Date	Performed By
Stabilities	80-000995-55	Clewiston Police Department	8/2/2022	DERR <i>[Signature]</i>

DERR

Digitally signed by DERR
Date: 2022.08.02 11:59:54
+0100

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L																																																																																																																																																
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083																																																																																																																																																
CLEWISTON P.D. Intoxilyzer - Alconol Analyzer Model: 8000 SN 80-000955 08/02/2022 Software: 8100.27 <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>07:53</td></tr><tr><td>Control Test</td><td>0.049</td><td>07:54</td></tr><tr><td>Air Blank</td><td>0.000</td><td>07:55</td></tr><tr><td>Control Test</td><td>0.049</td><td>07:55</td></tr><tr><td>Air Blank</td><td>0.000</td><td>07:56</td></tr><tr><td>Control Test</td><td>0.049</td><td>07:57</td></tr><tr><td>Air Blank</td><td>0.000</td><td>07:57</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0490</td><td></td></tr><tr><td>Std Dev</td><td>0.0000</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr></table>	Test	g/210L	Time	Air Blank	0.000	07:53	Control Test	0.049	07:54	Air Blank	0.000	07:55	Control Test	0.049	07:55	Air Blank	0.000	07:56	Control Test	0.049	07:57	Air Blank	0.000	07:57	Control Test Stats			Average	0.0490		Std Dev	0.0000		Rel Std Dev(%)	0.0000		CLEWISTON P.D. Intoxilyzer - Alconol Analyzer Model: 8000 SN 80-000955 08/02/2022 Software: 8100.27 <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>07:59</td></tr><tr><td>Control Test</td><td>0.079</td><td>07:59</td></tr><tr><td>Air Blank</td><td>0.000</td><td>08:00</td></tr><tr><td>Control Test</td><td>0.078</td><td>08:00</td></tr><tr><td>Air Blank</td><td>0.000</td><td>08:01</td></tr><tr><td>Control Test</td><td>0.078</td><td>08:02</td></tr><tr><td>Air Blank</td><td>0.000</td><td>08:02</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0783</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7370</td><td></td></tr></table>	Test	g/210L	Time	Air Blank	0.000	07:59	Control Test	0.079	07:59	Air Blank	0.000	08:00	Control Test	0.078	08:00	Air Blank	0.000	08:01	Control Test	0.078	08:02	Air Blank	0.000	08:02	Control Test Stats			Average	0.0783		Std Dev	0.0006		Rel Std Dev(%)	0.7370		CLEWISTON P.D. Intoxilyzer - Alconol Analyzer Model: 8000 SN 80-000955 08/02/2022 Software: 8100.27 <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>08:04</td></tr><tr><td>Control Test</td><td>0.201</td><td>08:05</td></tr><tr><td>Air Blank</td><td>0.000</td><td>08:05</td></tr><tr><td>Control Test</td><td>0.199</td><td>08:06</td></tr><tr><td>Air Blank</td><td>0.000</td><td>08:07</td></tr><tr><td>Control Test</td><td>0.200</td><td>08:07</td></tr><tr><td>Air Blank</td><td>0.000</td><td>08:08</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.2000</td><td></td></tr><tr><td>Std Dev</td><td>0.0010</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.5000</td><td></td></tr></table>	Test	g/210L	Time	Air Blank	0.000	08:04	Control Test	0.201	08:05	Air Blank	0.000	08:05	Control Test	0.199	08:06	Air Blank	0.000	08:07	Control Test	0.200	08:07	Air Blank	0.000	08:08	Control Test Stats			Average	0.2000		Std Dev	0.0010		Rel Std Dev(%)	0.5000		CLEWISTON P.D. Intoxilyzer - Alconol Analyzer Model: 8000 SN 80-000955 08/02/2022 Software: 8100.27 <table><tr><th>Test</th><th>g/210L</th><th>Time</th></tr><tr><td>Air Blank</td><td>0.000</td><td>08:09</td></tr><tr><td>Control Test</td><td>0.081</td><td>08:09</td></tr><tr><td>Air Blank</td><td>0.000</td><td>08:10</td></tr><tr><td>Control Test</td><td>0.080</td><td>08:10</td></tr><tr><td>Air Blank</td><td>0.000</td><td>08:11</td></tr><tr><td>Control Test</td><td>0.080</td><td>08:11</td></tr><tr><td>Air Blank</td><td>0.000</td><td>08:12</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0803</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7187</td><td></td></tr></table>	Test	g/210L	Time	Air Blank	0.000	08:09	Control Test	0.081	08:09	Air Blank	0.000	08:10	Control Test	0.080	08:10	Air Blank	0.000	08:11	Control Test	0.080	08:11	Air Blank	0.000	08:12	Control Test Stats			Average	0.0803		Std Dev	0.0006		Rel Std Dev(%)	0.7187	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	07:53																																																																																																																																																	
Control Test	0.049	07:54																																																																																																																																																	
Air Blank	0.000	07:55																																																																																																																																																	
Control Test	0.049	07:55																																																																																																																																																	
Air Blank	0.000	07:56																																																																																																																																																	
Control Test	0.049	07:57																																																																																																																																																	
Air Blank	0.000	07:57																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0490																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	07:59																																																																																																																																																	
Control Test	0.079	07:59																																																																																																																																																	
Air Blank	0.000	08:00																																																																																																																																																	
Control Test	0.078	08:00																																																																																																																																																	
Air Blank	0.000	08:01																																																																																																																																																	
Control Test	0.078	08:02																																																																																																																																																	
Air Blank	0.000	08:02																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0783																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7370																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	08:04																																																																																																																																																	
Control Test	0.201	08:05																																																																																																																																																	
Air Blank	0.000	08:05																																																																																																																																																	
Control Test	0.199	08:06																																																																																																																																																	
Air Blank	0.000	08:07																																																																																																																																																	
Control Test	0.200	08:07																																																																																																																																																	
Air Blank	0.000	08:08																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.2000																																																																																																																																																		
Std Dev	0.0010																																																																																																																																																		
Rel Std Dev(%)	0.5000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	08:09																																																																																																																																																	
Control Test	0.081	08:09																																																																																																																																																	
Air Blank	0.000	08:10																																																																																																																																																	
Control Test	0.080	08:10																																																																																																																																																	
Air Blank	0.000	08:11																																																																																																																																																	
Control Test	0.080	08:11																																																																																																																																																	
Air Blank	0.000	08:12																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0803																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.7187																																																																																																																																																		
 Operator's Signature	 Operator's Signature	 Operator's Signature	 Operator's Signature																																																																																																																																																

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: CLEWISTON P.D.
Time of Inspection: 07:05

Date of Inspection: 08/02/2022

Serial Number: 80-000955
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:


Number of Simulators Used: _____

Remarks:

COMPLIANCE NOT DETERMINED, AI NOT CONDUCTED.

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

 _____
Signature and Printed Name

DAVID E REYES-RIVERA

08/02/2022
Date