

INSTRUMENT PROCESSING SHEET

Agency Dade City PD

S/N 80-001073

Florida Department of
Law Enforcement

Date In 4/28/2022 DI Completion Date 5/3/2022

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By TDG	Quality Checks	By TDG	Date 5/3/2022	Flow Calibration	By	Date																																																											
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☒ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: _____ _____ _____ _____ _____ _____ _____ _____ _____ _____		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>236</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP104</u> 32 mm <u>0.152</u> (.139 - .169) 36 mm <u>0.167</u> (.156 - .190) 53 mm <u>0.246</u> (.228 - .278) 103 mm <u>0.503</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>68639</u> ☒ Stability Checks			Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																																																													
		<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #/Exp</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td rowspan="2">MP5092</td> <td>202010A 10/05/2022</td> </tr> <tr> <td>0.080</td> <td>202010B 10/05/2022</td> </tr> <tr> <td>0.200</td> <td rowspan="2">MP5094</td> <td>202010D 10/06/2022</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A AG115904 06/08/2023</td> </tr> </tbody> </table>			Simulator	Serial #	Lot #/Exp	0.050	MP5092	202010A 10/05/2022	0.080	202010B 10/05/2022	0.200	MP5094	202010D 10/06/2022	0.080 DGS	N/A AG115904 06/08/2023	Maintenance By _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____ _____ _____ _____ _____																																																
Simulator	Serial #	Lot #/Exp																																																																
0.050	MP5092	202010A 10/05/2022																																																																
0.080		202010B 10/05/2022																																																																
0.200	MP5094	202010D 10/06/2022																																																																
0.080 DGS		N/A AG115904 06/08/2023																																																																
Calibration Adjustment By _____ Barometric Pressure Gauge _____ ID # _____ <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> ☐ Post Calibration Adjustment Stability Checks <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> </tbody> </table> Notes/Suggested Service: _____ _____ _____ _____ _____ _____ _____ _____ _____		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			Department Inspection By TDG _____ Barometric Pressure ID# <u>28199</u> Gauge <u>1017</u> Instrument <u>1011</u> Mouth Alcohol Solution Lot # <u>2021-B</u> Acetone Stock Solution Lot # <u>2021-C</u> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial Number</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>SD3963</td> </tr> <tr> <td>Interferent</td> <td>SD3966</td> </tr> <tr> <td>0.050</td> <td>MP5092</td> </tr> <tr> <td>0.080</td> <td>MP5093</td> </tr> <tr> <td>0.200</td> <td>MP5094</td> </tr> </tbody> </table> Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Calibration Adjustment ☐ Post-Stability Checks ☐ Flow Calibration ☐ Form 40 ☐ Other _____ ☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use 2022.05 05 13:34:4 2024.06 Israel Soto Tech Review / Date Admin Review / Date					Simulator	Serial Number	0.000	SD3963	Interferent	SD3966	0.050	MP5092	0.080	MP5093	0.200	MP5094
Simulator	Serial #	Lot #	Expiration																																																															
0.000		N/A	N/A																																																															
0.040																																																																		
0.100																																																																		
0.200																																																																		
0.300																																																																		
0.080 DGS	N/A																																																																	
Simulator	Serial #	Lot #	Expiration																																																															
0.050																																																																		
0.080																																																																		
0.200																																																																		
0.080 DGS	N/A																																																																	
Simulator	Serial Number																																																																	
0.000	SD3963																																																																	
Interferent	SD3966																																																																	
0.050	MP5092																																																																	
0.080	MP5093																																																																	
0.200	MP5094																																																																	

Type of Test	Serial Number	Agency	Date	Performed By
Stabilities	80-001073	Dade City PD	05/03/2022	TDG <i>ML</i>

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L																																																																																																																																				
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083																																																																																																																																				
DADE CITY PD Intoxilyzer - Alcohol Analyzer Model: 8000 05/03/2022 Software: 8100.27 Test: 9/210L Time <table><tr><td>Air Blank</td><td>0.000</td><td>12:51</td></tr><tr><td>Control Test</td><td>0.049</td><td>12:52</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:53</td></tr><tr><td>Control Test</td><td>0.049</td><td>12:53</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:54</td></tr><tr><td>Control Test</td><td>0.049</td><td>12:55</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:55</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0490</td><td></td></tr><tr><td>Std Dev</td><td>0.0000</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr></table> Operator's Signature ML	Air Blank	0.000	12:51	Control Test	0.049	12:52	Air Blank	0.000	12:53	Control Test	0.049	12:53	Air Blank	0.000	12:54	Control Test	0.049	12:55	Air Blank	0.000	12:55	Control Test Stats			Average	0.0490		Std Dev	0.0000		Rel Std Dev(%)	0.0000		DADE CITY PD Intoxilyzer - Alcohol Analyzer Model: 8000 05/03/2022 Software: 8100.27 Test: 9/210L Time <table><tr><td>Air Blank</td><td>0.000</td><td>12:41</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:42</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:43</td></tr><tr><td>Control Test</td><td>0.078</td><td>12:43</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:44</td></tr><tr><td>Control Test</td><td>0.078</td><td>12:45</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:45</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0783</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7370</td><td></td></tr></table> Operator's Signature ML	Air Blank	0.000	12:41	Control Test	0.079	12:42	Air Blank	0.000	12:43	Control Test	0.078	12:43	Air Blank	0.000	12:44	Control Test	0.078	12:45	Air Blank	0.000	12:45	Control Test Stats			Average	0.0783		Std Dev	0.0006		Rel Std Dev(%)	0.7370		DADE CITY PD Intoxilyzer - Alcohol Analyzer Model: 8000 05/03/2022 Software: 8100.27 Test: 9/210L Time <table><tr><td>Air Blank</td><td>0.000</td><td>12:34</td></tr><tr><td>Control Test</td><td>0.197</td><td>12:34</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:35</td></tr><tr><td>Control Test</td><td>0.198</td><td>12:36</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:36</td></tr><tr><td>Control Test</td><td>0.197</td><td>12:37</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:37</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.1973</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.2926</td><td></td></tr></table> Operator's Signature ML	Air Blank	0.000	12:34	Control Test	0.197	12:34	Air Blank	0.000	12:35	Control Test	0.198	12:36	Air Blank	0.000	12:36	Control Test	0.197	12:37	Air Blank	0.000	12:37	Control Test Stats			Average	0.1973		Std Dev	0.0006		Rel Std Dev(%)	0.2926		DADE CITY PD Intoxilyzer - Alcohol Analyzer Model: 8000 05/03/2022 Software: 8100.27 Test: 9/210L Time <table><tr><td>Air Blank</td><td>0.000</td><td>12:57</td></tr><tr><td>Control Test</td><td>0.060</td><td>12:57</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:57</td></tr><tr><td>Control Test</td><td>0.079</td><td>12:58</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:58</td></tr><tr><td>Control Test</td><td>0.080</td><td>12:59</td></tr><tr><td>Air Blank</td><td>0.000</td><td>12:59</td></tr><tr><td>Control Test Stats</td><td></td><td></td></tr><tr><td>Average</td><td>0.0797</td><td></td></tr><tr><td>Std Dev</td><td>0.0006</td><td></td></tr><tr><td>Rel Std Dev(%)</td><td>0.7247</td><td></td></tr></table> Operator's Signature ML	Air Blank	0.000	12:57	Control Test	0.060	12:57	Air Blank	0.000	12:57	Control Test	0.079	12:58	Air Blank	0.000	12:58	Control Test	0.080	12:59	Air Blank	0.000	12:59	Control Test Stats			Average	0.0797		Std Dev	0.0006		Rel Std Dev(%)	0.7247	
Air Blank	0.000	12:51																																																																																																																																					
Control Test	0.049	12:52																																																																																																																																					
Air Blank	0.000	12:53																																																																																																																																					
Control Test	0.049	12:53																																																																																																																																					
Air Blank	0.000	12:54																																																																																																																																					
Control Test	0.049	12:55																																																																																																																																					
Air Blank	0.000	12:55																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.0490																																																																																																																																						
Std Dev	0.0000																																																																																																																																						
Rel Std Dev(%)	0.0000																																																																																																																																						
Air Blank	0.000	12:41																																																																																																																																					
Control Test	0.079	12:42																																																																																																																																					
Air Blank	0.000	12:43																																																																																																																																					
Control Test	0.078	12:43																																																																																																																																					
Air Blank	0.000	12:44																																																																																																																																					
Control Test	0.078	12:45																																																																																																																																					
Air Blank	0.000	12:45																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.0783																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.7370																																																																																																																																						
Air Blank	0.000	12:34																																																																																																																																					
Control Test	0.197	12:34																																																																																																																																					
Air Blank	0.000	12:35																																																																																																																																					
Control Test	0.198	12:36																																																																																																																																					
Air Blank	0.000	12:36																																																																																																																																					
Control Test	0.197	12:37																																																																																																																																					
Air Blank	0.000	12:37																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.1973																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.2926																																																																																																																																						
Air Blank	0.000	12:57																																																																																																																																					
Control Test	0.060	12:57																																																																																																																																					
Air Blank	0.000	12:57																																																																																																																																					
Control Test	0.079	12:58																																																																																																																																					
Air Blank	0.000	12:58																																																																																																																																					
Control Test	0.080	12:59																																																																																																																																					
Air Blank	0.000	12:59																																																																																																																																					
Control Test Stats																																																																																																																																							
Average	0.0797																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.7247																																																																																																																																						

Comments:

Florida Department of Law Enforcement

Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: DADE CITY PD
Time of Inspection: 15:32

Date of Inspection: 05/03/2022

Serial Number: 80-001073
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202010A Exp: 10/05/2022	0.08g/210L Test (g/210L) Lot#:202010B Exp: 10/05/2022	0.20g/210L Test (g/210L) Lot#:202010D Exp: 10/06/2022	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG115904 Exp: 06/08/2023
0.000	0.049	0.078	0.181 / 0.196	0.079
0.000	0.048	0.078	INT / 0.196	0.080
0.000	0.049	0.078	/ 0.196	0.080
0.000	0.049	0.078	/ 0.196	0.080
0.000	0.048	0.078	/ 0.195	0.080
0.000	0.049	0.078	/ 0.195	0.080
0.000	0.050	0.078	/ 0.196	0.080
0.000	0.050	0.078	/ 0.195	0.080
0.000	0.049	0.078	/ 0.195	0.080
0.000	0.050	0.078	/ 0.195	0.080

Standard Deviations	0.0007	0.0000	/ 0.0005	0.0003
---------------------	--------	--------	----------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

20: Interferent Detect..

Unknown interferent detected during 0.20 Test. Allowed room to air out, then repeated the test. No interferent detected during repeated test.

TLB
05/03/2022

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Taylor D Gutschow

TAYLOR D GUTSCHOW

Signature and Printed Name

05/03/2022
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001073, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001073</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>DADE CITY PD</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>05/03/2022</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>15:32</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.
*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.
This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.


Date 05/03/2022
Taylor D GUTSCHOW,
Department Inspector

FDLE/ATP Form 69 December 2021
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality