

INSTRUMENT PROCESSING SHEET

S/N 80-006240

aw Enforce	By TC)G	Quality Che	cks By TDG	Date 01/10/2023	Flow Calib	ration By Date
Annual	-1-			ube Screen			nn #
☐ Registrati	on		Replace External O-Rings			☐ 5L/min – 17mm	
	om CMI / EE		■ Instrument Set Up Verified			☐ 15L/min – 53mm	
			R-Value		VXX 575040	□ 30L/min – 103mm	
Visual Inspection:				rification (L/s)		□ R-Value	
Case Handle			Flow Column # ATP104			☐ Post Calibration Verification (L/s)	
Keyboard Dry Gas Shelf					(.139169)	Flow Colun	nn #
Feet Breath Tube			36 mm 0).167	(.156190)	32 mm	(.139169)
Ports Screws Tight			53 mm 0.234				(.156190)
	ment/ Accessories:		103 mm	0.507	(.447547)	53 mm	(.228278)
Power co	rd 🔲 Printer Cab	le	■ Baromet	ric Pressure Ch	eck		(.447547)
☐ Static Bag	g 🔲 12V DC Cab	le	Gauge ID #	68639			
Notes:			■ Stability				
			Simulator		Lot #/Exp	Maintenar	ce By DERR
			0.050		2022010		Replacement
	V.		0.050	MP6286	202201C	The state of the s	Regulator Replacement
			0.000		01/11/2024		ube Replacement
			0.080	MP6287	202201D		Replaced int. printer paper
				7,000,000	01/18/2024		
			0.200	MP5094	202201E		
el				WII 000 T	01/18/2024		
			0.080 DGS	N/A	00521080A2		
				1'	02/05/2023		
Calibration	Adjustment	3.3.3	Same 5	Ву	Department Inspec		By TDG
Barometric	Pressure Gauge		ID #	- 7	Barometric Pressur		
Simulator	Serial #	Lot#		Expiration	Gauge <u>1020</u>		strument 1020
0.000			N/A	N/A	Mouth Alcohol Solu		
0.040				.3	Acetone Stock Solu		
0.100					Simulator .	emer definit	Serial Number
0.200					0.000		MP6284
1 0.200		1					
0.200					Interferent		MP6285
0.300	N/A)	Interferent 0.050		MP6285 MP6286
0.300 0.080 DGS	N/A			,)	Interferent		MP6285
0.300 0.080 DGS	bration Adjustment	Stabilit	y Checks		Interferent 0.050 0.080 0.200		MP6285 MP6286 MP6287
0.300 0.080 DGS Post Calil Simulator		Stabilit	y Checks	Expiration	Interferent 0.050 0.080 0.200 Attachments		MP6285 MP6286 MP6287 MP5094
0.300 0.080 DGS Post Calil Simulator 0.050	bration Adjustment	Stabilit Lot#	y Checks		Interferent 0.050 0.080 0.200 Attachments ■ Form 41		MP6285 MP6286 MP6287 MP5094
0.300 0.080 DGS Post Calil Simulator	bration Adjustment	Stabilit Lot#	y Checks		Interferent 0.050 0.080 0.200 Attachments Form 41 Stability Checks		MP6285 MP6286 MP6287 MP5094 Post-Stability Checks Flow Calibration
0.300 0.080 DGS Post Calil Simulator 0.050	bration Adjustment	Stabilit Lot #	y Checks		Interferent 0.050 0.080 0.200 Attachments Form 41 Stability Checks Calibration Cer	tificate	MP6285 MP6286 MP6287 MP5094 Post-Stability Checks Flow Calibration Form 40
0.300 0.080 DGS Post Calil Simulator 0.050 0.080	bration Adjustment	Stabilit Lot#	y Checks		Interferent 0.050 0.080 0.200 Attachments Form 41 Stability Checks	tificate	MP6285 MP6286 MP6287 MP5094 Post-Stability Checks Flow Calibration
0.300 0.080 DGS Post Calil Simulator 0.050 0.080 0.200 0.080 DGS Notes/Sugg	bration Adjustment Serial #	Lot#	aced the ba	Expiration attery and	Interferent 0.050 0.080 0.200 Attachments Form 41 Stability Checks Calibration Cer Calibration Adj Instrument Co Instrument Do Return to/Place	tificate ustment mplies with es Not Comp e into Evide	MP6285 MP6286 MP6287 MP5094 Post-Stability Checks Flow Calibration Other Other Chapter 11D-8, FAC Dly with Chapter 11D-8, FAC
0.300 0.080 DGS Post Calil Simulator 0.050 0.080 0.200 0.080 DGS Notes/Sugg	bration Adjustment Serial #	Lot#	aced the ba	Expiration attery and	Interferent 0.050 0.080 0.200 Attachments Form 41 Stability Checks Calibration Cer Calibration Adj Instrument Co Instrument Do Return to/Plac	tificate ustment mplies with es Not Complee into Evidentiary	MP6285 MP6286 MP6287 MP5094 Post-Stability Checks Flow Calibration Other Chapter 11D-8, FAC Oly with Chapter 11D-8, FAC
0.300 0.080 DGS Post Calil Simulator 0.050 0.080 0.200 0.080 DGS Notes/Sugg	bration Adjustment Serial #	Lot#	aced the ba	Expiration attery and	Interferent 0.050 0.080 0.200 Attachments Form 41 Stability Checks Calibration Cer Calibration Adj Instrument Co Instrument Do Return to/Plac	tificate ustment mplies with es Not Complee into Evidentiary	MP6285 MP6286 MP6287 MP5094 Post-Stability Checks Flow Calibration Other Chapter 11D-8, FAC Oly with Chapter 11D-8, FAC Intiary Use Use

Type of Test	Serial Number	Agency	Date ,	Performed By
Stabilities	80-00 6240	Palm Beach OD	01/10/2023	TDG MG-

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L	
0.047 to 0.053	0.077 to 0.083	0.194 to 0.206	0.077 to 0.083 ✓ ≤0.003 of Wet	
PALM BEACH CO SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006240 01/10/2023 Software: 8100.27 Test g/210L Time	PALM BEACH CO SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006240 01/10/2023 Software: 8100.27 Test g/210L Time Air Blank 0.000 19:38 Control Test 0.079 10:39 Control Test 0.079 10:40 Air Blank 0.000 10:40 Control Test 0.079 10:40 Control Test 0.079 10:41 Control Test 0.079 10:41 Air Blank 0.000 10:41 Control Test 0.079 10:41 Control Test Stats Average 0.0790 Std Dev 0.0000 Rel Std Dev(%) 0.0000	PALM BEACH CO SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-006240 D1/10/2023 Software: 8100.27 Test g/210L Time	PALM BEACH CO SO Intoxilyzer - Alconol Analyzer Model 8000\ SN 80-006240 01/10/2023\ Software: 8100.27 Test g/210L Time Air Blank 0.000 10:02 Control Test 0.082 10:03 Air Blank 0.000 10:03 Control Test 0.081 10:04 Air Blank 0.000 10:04 Control Test 0.090 10:04 Control Test 0.090 10:04 Control Test 0.090 10:05 Control Test 0.000 10:05 Control Test Stats Average 0.0810 Std Dev 0.0010 Rel Std Dev(%) 1.2346	
Operator's Signature	M C- Operator's Signature	Operator's Signature	Operator's Signature	

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: PALM BEACH CO SO Time of Inspection: 14:25

Date of Inspection: 01/10/2023

Serial Number: 80-006240

Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check			Date and/or Time Adjusted		
(Pre-Inspection): OK	Yes				No
Minimum Sample Volume			Barometric Pressure Sensor		
Check: OK	Yes		Check: OK	Yes	
Alcohol Free Subject			Mouth Alcohol Test:		
Test: 0.000	Yes		Slope Not Met	Yes	
Interferent Detect Test:			Diagnostic Check		
Interferent Detect	Yes		(Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202201C Exp: 01/11/2024	0.08g/210L Test (g/210L) Lot#:202201D Exp: 01/18/2024	0.20g/210L Test (g/210L) Lot#:202201E Exp: 01/18/2024	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:00521080A2 Exp: 02/05/2023
0.000	0.050	0.080	0.202	0.081
0.000	0.050	0.081	0.201	0.081
0.000	0.050	0.080	0.201	0.080
0.000	0.050	0.080	0.201	0.081
0.000	0.050	0.080	0.201	0.081
0.000	0.050	0.080	0.201	0.081
0.000	0.050	0.079	0.201	0.080
0.000	0.050	0.080	0.201	0.081
0.000	0.050	0.080	0.200	0.081
0.000	0.050	0.080	0:200	0.081
			3	
Standard Deviations	0 0000	0.0004	0.0005	0.0004

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

The above instrument complies (X) does not comply () with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

TAYLOR D GUTSCHOW

Signature and Printed Name

01/10/2023 Date



Calibration Certificate

Florida Department of Law Enforcement Alcohol Testing Program 4700 Terminal Drive, Suite 1 Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-006240, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	80-006240	UNCERTAINTY* ±	
Owning Agency:	PALM BEACH CO SO	0.050 g/ 210 L	0.004
Calibration Date:	01/10/2023	0.080 g/ 210 L	0.004
Calibration Time:	14:25	0.200 g/ 210 L	0.007
		0.080 g/210 L Dry Gas Control	0.005

All results are reported in g/210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration. *Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full,

without written approval of the Florida Department of

Law Enforcement Alcohol Testing Program.

01/10/2023

Date

OR D GUTSCHOW.

Department Inspector

FDLE/ATP Form 69 December 2021 Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality

Page 1 of 1