



INSTRUMENT PROCESSING SHEET

Agency Pinellas County SOS/N 80-005290Florida Department of
Law EnforcementDate In 04/22/2024 DI Completion Date 5/16/2024☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By <u>ALL</u> Date <u>04/22/2024</u> ☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☐ Handle ☒ Keyboard ☐ Dry Gas Shelf ☒ Feet ☐ Breath Tube ☒ Ports ☐ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes: <u>missing left peg holding keyboard</u>	Quality Checks By <u>ALL</u> Date <u>04/23/2024</u> ☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>140</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP-105</u> 32 mm <u>.152</u> (.139 - .169) 36 mm <u>.171</u> (.156 - .190) 53 mm <u>.238</u> (.228 - .278) 103 mm <u>.515</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28421</u> ☒ Stability Checks <table border="1" style="width:100%"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP6291</td><td>202303K 03/29/2025</td></tr><tr><td>0.080</td><td>MP6292</td><td>202303L 03/29/2025</td></tr><tr><td>0.200</td><td>MP6293</td><td>202304C 04/05/2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG325603 09/13/2025</td></tr></tbody></table>	Simulator	Serial #	Lot #/Exp	0.050	MP6291	202303K 03/29/2025	0.080	MP6292	202303L 03/29/2025	0.200	MP6293	202304C 04/05/2025	0.080 DGS	N/A	AG325603 09/13/2025	Flow Calibration By _____ Date _____ Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547) Maintenance By _____ Date _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____																																	
Simulator	Serial #	Lot #/Exp																																																
0.050	MP6291	202303K 03/29/2025																																																
0.080	MP6292	202303L 03/29/2025																																																
0.200	MP6293	202304C 04/05/2025																																																
0.080 DGS	N/A	AG325603 09/13/2025																																																
Calibration Adjustment By <u>PN</u> Barometric Pressure Gauge <u>1008</u> ID # <u>28427</u> <table border="1" style="width:100%"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td>MP6294</td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td>MP6295</td><td>23400</td><td>10/24/2025</td></tr><tr><td>0.100</td><td>MP6296</td><td>23450</td><td>12/5/2025</td></tr><tr><td>0.200</td><td>MP6297</td><td>24080</td><td>2/13/2026</td></tr><tr><td>0.300</td><td>MP6298</td><td>23410</td><td>11/1/2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG310901</td><td>04/19/2025</td></tr></tbody></table> ☒ Post Calibration Adjustment Stability Checks <table border="1" style="width:100%"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td>MP5088</td><td>202303K</td><td>03/29/2025</td></tr><tr><td>0.080</td><td>MP5089</td><td>202303L</td><td>03/29/2025</td></tr><tr><td>0.200</td><td>MP5090</td><td>202304C</td><td>04/05/2025</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>06723080A5</td><td>04/05/2025</td></tr></tbody></table> Notes/Suggested Service: <u>ran stability check with .150g/210L DGS 04/23/2024 ALL</u> <u>stability for .080g/210L DGS repeated, values outside tolerance 04/23/2024 ALL</u>			Simulator	Serial #	Lot #	Expiration	0.000	MP6294	N/A	N/A	0.040	MP6295	23400	10/24/2025	0.100	MP6296	23450	12/5/2025	0.200	MP6297	24080	2/13/2026	0.300	MP6298	23410	11/1/2025	0.080 DGS	N/A	AG310901	04/19/2025	Simulator	Serial #	Lot #	Expiration	0.050	MP5088	202303K	03/29/2025	0.080	MP5089	202303L	03/29/2025	0.200	MP5090	202304C	04/05/2025	0.080 DGS	N/A	06723080A5	04/05/2025
Simulator	Serial #	Lot #	Expiration																																															
0.000	MP6294	N/A	N/A																																															
0.040	MP6295	23400	10/24/2025																																															
0.100	MP6296	23450	12/5/2025																																															
0.200	MP6297	24080	2/13/2026																																															
0.300	MP6298	23410	11/1/2025																																															
0.080 DGS	N/A	AG310901	04/19/2025																																															
Simulator	Serial #	Lot #	Expiration																																															
0.050	MP5088	202303K	03/29/2025																																															
0.080	MP5089	202303L	03/29/2025																																															
0.200	MP5090	202304C	04/05/2025																																															
0.080 DGS	N/A	06723080A5	04/05/2025																																															
Department Inspection By <u>PN</u> Barometric Pressure ID# <u>28421</u> Gauge <u>1009</u> Instrument <u>1007</u> Mouth Alcohol Solution Lot # <u>2024-A</u> Acetone Stock Solution Lot # <u>2023-B</u> <table border="1" style="width:100%"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5086</td></tr><tr><td>Interferent</td><td>MP5087</td></tr><tr><td>0.050</td><td>MP5088</td></tr><tr><td>0.080</td><td>MP5089</td></tr><tr><td>0.200</td><td>MP5090</td></tr></tbody></table> Attachments <table border="1" style="width:100%"><tbody><tr><td>☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☒ Calibration Adjustment</td><td>☒ Post-Stability Checks ☐ Flow Calibration ☒ Form 40 ☐ Other _____</td></tr></tbody></table> ☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use Digitally signed by Taylor Gutschow Date: 2024.05.29 12:21:30 -0400 Shayla Platt Date: 2024.06.03 13:38:35 -0400 Tech Review / Date _____ Admin Review / Date _____			Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090	☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☒ Calibration Adjustment	☒ Post-Stability Checks ☐ Flow Calibration ☒ Form 40 ☐ Other _____																																		
Simulator	Serial Number																																																	
0.000	MP5086																																																	
Interferent	MP5087																																																	
0.050	MP5088																																																	
0.080	MP5089																																																	
0.200	MP5090																																																	
☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☒ Calibration Adjustment	☒ Post-Stability Checks ☐ Flow Calibration ☒ Form 40 ☐ Other _____																																																	

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: PINELLAS COUNTY SO
Time of Inspection: 09:55

Date of Inspection: 04/23/2024

Serial Number: 80-005290
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK	Yes	
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#: Exp:	0.08g/210L Test (g/210L) Lot#: Exp:	0.20g/210L Test (g/210L) Lot#: Exp:	0.08 g/210L Dry Gas Std Test (g/210L) Lot#: Exp:


Number of Simulators Used: 0

Remarks:

BYPASS AI FOR INSPECTION COMPLIANCE NOT DETERMINED

The above instrument complies (04/23/24 AL) does not comply () with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



ALFONSO L LOWRY

Signature and Printed Name

04/23/2024
Date

stability checks 80-005290 04/23/24

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
04/23/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:43
Control Test	0.050	11:44
Air Blank	0.000	11:44
Control Test	0.050	11:45
Air Blank	0.000	11:45
Control Test	0.050	11:46
Air Blank	0.000	11:46
Control Test Stats		
Average	0.0500	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

AZ
Operator's Signature

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
04/23/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:51
Control Test	0.080	11:51
Air Blank	0.000	11:52
Control Test	0.080	11:52
Air Blank	0.000	11:53
Control Test	0.079	11:54
Air Blank	0.000	11:54
Control Test Stats		
Average	0.0797	
Std Dev	0.0006	
Rel Std Dev(%)	0.7247	

Wet
AZ
Operator's Signature

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
04/23/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:59
Control Test	0.199	12:00
Air Blank	0.000	12:00
Control Test	0.198	12:01
Air Blank	0.000	12:01
Control Test	0.198	12:02
Air Blank	0.000	12:03
Control Test Stats		
Average	0.1983	
Std Dev	0.0006	
Rel Std Dev(%)	0.2911	

AZ
Operator's Signature

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
04/23/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:33
Control Test	0.085	12:33
Air Blank	0.000	12:34
Control Test	0.084	12:34
Air Blank	0.000	12:34
Control Test	0.084	12:35
Air Blank	0.000	12:35
Control Test Stats		
Average	0.0843	
Std Dev	0.0006	
Rel Std Dev(%)	0.6846	

D65
AZ
Operator's Signature

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
04/23/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	12:37
Control Test	0.084	12:38
Air Blank	0.000	12:38
Control Test	0.083	12:38
Air Blank	0.000	12:39
Control Test	0.083	12:39
Air Blank	0.000	12:40
Control Test Stats		
Average	0.0833	
Std Dev	0.0006	
Rel Std Dev(%)	0.6928	

D65
AZ
Operator's Signature

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 00-005290
04/23/2024
Software: 8100.27


Test	g/210L	Time
Air Blank	0.000	12:06
Control Test	0.153	12:06
Air Blank	0.000	12:07
Control Test	0.153	12:07
Air Blank	0.000	12:07
Control Test	0.153	12:08
Air Blank	0.000	12:08
Control Test Stats		
Average	0.1530	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

AZ
Operator's Signature

can stability with
150g/210L OGS
04/23/24
AZ

CALIBRATION

ADJUSTMENT



PINELLAS COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-005290
 05/16/2024 09:56:26

Auto Calibration

pg 1 of 2

<<<< 3um >>>>			<<<< 9um >>>>		

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	0.0370	(-0.0060)	0.1250	(-0.0070)	
Sample #2	0.0280	(0.0350)	0.1200	(0.0050)	
Sample #3	0.0550	(0.0630)	0.1400	(0.0070)	
Sample #4	0.0660	(0.0760)	0.1350	(0.0120)	
Avg % Abs	0.0497	(0.0580)	0.1317	(0.0080)	
STD DEV	0.0196	(0.0210)	0.0104	(0.0036)	
REL STD DEV	39.369	(36.125)	7.905	(45.069)	

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	0.7650	(-0.0040)	1.5340	(-0.0100)	
Sample #2	0.7640	(0.0100)	1.5330	(-0.0110)	
Sample #3	0.7530	(0.0240)	1.5410	(-0.0110)	
Sample #4	0.7610	(0.0320)	1.5320	(-0.0020)	
Avg % Abs	0.7593	(0.0220)	1.5353	(-0.0080)	
STD DEV	0.0057	(0.0111)	0.0049	(0.0052)	
REL STD DEV	0.749	(50.616)	0.321	(64.952)	

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	1.8810	(-0.0140)	3.6440	(-0.0070)	
Sample #2	1.9090	(-0.0040)	3.6490	(0.0020)	
Sample #3	1.9000	(0.0050)	3.6520	(0.0090)	
Sample #4	1.8880	(0.0170)	3.6350	(0.0110)	
Avg % Abs	1.8990	(0.0060)	3.6453	(0.0073)	
STD DEV	0.0105	(0.0105)	0.0091	(0.0047)	
REL STD DEV	0.555	(175.594)	0.249	(64.443)	

Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	3.5850	(-0.0060)	6.8370	(0.0170)	
Sample #2	3.6390	(-0.0180)	6.8710	(-0.0020)	
Sample #3	3.5990	(-0.0090)	6.8540	(0.0070)	
Sample #4	3.6190	(0.0020)	6.8760	(0.0080)	
Avg % Abs	3.6190	(-0.0083)	6.8670	(0.0043)	
STD DEV	0.0200	(0.0100)	0.0115	(0.0055)	
REL STD DEV	0.553	(120.200)	0.168	(127.098)	

Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1					
Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)	
Sample #1	5.2850	(-0.0100)	9.9910	(0.0000)	
Sample #2	5.3290	(-0.0140)	9.9840	(0.0090)	
Sample #3	5.3160	(0.0070)	9.9970	(0.0090)	
Sample #4	5.2940	(0.0130)	9.9940	(0.0200)	
Avg % Abs	5.3130	(0.0020)	9.9917	(0.0127)	
STD DEV	0.0177	(0.0142)	0.0068	(0.0064)	
REL STD DEV	0.333	(708.872)	0.068	(50.138)	

CALIBRATION

ADJUSTMENT

(P)

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
05/16/2024 09:56:26

Auto Calibration

pg 2 of 2

<<<< 3um >>>>

Zero Order Coef -99.50
First Order Coef 2527.49
Second Order Coef 34.13

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.001	-0.0005
0.040	0.039	0.0014
0.100	0.101	-0.0013
0.200	0.199	0.0006
0.300	0.300	-0.0001

<<<< 9um >>>>

-171.43
1317.44
13.06

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0000
0.040	0.040	0.0005
0.100	0.101	-0.0009
0.200	0.199	0.0007
0.300	0.300	-0.0002

<<<< 3um >>>> <<<< 9um >>>>

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
Sample
Sample #1 3421.00 3333.00
Sample #2 3332.00 3304.00
Sample #3 3314.00 3295.00
Sample #4 3255.00 3267.00
Avg 3300.3333 3288.6667
STD DEV 40.2782 19.2959
REL STD DEV 1.220 0.587
H2O adjust (mg/l*10k) 509 521

Barometric Pressure = 1008

*****CALIBRATION SUCCESSFUL*****

POST STABILITY CHECKS

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
05/16/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:43
Control Test	0.050	11:43
Air Blank	0.000	11:44
Control Test	0.050	11:45
Air Blank	0.000	11:45
Control Test	0.050	11:46
Air Blank	0.000	11:46
Control Test Stats		
Average	0.0500	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	


Operator's Signature

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
05/16/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:47
Control Test	0.079	11:48
Air Blank	0.000	11:49
Control Test	0.079	11:49
Air Blank	0.000	11:50
Control Test	0.079	11:51
Air Blank	0.000	11:51
Control Test Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	


Operator's Signature


PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
05/16/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:52
Control Test	0.198	11:53
Air Blank	0.000	11:54
Control Test	0.197	11:54
Air Blank	0.000	11:55
Control Test	0.197	11:56
Air Blank	0.000	11:56
Control Test Stats		
Average	0.1973	
Std Dev	0.0006	
Rel Std Dev(%)	0.2926	


Operator's Signature

PINELLAS COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-005290
05/16/2024
Software: 8100.27

Test	g/210L	Time
Air Blank	0.000	11:57
Control Test	0.080	11:57
Air Blank	0.000	11:58
Control Test	0.081	11:58
Air Blank	0.000	11:59
Control Test	0.080	11:59
Air Blank	0.000	11:59
Control Test Stats		
Average	0.0803	
Std Dev	0.0006	
Rel Std Dev(%)	0.7187	


Operator's Signature

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: PINELLAS COUNTY SO
Time of Inspection: 14:20

Date of Inspection: 05/16/2024

Serial Number: 80-005290
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202303K Exp: 03/29/2025	0.08g/210L Test (g/210L) Lot#:202303L Exp: 03/29/2025	0.20g/210L Test (g/210L) Lot#:202304C Exp: 04/05/2025	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:06723080A5 Exp: 04/05/2025
0.000	0.050	0.080	0.198	0.081
0.000	0.050	0.080	0.199	0.080
0.000	0.050	0.080	0.199	0.080
0.000	0.050	0.080	0.199	0.080
0.000	0.050	0.080	0.199	0.080
0.000	0.051	0.080	0.199	0.080
0.000	0.051	0.080	0.199	0.080
0.000	0.051	0.080	0.199	0.080
0.000	0.051	0.080	0.199	0.080
0.000	0.051	0.080	0.199	0.080

Standard Deviations	0.0005	0.0000	0.0003	0.0003
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0002 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



PHIL NICODEMO
Signature and Printed Name

05/16/2024
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road.
Suite B1032
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-005290, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-005290</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>PINELLAS COUNTY SO</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>05/16/2024</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>14:20</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

05/16/2024

Date

PHIL NICODEMO,

Department Inspector

FDLE/ATP Form 69 March 2022

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality