



INSTRUMENT PROCESSING SHEET

Agency Hillsborough County Sheriff's Office

S/N 80-007483

Florida Department of
Law Enforcement

Date In 04/25/2025 DI Completion Date 05/20/2025

☐ Ship ☒ P/U ☐ H/D ☐ CMI ☐ EE

Intake	By TDG	Date 05/19/2025	Quality Checks	By TDG	Date 05/19/2025	Flow Calibration	By	Date																																																																																															
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes: Dropped off with box.			☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value 124 ☒ Flow Verification (L/s) Flow Column # ATP104 32 mm 0.144 (.139 - .169) 36 mm 0.160 (.156 - .190) 53 mm 0.238 (.228 - .278) 103 mm 0.511 (.447 - .547) ☒ Barometric Pressure Check Gauge ID # 28199 ☒ Stability Checks			Flow Column # ☐ 5L/min - 17mm ☐ 15L/min - 53mm ☐ 30L/min - 103mm ☐ R-Value ☐ Post Calibration Verification (L/s) Flow Column # 32 mm (.139 - .169) 36 mm (.156 - .190) 53 mm (.228 - .278) 103 mm (.447 - .547)																																																																																																	
			<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr></thead><tbody><tr><td>0.050</td><td>MP6286</td><td>202406K 06/19/2026</td></tr><tr><td>0.080</td><td>MP6287</td><td>202406L 06/19/2026</td></tr><tr><td>0.200</td><td>MP6288</td><td>202406N 06/20/2026</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG429602 10/22/2026</td></tr></tbody></table>			Simulator	Serial #	Lot #/Exp	0.050	MP6286	202406K 06/19/2026	0.080	MP6287	202406L 06/19/2026	0.200	MP6288	202406N 06/20/2026	0.080 DGS	N/A	AG429602 10/22/2026	<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>			Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>			Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP6284</td></tr><tr><td>Interferent</td><td>MP6285</td></tr><tr><td>0.050</td><td>MP6286</td></tr><tr><td>0.080</td><td>MP6287</td></tr><tr><td>0.200</td><td>MP6288</td></tr></tbody></table>			Simulator	Serial Number	0.000	MP6284	Interferent	MP6285	0.050	MP6286	0.080	MP6287	0.200	MP6288	<table border="1"><thead><tr><th colspan="2">Attachments</th></tr></thead><tbody><tr><td>☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Calibration Adjustment</td><td>☐ Post-Stability Checks ☐ Flow Calibration ☐ Form 40 ☐ Other</td></tr></tbody></table>			Attachments		☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Calibration Adjustment	☐ Post-Stability Checks ☐ Flow Calibration ☐ Form 40 ☐ Other	<table border="1"><thead><tr><th colspan="2">Maintenance</th></tr></thead><tbody><tr><td>☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other</td><td></td></tr></tbody></table>			Maintenance		☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other	
Simulator	Serial #	Lot #/Exp																																																																																																					
0.050	MP6286	202406K 06/19/2026																																																																																																					
0.080	MP6287	202406L 06/19/2026																																																																																																					
0.200	MP6288	202406N 06/20/2026																																																																																																					
0.080 DGS	N/A	AG429602 10/22/2026																																																																																																					
Simulator	Serial #	Lot #	Expiration																																																																																																				
0.000		N/A	N/A																																																																																																				
0.040																																																																																																							
0.100																																																																																																							
0.200																																																																																																							
0.300																																																																																																							
0.080 DGS	N/A																																																																																																						
Simulator	Serial #	Lot #	Expiration																																																																																																				
0.050																																																																																																							
0.080																																																																																																							
0.200																																																																																																							
0.080 DGS	N/A																																																																																																						
Simulator	Serial Number																																																																																																						
0.000	MP6284																																																																																																						
Interferent	MP6285																																																																																																						
0.050	MP6286																																																																																																						
0.080	MP6287																																																																																																						
0.200	MP6288																																																																																																						
Attachments																																																																																																							
☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Calibration Adjustment	☐ Post-Stability Checks ☐ Flow Calibration ☐ Form 40 ☐ Other																																																																																																						
Maintenance																																																																																																							
☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other																																																																																																							
Calibration Adjustment By			Department Inspection By TDG																																																																																																				
Barometric Pressure Gauge ID #			Barometric Pressure ID# 28199																																																																																																				
Gauge 1016 Instrument 1013			Mouth Alcohol Solution Lot # 2024-A																																																																																																				
Acetone Stock Solution Lot # 2023-B																																																																																																							
Notes/Suggested Service:			☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																																																																																																				
			Digitally signed by Phil Nicodemo Date: 2025.05.29 13:38:55 -0400																																																																																																				
			Shayla Platt Date: 2025.05.31 15:01:33 -0400																																																																																																				
			Tech Review / Date Admin Review / Date																																																																																																				

Stability Checks

0.05g/210L	0.08g/210L	0.20g/210L	DGS 0.08g/210L																																																																																																																																				
✓ 0.047 to 0.053	✓ 0.077 to 0.083	✓ 0.194 to 0.206	✓ 0.077 to 0.083 ✓ ≤0.003 of Wet ✓																																																																																																																																				
HILLSBOROUGH CSO Intoxilyzer - Alcohol Analyzer Model 8000 05/19/2025 Software: 8100.27 SN 80-007483 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>13:59</td></tr> <tr><td>Control Test</td><td>0.048</td><td>14:00</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:00</td></tr> <tr><td>Control Test</td><td>0.048</td><td>14:01</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:01</td></tr> <tr><td>Control Test</td><td>0.048</td><td>14:02</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:02</td></tr> <tr><td>Control Test</td><td>0.048</td><td>14:02</td></tr> <tr><td>Average</td><td>0.0480</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </table>	Air Blank	0.000	13:59	Control Test	0.048	14:00	Air Blank	0.000	14:00	Control Test	0.048	14:01	Air Blank	0.000	14:01	Control Test	0.048	14:02	Air Blank	0.000	14:02	Control Test	0.048	14:02	Average	0.0480		Std Dev	0.0000		Rel Std Dev(%)	0.0000		HILLSBOROUGH CSO Intoxilyzer - Alcohol Analyzer Model 8000 05/19/2025 Software: 8100.27 SN 80-007483 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>13:38</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:39</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:39</td></tr> <tr><td>Control Test</td><td>0.079</td><td>13:40</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:41</td></tr> <tr><td>Control Test</td><td>0.078</td><td>13:41</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:42</td></tr> <tr><td>Control Test</td><td>0.000</td><td>13:42</td></tr> <tr><td>Average</td><td>0.0787</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7339</td><td></td></tr> </table>	Air Blank	0.000	13:38	Control Test	0.079	13:39	Air Blank	0.000	13:39	Control Test	0.079	13:40	Air Blank	0.000	13:41	Control Test	0.078	13:41	Air Blank	0.000	13:42	Control Test	0.000	13:42	Average	0.0787		Std Dev	0.0006		Rel Std Dev(%)	0.7339		HILLSBOROUGH CSO Intoxilyzer - Alcohol Analyzer Model 8000 05/19/2025 Software: 8100.27 SN 80-007483 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>13:47</td></tr> <tr><td>Control Test</td><td>0.199</td><td>13:48</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:48</td></tr> <tr><td>Control Test</td><td>0.198</td><td>13:49</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:50</td></tr> <tr><td>Control Test</td><td>0.198</td><td>13:50</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>13:51</td></tr> <tr><td>Control Test</td><td>0.000</td><td>13:51</td></tr> <tr><td>Average</td><td>0.1983</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2911</td><td></td></tr> </table>	Air Blank	0.000	13:47	Control Test	0.199	13:48	Air Blank	0.000	13:48	Control Test	0.198	13:49	Air Blank	0.000	13:50	Control Test	0.198	13:50	Air Blank	0.000	13:51	Control Test	0.000	13:51	Average	0.1983		Std Dev	0.0006		Rel Std Dev(%)	0.2911		HILLSBOROUGH CSO Intoxilyzer - Alcohol Analyzer Model 8000 05/19/2025 Software: 8100.27 SN 80-007483 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>14:05</td></tr> <tr><td>Control Test</td><td>0.077</td><td>14:05</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:06</td></tr> <tr><td>Control Test</td><td>0.077</td><td>14:06</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:06</td></tr> <tr><td>Control Test</td><td>0.077</td><td>14:07</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>14:07</td></tr> <tr><td>Control Test</td><td>0.000</td><td>14:07</td></tr> <tr><td>Average</td><td>0.0770</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </table>	Air Blank	0.000	14:05	Control Test	0.077	14:05	Air Blank	0.000	14:06	Control Test	0.077	14:06	Air Blank	0.000	14:06	Control Test	0.077	14:07	Air Blank	0.000	14:07	Control Test	0.000	14:07	Average	0.0770		Std Dev	0.0000		Rel Std Dev(%)	0.0000	
Air Blank	0.000	13:59																																																																																																																																					
Control Test	0.048	14:00																																																																																																																																					
Air Blank	0.000	14:00																																																																																																																																					
Control Test	0.048	14:01																																																																																																																																					
Air Blank	0.000	14:01																																																																																																																																					
Control Test	0.048	14:02																																																																																																																																					
Air Blank	0.000	14:02																																																																																																																																					
Control Test	0.048	14:02																																																																																																																																					
Average	0.0480																																																																																																																																						
Std Dev	0.0000																																																																																																																																						
Rel Std Dev(%)	0.0000																																																																																																																																						
Air Blank	0.000	13:38																																																																																																																																					
Control Test	0.079	13:39																																																																																																																																					
Air Blank	0.000	13:39																																																																																																																																					
Control Test	0.079	13:40																																																																																																																																					
Air Blank	0.000	13:41																																																																																																																																					
Control Test	0.078	13:41																																																																																																																																					
Air Blank	0.000	13:42																																																																																																																																					
Control Test	0.000	13:42																																																																																																																																					
Average	0.0787																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.7339																																																																																																																																						
Air Blank	0.000	13:47																																																																																																																																					
Control Test	0.199	13:48																																																																																																																																					
Air Blank	0.000	13:48																																																																																																																																					
Control Test	0.198	13:49																																																																																																																																					
Air Blank	0.000	13:50																																																																																																																																					
Control Test	0.198	13:50																																																																																																																																					
Air Blank	0.000	13:51																																																																																																																																					
Control Test	0.000	13:51																																																																																																																																					
Average	0.1983																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.2911																																																																																																																																						
Air Blank	0.000	14:05																																																																																																																																					
Control Test	0.077	14:05																																																																																																																																					
Air Blank	0.000	14:06																																																																																																																																					
Control Test	0.077	14:06																																																																																																																																					
Air Blank	0.000	14:06																																																																																																																																					
Control Test	0.077	14:07																																																																																																																																					
Air Blank	0.000	14:07																																																																																																																																					
Control Test	0.000	14:07																																																																																																																																					
Average	0.0770																																																																																																																																						
Std Dev	0.0000																																																																																																																																						
Rel Std Dev(%)	0.0000																																																																																																																																						
Operator's Signature 	Operator's Signature 	Operator's Signature 	Operator's Signature 																																																																																																																																				

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: HILLSBOROUGH CSO
Time of Inspection: 13:54

Date of Inspection: 05/20/2025

Serial Number: 80-007483
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	


Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202406K Exp: 06/19/2026	0.08g/210L Test (g/210L) Lot#:202406L Exp: 06/19/2026	0.20g/210L Test (g/210L) Lot#:202406N Exp: 06/20/2026	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG429602 Exp: 10/22/2026
0.000	0.047	0.078	0.196	0.078
0.000	0.047	0.078	0.197	0.078
0.000	0.048	0.078	0.197	0.078
0.000	0.047	0.078	0.197	0.078
0.000	0.047	0.078	0.197	0.077
0.000	0.048	0.078	0.197	0.077
0.000	0.048	0.078	0.197	0.078
0.000	0.048	0.078	0.197	0.077
0.000	0.047	0.078	0.197	0.078
0.000	0.047	0.078	0.197	0.077
Standard Deviations	0.0005	0.0000	0.0003	0.0005

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.


 TAYLOR D GUTSCHOW
 Signature and Printed Name

05/20/2025
 Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
4700 Terminal Drive, Suite 1
Ft. Myers, FL 33907

This is to certify the calibration of Intoxilyzer 8000 serial number 80-007483, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-007483</u>	UNCERTAINTY* $\pm$
Owning Agency:	<u>HILLSBOROUGH CSO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>05/20/2025</u>	0.080 g/ 210 L 0.004
Calibration Time:	<u>13:54</u>	0.200 g/ 210 L 0.007
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.
*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.
This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Taylor
Gutschow
Digitally signed by Taylor
Gutschow
Date: 2025.05.23 11:31:28
-04'00'

05/20/2025

Date

TAYLOR D GUTSCHOW,
Department Inspector

FDLE/ATP Form 69 December 2021

Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality