

Alcohol Testing Program - Instrument Processing Sheet

Agency: PUTNAM COUNTY SO

Instrument Serial Number: 80-001128

Date In: 7/1/2026

DI Completion Date: 7/9/2026

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: <u>WKP</u> Date: <u>7/1/2026</u>		Quality Checks By: <u>SLH</u> Date: <u>7/1/2026</u>		Flow Adjustment By: _____																															
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>212</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP 102</u> 32 mm 0.148 (.139-.169) 36 mm 0.160 (.156-.190) 53 mm 0.230 (.228-.278) 103 mm 0.476 (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>28662</u> Gauge: <u>1013</u> Instrument: <u>1010</u> ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm (.139-.169) 36 mm (.156-.190) 53 mm (.228-.278) 103 mm (.447-.547)																															
						<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 20%;">Serial #</th> <th style="width: 25%;">Lot#/Exp</th> <th style="width: 40%;">Maintenance By: _____ Date: _____</th> </tr> </thead> <tbody> <tr> <td rowspan="2">0.050</td> <td rowspan="2">MP5088</td> <td>202603A</td> <td>☐ Battery Replacement</td> </tr> <tr> <td>3/3/2028</td> <td>☐ Dry Gas Regulator Replacement</td> </tr> <tr> <td rowspan="2">0.080</td> <td rowspan="2">MP5089</td> <td>202603B</td> <td>☐ Tank Sensor Tare</td> </tr> <tr> <td>3/3/2028</td> <td>☐ Breath Tube Replacement</td> </tr> <tr> <td rowspan="2">0.200</td> <td rowspan="2">MP5090</td> <td>202603D</td> <td>☐ Other: _____</td> </tr> <tr> <td>3/4/2026- 28 <i>SLH</i></td> <td></td> </tr> <tr> <td rowspan="2">0.080 DGS</td> <td rowspan="2">N/A</td> <td>23825080A1</td> <td></td> </tr> <tr> <td>10/5/2027</td> <td></td> </tr> </tbody> </table>		Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____	0.050	MP5088	202603A	☐ Battery Replacement	3/3/2028	☐ Dry Gas Regulator Replacement	0.080	MP5089	202603B	☐ Tank Sensor Tare	3/3/2028	☐ Breath Tube Replacement	0.200	MP5090	202603D	☐ Other: _____	3/4/2026- 28 <i>SLH</i>		0.080 DGS	N/A	23825080A1		10/5/2027	
						Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____																										
						0.050	MP5088	202603A	☐ Battery Replacement																										
								3/3/2028	☐ Dry Gas Regulator Replacement																										
0.080	MP5089	202603B	☐ Tank Sensor Tare																																
		3/3/2028	☐ Breath Tube Replacement																																
0.200	MP5090	202603D	☐ Other: _____																																
		3/4/2026- 28 <i>SLH</i>																																	
0.080 DGS	N/A	23825080A1																																	
		10/5/2027																																	
Optical Bench Adjustment By: <u>SLH</u>		Department Inspection By: <u>SLH</u>																																	
Barometric Pressure Gauge: <u>1014</u> ID#: <u>34419</u>		Barometric Pressure ID#: <u>28427</u>																																	
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td>MP6294</td> <td>N/A</td> <td>N/A</td> </tr> <tr> <td>0.040</td> <td>MP6295</td> <td>25090</td> <td>3/11/2027</td> </tr> <tr> <td>0.100</td> <td>MP6296</td> <td>24390</td> <td>10/29/2026</td> </tr> <tr> <td>0.200</td> <td>MP6297</td> <td>25020</td> <td>1/14/2027</td> </tr> <tr> <td>0.300</td> <td>MP6298</td> <td>25150</td> <td>5/6/2027</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>AG526603</td> <td>09/23/2027</td> </tr> </tbody> </table>		Simulator	Serial #	Lot #	Expiration	0.000	MP6294	N/A	N/A	0.040	MP6295	25090	3/11/2027	0.100	MP6296	24390	10/29/2026	0.200	MP6297	25020	1/14/2027	0.300	MP6298	25150	5/6/2027	0.080 DGS	N/A	AG526603	09/23/2027	Gauge: <u>1015</u> Instrument: <u>1015</u> Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>09/25/2027</u> Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>09/22/2027</u>					
Simulator	Serial #	Lot #	Expiration																																
0.000	MP6294	N/A	N/A																																
0.040	MP6295	25090	3/11/2027																																
0.100	MP6296	24390	10/29/2026																																
0.200	MP6297	25020	1/14/2027																																
0.300	MP6298	25150	5/6/2027																																
0.080 DGS	N/A	AG526603	09/23/2027																																
☒ Post Optical Bench Adjustment Stability Checks		<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot #</th> <th>Expiration</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td>MP5088</td> <td>202603A</td> <td>03/03/2028</td> </tr> <tr> <td>0.080</td> <td>MP5089</td> <td>202603B</td> <td>03/03/2028</td> </tr> <tr> <td>0.200</td> <td>MP5090</td> <td>202603D</td> <td>03/04/2028</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td>23825080A1</td> <td>10/5/2027</td> </tr> </tbody> </table>		Simulator	Serial #	Lot #	Expiration	0.050	MP5088	202603A	03/03/2028	0.080	MP5089	202603B	03/03/2028	0.200	MP5090	202603D	03/04/2028	0.080 DGS	N/A	23825080A1	10/5/2027												
Simulator	Serial #	Lot #	Expiration																																
0.050	MP5088	202603A	03/03/2028																																
0.080	MP5089	202603B	03/03/2028																																
0.200	MP5090	202603D	03/04/2028																																
0.080 DGS	N/A	23825080A1	10/5/2027																																
Gauge ID #: <u>28662</u> Gauge: <u>1014</u> Instrument: <u>1014</u>		Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☒ Optical Bench Adjustment																																	
Notes/Suggested Service: Flow verification performed 7/7/2026 and post stabilities performed 2 separate days. SLH 7/7/2026 Corrected expiration date for the 0.200 g/210L stability. SLH 7/9/2026		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC																																	
		☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use																																	
		☒ Conduct an Agency Inspection Before Evidentiary Use																																	
Taylor Gutschow Digitally signed by Taylor Gutschow Date: 2026.07.16 10:50:39 -04'00'		Kaitlyn Spearin Digitally signed by Kaitlyn Spearin Date: 2026.07.21 10:20:26 -04'00'																																	
Tech Review		Admin Review																																	

Florida Department of Law Enforcement Alcohol Testing Program

AGENCY INSPECTION REPORT - INTOXILYZER 8000

Agency: PUTNAM COUNTY SO
Time of Inspection: 16:13

Date of Inspection: 07/01/2026

Serial Number: 80-001128
Software: 8100.27

Check or Test	YES	NO
Date and/or Time Adjusted		No
Diagnostic Check (Pre-Inspection): OK		No
Alcohol Free Subject Test: 0.000		No
Mouth Alcohol Test: Slope Not Met		No
Interferent Detect Test: Interferent Detect		No
Diagnostic Check (Post-Inspection): OK		No

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:_____ Exp:_____	0.08g/210L Test (g/210L) Lot#:_____ Exp:_____	0.20g/210L Test (g/210L) Lot#:_____ Exp:_____	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:_____ Exp:_____


Number of Simulators Used: _____

Remarks:

BYPASS AI FOR OPERATION, COMPLIANCE UNDETERMINED

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I hold a valid Florida Department of Law Enforcement Agency Inspector Permit and that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

 LEANDRA HIGGINBOTHAM
Signature and Printed Name

07/01/2026
Date

Stability Checks

80-00 1128

0.050 g/210L 0.047 to 0.053 g/210L	0.080 g/210L 0.077 to 0.083 g/210L	0.200 g/210L 0.194 to 0.206 g/210L	DGS 0.080 g/210L 0.077 to 0.083 g/210L ≤0.003 g/210L of Wet
PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001128 07/01/2026 Software: 8100.27 Test g/210L Time Air Blank 0.000 16:19 Control Test 0.049 16:20 Air Blank 0.000 16:21 Control Test 0.048 16:21 Air Blank 0.000 16:22 Control Test 0.049 16:23 Air Blank 0.000 16:23 Control Test Stats Average 0.0487 Std Dev 0.0006 Rel Std Dev(%) 1.1863 <i>Miguel</i> Operator's Signature	PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001128 07/01/2026 Software: 8100.27 Test g/210L Time Air Blank 0.000 16:24 Control Test 0.078 16:25 Air Blank 0.000 16:25 Control Test 0.078 16:26 Air Blank 0.000 16:27 Control Test 0.078 16:28 Air Blank 0.000 16:28 Control Test Stats Average 0.0780 Std Dev 0.0000 Rel Std Dev(%) 0.0000 <i>Miguel</i> Operator's Signature	PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 3000 SN 80-001128 07/01/2026 Software: 8100.27 Test g/210L Time Air Blank 0.000 16:15 Control Test 0.199 16:15 Air Blank 0.000 16:16 Control Test 0.197 16:17 Air Blank 0.000 16:17 Control Test 0.197 16:18 Air Blank 0.000 16:19 Control Test Stats Average 0.1977 Std Dev 0.0012 Rel Std Dev(%) 0.5842 <i>Miguel</i> Operator's Signature	PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001128 07/01/2026 Software: 8100.27 Test g/210L Time Air Blank 0.000 16:29 Control Test 0.082 16:29 Air Blank 0.000 16:30 Control Test 0.082 16:30 Air Blank 0.000 16:31 Control Test 0.082 16:31 Air Blank 0.000 16:32 Control Test Stats Average 0.0820 Std Dev 0.0000 Rel Std Dev(%) 0.0000 <i>Miguel</i> Operator's Signature

OBA=optical bench adjustment

needs OBA

SUB

Verified instrument was able to read gas pressure.

PUTNAM COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001128
07/07/2026 13:56:18

SW

Auto Calibration

pg 1 of 2

<<<<< 3um >>>>> <<<<< 9um >>>>>

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1
Sample % Abs (% Abs Ref) % Abs (% Abs Ref)
Sample #1 0.1130 (0.0000) 0.0920 (-0.0210)
Sample #2 0.1100 (0.0140) 0.1040 (-0.0190)
Sample #3 0.1130 (0.0260) 0.1030 (-0.0160)
Sample #4 0.1130 (0.0340) 0.1400 (-0.0260)
Avg % Abs 0.1120 (0.0247) 0.1157 (-0.0203)
STD DEV 0.0017 (0.0101) 0.0211 (0.0051)
REL STD DEV 1.546 (40.810) 18.224 (25.237)

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1
Sample % Abs (% Abs Ref) % Abs (% Abs Ref)
Sample #1 0.8550 (-0.0370) 1.5620 (-0.0080)
Sample #2 0.8490 (-0.0170) 1.5610 (0.0070)
Sample #3 0.8330 (-0.0010) 1.5430 (0.0080)
Sample #4 0.8370 (0.0180) 1.5430 (0.0250)
Avg % Abs 0.8397 (-0.0000) 1.5490 (0.0133)
STD DEV 0.0083 (0.0175) 0.0104 (0.0101)
REL STD DEV 0.992 (2822021376.000) 0.671 (75.870)

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1
Sample % Abs (% Abs Ref) % Abs (% Abs Ref)
Sample #1 1.9140 (-0.0060) 3.6640 (-0.0050)
Sample #2 1.9210 (0.0010) 3.6630 (0.0180)
Sample #3 1.9090 (0.0280) 3.6700 (0.0180)
Sample #4 1.9310 (0.0290) 3.6450 (0.0380)
Avg % Abs 1.9203 (0.0193) 3.6593 (0.0247)
STD DEV 0.0110 (0.0159) 0.0129 (0.0115)
REL STD DEV 0.574 (82.164) 0.352 (46.812)

Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1
Sample % Abs (% Abs Ref) % Abs (% Abs Ref)
Sample #1 3.7190 (-0.0060) 7.0380 (0.0070)
Sample #2 3.7100 (0.0190) 7.0310 (0.0410)
Sample #3 3.7030 (0.0260) 7.0260 (0.0320)
Sample #4 3.7100 (0.0340) 7.0370 (0.0390)
Avg % Abs 3.7077 (0.0263) 7.0313 (0.0373)
STD DEV 0.0040 (0.0075) 0.0055 (0.0047)
REL STD DEV 0.109 (28.502) 0.078 (12.658)

Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1
Sample % Abs (% Abs Ref) % Abs (% Abs Ref)
Sample #1 5.4260 (-0.0120) 10.2570 (-0.0160)
Sample #2 5.3790 (0.0330) 10.1690 (0.0600)
Sample #3 5.3720 (0.0710) 10.1770 (0.0810)
Sample #4 5.3510 (0.0610) 10.1300 (0.0800)
Avg % Abs 5.3673 (0.0550) 10.1587 (0.0737)
STD DEV 0.0146 (0.0197) 0.0251 (0.0118)
REL STD DEV 0.271 (35.814) 0.248 (16.081)

PUTNAM COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001128
07/07/2026 13:56:18

Auto Calibration

pg 2 of 2

<<<<< 3um >>>>>

Zero Order Coef -259.89
First Order Coef 2531.59
Second Order Coef 32.59

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.001	-0.0005
0.040	0.040	0.0003
0.100	0.099	0.0008
0.200	0.201	-0.0011
0.300	0.300	0.0004

<<<<< 9um >>>>>

-132.71
1284.84
13.13

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0003
0.040	0.040	0.0003
0.100	0.100	0.0004
0.200	0.201	-0.0006
0.300	0.300	0.0002

<<<<< 3um >>>>> <<<<< 9um >>>>>

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
Sample
Sample #1 3047.00 3416.00
Sample #2 3011.00 3438.00
Sample #3 3032.00 3405.00
Sample #4 3049.00 3408.00
Avg 3030.6667 3417.0000
STD DEV 19.0351 18.2483
REL STD DEV 0.628 0.534
H2O adjust (mg/l*10k) 779 392


Barometric Pressure = 1014

*****CALIBRATION SUCCESSFUL*****

SWA

Post Optical Bench Adjustment Stability Checks

80-001128
Sub

0.050 g/210L	0.080 g/210L	0.200 g/210L	DGS 0.080 g/210L																																																																																																																																																
0.047 to 0.053 g/210L ✓	0.077 to 0.083 g/210L ✓	0.194 to 0.206 g/210L ✓	0.077 to 0.083 g/210L 0.003 g/210L of Wet ✓																																																																																																																																																
PUTNAM COUNTY SC Intoxilyzer - Alcohol Analyzer Model 3000 SN 80-001128 07/09/2026 Software: 8100.27	PUTNAM COUNTY SC Intoxilyzer - Alcohol Analyzer Model 3000 SN 80-001128 07/09/2026 Software: 8100.27	PUTNAM COUNTY SC Intoxilyzer - Alcohol Analyzer Model 3000 SN 80-001128 07/09/2026 Software: 8100.27	PUTNAM COUNTY SC Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001128 07/07/2026 Software: 8100.27																																																																																																																																																
<table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>09:14</td></tr> <tr><td>Control Test</td><td>0.050</td><td>09:15</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:15</td></tr> <tr><td>Control Test</td><td>0.050</td><td>09:16</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:17</td></tr> <tr><td>Control Test</td><td>0.050</td><td>09:17</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:18</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0500</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	09:14	Control Test	0.050	09:15	Air Blank	0.000	09:15	Control Test	0.050	09:16	Air Blank	0.000	09:17	Control Test	0.050	09:17	Air Blank	0.000	09:18	Control Test Stats			Average	0.0500		Std Dev	0.0000		Rel Std Dev(%)	0.0000		<table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>15:31</td></tr> <tr><td>Control Test</td><td>0.079</td><td>15:31</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:32</td></tr> <tr><td>Control Test</td><td>0.079</td><td>15:33</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:33</td></tr> <tr><td>Control Test</td><td>0.079</td><td>15:34</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:35</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0790</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	15:31	Control Test	0.079	15:31	Air Blank	0.000	15:32	Control Test	0.079	15:33	Air Blank	0.000	15:33	Control Test	0.079	15:34	Air Blank	0.000	15:35	Control Test Stats			Average	0.0790		Std Dev	0.0000		Rel Std Dev(%)	0.0000		<table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>09:07</td></tr> <tr><td>Control Test</td><td>0.199</td><td>09:08</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:09</td></tr> <tr><td>Control Test</td><td>0.198</td><td>09:09</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:10</td></tr> <tr><td>Control Test</td><td>0.198</td><td>09:11</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>09:11</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.1983</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2911</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	09:07	Control Test	0.199	09:08	Air Blank	0.000	09:09	Control Test	0.198	09:09	Air Blank	0.000	09:10	Control Test	0.198	09:11	Air Blank	0.000	09:11	Control Test Stats			Average	0.1983		Std Dev	0.0006		Rel Std Dev(%)	0.2911		<table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>15:37</td></tr> <tr><td>Control Test</td><td>0.081</td><td>15:38</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:38</td></tr> <tr><td>Control Test</td><td>0.081</td><td>15:38</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:39</td></tr> <tr><td>Control Test</td><td>0.081</td><td>15:39</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:40</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0810</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table>	Test	g/210L	Time	Air Blank	0.000	15:37	Control Test	0.081	15:38	Air Blank	0.000	15:38	Control Test	0.081	15:38	Air Blank	0.000	15:39	Control Test	0.081	15:39	Air Blank	0.000	15:40	Control Test Stats			Average	0.0810		Std Dev	0.0000		Rel Std Dev(%)	0.0000	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:14																																																																																																																																																	
Control Test	0.050	09:15																																																																																																																																																	
Air Blank	0.000	09:15																																																																																																																																																	
Control Test	0.050	09:16																																																																																																																																																	
Air Blank	0.000	09:17																																																																																																																																																	
Control Test	0.050	09:17																																																																																																																																																	
Air Blank	0.000	09:18																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0500																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	15:31																																																																																																																																																	
Control Test	0.079	15:31																																																																																																																																																	
Air Blank	0.000	15:32																																																																																																																																																	
Control Test	0.079	15:33																																																																																																																																																	
Air Blank	0.000	15:33																																																																																																																																																	
Control Test	0.079	15:34																																																																																																																																																	
Air Blank	0.000	15:35																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0790																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	09:07																																																																																																																																																	
Control Test	0.199	09:08																																																																																																																																																	
Air Blank	0.000	09:09																																																																																																																																																	
Control Test	0.198	09:09																																																																																																																																																	
Air Blank	0.000	09:10																																																																																																																																																	
Control Test	0.198	09:11																																																																																																																																																	
Air Blank	0.000	09:11																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.1983																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.2911																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	15:37																																																																																																																																																	
Control Test	0.081	15:38																																																																																																																																																	
Air Blank	0.000	15:38																																																																																																																																																	
Control Test	0.081	15:38																																																																																																																																																	
Air Blank	0.000	15:39																																																																																																																																																	
Control Test	0.081	15:39																																																																																																																																																	
Air Blank	0.000	15:40																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0810																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
 Operator's Signature	 Operator's Signature	 Operator's Signature	 Operator's Signature																																																																																																																																																

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: PUTNAM COUNTY SO
Time of Inspection: 11:45

Date of Inspection: 07/09/2026

Serial Number: 80-001128
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23825080A1 Exp: 10/05/2027
0.000	0.049	0.079	0.198	0.080
0.000	0.049	0.079	0.198	0.080
0.000	0.050	0.079	0.198	0.080
0.000	0.050	0.079	0.198	0.080
0.000	0.050	0.078	0.198	0.080
0.000	0.050	0.079	0.198	0.080
0.000	0.050	0.079	0.198	0.080
0.000	0.050	0.079	0.199	0.080
0.000	0.050	0.079	0.198	0.080
0.000	0.050	0.079	0.199	0.080


Standard Deviations	0.0004	0.0003	0.0004	0.0000
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0002 Number of Simulators Used: 5


Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.


LEANDRA HIGGINBOTHAM
Signature and Printed Name

07/09/2026
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001128, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001128</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>PUTNAM COUNTY SO</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>07/09/2026</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>11:45</u>	0.200 g/ 210 L	0.008
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.



07/09/2026

Date

LEANDRA HIGGINBOTHAM,

Department Inspector

FDLE/ATP Form 69 January 2026
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality