



INSTRUMENT PROCESSING SHEET

Agency TALLAHASSEE PDS/N 80-002989Florida Department of
Law EnforcementDate In 9/10/25DI Completion Date 09/16/2025☐ Ship ☒ P/U ☐ H/D ☐ CMI ☐ EE

Intake By <u>KTS</u> Date <u>9/10/25</u>		Quality Checks By <u>SLH</u> Date <u>9/15/2025</u>		Flow Calibration By _____ Date _____																													
☒ Annual ☐ Registration ☐ Return from CMI / EE Visual Inspection: ☒ Case ☐ Handle ☒ Keyboard ☐ Dry Gas Shelf ☒ Feet ☐ Breath Tube ☒ Ports ☐ Screws Tight Other Equipment/ Accessories: ☐ Power cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes: <u>Dropped by Agency. SLH</u> <u>9/15/25</u>		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value <u>219</u> ☒ Flow Verification (L/s) Flow Column # <u>ATP103</u> 32 mm <u>0.148</u> (.139 - .169) 36 mm <u>0.160</u> (.156 - .190) 53 mm <u>0.230</u> (.228 - .278) 103 mm <u>0.492</u> (.447 - .547) ☒ Barometric Pressure Check Gauge ID # <u>28421</u> ☒ Stability Checks		☐ Flow Column # _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value _____ ☐ Post Calibration Verification (L/s) Flow Column # _____ 32 mm _____ (.139 - .169) 36 mm _____ (.156 - .190) 53 mm _____ (.228 - .278) 103 mm _____ (.447 - .547)																													
		<table border="1" style="width:100%"><tr><th>Simulator</th><th>Serial #</th><th>Lot #/Exp</th></tr><tr><td>0.050</td><td>MP6291</td><td>202406K 06/19/2026</td></tr><tr><td>0.080</td><td>MP6292</td><td>202406L 06/19/2026</td></tr><tr><td>0.200</td><td>MP6293</td><td>202406N 06/20/2026</td></tr><tr><td>0.080 DGS</td><td>N/A</td><td>AG429602 10/22/2026</td></tr></table>		Simulator	Serial #	Lot #/Exp	0.050	MP6291	202406K 06/19/2026	0.080	MP6292	202406L 06/19/2026	0.200	MP6293	202406N 06/20/2026	0.080 DGS	N/A	AG429602 10/22/2026	Maintenance By _____ Date _____ ☐ Battery Replacement ☐ Dry Gas Regulator Replacement ☐ Breath Tube Replacement ☐ Other _____														
Simulator	Serial #	Lot #/Exp																															
0.050	MP6291	202406K 06/19/2026																															
0.080	MP6292	202406L 06/19/2026																															
0.200	MP6293	202406N 06/20/2026																															
0.080 DGS	N/A	AG429602 10/22/2026																															
Calibration Adjustment By _____		Department Inspection By <u>SLH</u>																															
Barometric Pressure Gauge _____ ID # _____		Barometric Pressure ID# <u>28421</u>																															
<table border="1" style="width:100%"><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Gauge <u>1015</u> Instrument <u>1014</u> Mouth Alcohol Solution Lot # <u>2025-A</u> Acetone Stock Solution Lot # <u>2024-B</u>			
Simulator	Serial #	Lot #	Expiration																														
0.000		N/A	N/A																														
0.040																																	
0.100																																	
0.200																																	
0.300																																	
0.080 DGS	N/A																																
☐ Post Calibration Adjustment Stability Checks		<table border="1" style="width:100%"><tr><th>Simulator</th><th>Serial Number</th></tr><tr><td>0.000</td><td>MP6294</td></tr><tr><td>Interferent</td><td>MP6290</td></tr><tr><td>0.050</td><td>MP5088</td></tr><tr><td>0.080</td><td>MP6292</td></tr><tr><td>0.200</td><td>MP6293</td></tr></table>				Simulator	Serial Number	0.000	MP6294	Interferent	MP6290	0.050	MP5088	0.080	MP6292	0.200	MP6293																
Simulator	Serial Number																																
0.000	MP6294																																
Interferent	MP6290																																
0.050	MP5088																																
0.080	MP6292																																
0.200	MP6293																																
<table border="1" style="width:100%"><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></table>		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			Attachments ☒ Form 41 ☐ Post-Stability Checks ☒ Stability Checks ☐ Flow Calibration ☒ Calibration Certificate ☐ Form 40 ☐ Calibration Adjustment ☐ Other _____											
Simulator	Serial #	Lot #	Expiration																														
0.050																																	
0.080																																	
0.200																																	
0.080 DGS	N/A																																
Notes/Suggested Service: _____		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																															
		Shayla Platt Digitally signed by Shayla Platt Date: 2025.09.17 13:27:20 -04'00'																															
		Taylor Gutschow Digitally signed by Taylor Gutschow Date: 2025.09.18 09:31:08 -04'00'																															
		Tech Review / Date _____ Admin Review / Date _____																															

Stability Checks

80-002989 sub
9/15/25

0.050 g/210L	0.080 g/210L	0.200 g/210L	DGS 0.080 g/210L																																																																																																																																				
0.047 to 0.053 g/210L ☒ 0.077 to 0.083 g/210L ☒ 0.194 to 0.206 g/210L ☒ 0.077 to 0.083 g/210L ☒ 50.003 g/210L of Wet ☒	0.077 to 0.083 g/210L ☒ 0.194 to 0.206 g/210L ☒ 0.077 to 0.083 g/210L ☒ 50.003 g/210L of Wet ☒	0.077 to 0.083 g/210L ☒ 0.194 to 0.206 g/210L ☒ 0.077 to 0.083 g/210L ☒ 50.003 g/210L of Wet ☒	0.077 to 0.083 g/210L ☒ 0.194 to 0.206 g/210L ☒ 0.077 to 0.083 g/210L ☒ 50.003 g/210L of Wet ☒																																																																																																																																				
Performed Root Case Analysis	Performed Root Case Analysis	Performed Root Case Analysis	Performed Root Case Analysis																																																																																																																																				
TALLAHASSEE PD Intoxilyzer - Alcohol Analyzer Model 8000 09/15/2025 Software: 8100.27 SN 80-002989 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>11:30</td></tr> <tr><td>Control Test</td><td>0.049</td><td>11:30</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:31</td></tr> <tr><td>Control Test</td><td>0.049</td><td>11:32</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:32</td></tr> <tr><td>Control Test</td><td>0.049</td><td>11:33</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:33</td></tr> <tr><td>Control Test Status</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0490</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </table>	Air Blank	0.000	11:30	Control Test	0.049	11:30	Air Blank	0.000	11:31	Control Test	0.049	11:32	Air Blank	0.000	11:32	Control Test	0.049	11:33	Air Blank	0.000	11:33	Control Test Status			Average	0.0490		Std Dev	0.0000		Rel Std Dev(%)	0.0000		TALLAHASSEE PD Intoxilyzer - Alcohol Analyzer Model 8000 09/15/2025 Software: 8100.27 SN 80-002989 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>11:35</td></tr> <tr><td>Control Test</td><td>0.079</td><td>11:36</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:37</td></tr> <tr><td>Control Test</td><td>0.079</td><td>11:38</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:39</td></tr> <tr><td>Control Test</td><td>0.079</td><td>11:39</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:39</td></tr> <tr><td>Control Test Status</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0790</td><td></td></tr> <tr><td>Std Dev</td><td>0.0010</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0010</td><td></td></tr> </table>	Air Blank	0.000	11:35	Control Test	0.079	11:36	Air Blank	0.000	11:37	Control Test	0.079	11:38	Air Blank	0.000	11:39	Control Test	0.079	11:39	Air Blank	0.000	11:39	Control Test Status			Average	0.0790		Std Dev	0.0010		Rel Std Dev(%)	0.0010		TALLAHASSEE PD Intoxilyzer - Alcohol Analyzer Model 8000 09/15/2025 Software: 8100.27 SN 80-002989 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>11:44</td></tr> <tr><td>Control Test</td><td>0.199</td><td>11:45</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:45</td></tr> <tr><td>Control Test</td><td>0.198</td><td>11:46</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:47</td></tr> <tr><td>Control Test</td><td>0.199</td><td>11:47</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>11:48</td></tr> <tr><td>Control Test Status</td><td></td><td></td></tr> <tr><td>Average</td><td>0.1987</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2906</td><td></td></tr> </table>	Air Blank	0.000	11:44	Control Test	0.199	11:45	Air Blank	0.000	11:45	Control Test	0.198	11:46	Air Blank	0.000	11:47	Control Test	0.199	11:47	Air Blank	0.000	11:48	Control Test Status			Average	0.1987		Std Dev	0.0006		Rel Std Dev(%)	0.2906		TALLAHASSEE PD Intoxilyzer - Alcohol Analyzer Model 8000 09/15/2025 Software: 8100.27 SN 80-002989 Test g/210L Time <table border="1"> <tr><td>Air Blank</td><td>0.000</td><td>12:01</td></tr> <tr><td>Control Test</td><td>0.078</td><td>12:01</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:02</td></tr> <tr><td>Control Test</td><td>0.079</td><td>12:02</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:03</td></tr> <tr><td>Control Test</td><td>0.080</td><td>12:03</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:03</td></tr> <tr><td>Control Test Status</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0790</td><td></td></tr> <tr><td>Std Dev</td><td>0.0010</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.2658</td><td></td></tr> </table>	Air Blank	0.000	12:01	Control Test	0.078	12:01	Air Blank	0.000	12:02	Control Test	0.079	12:02	Air Blank	0.000	12:03	Control Test	0.080	12:03	Air Blank	0.000	12:03	Control Test Status			Average	0.0790		Std Dev	0.0010		Rel Std Dev(%)	1.2658	
Air Blank	0.000	11:30																																																																																																																																					
Control Test	0.049	11:30																																																																																																																																					
Air Blank	0.000	11:31																																																																																																																																					
Control Test	0.049	11:32																																																																																																																																					
Air Blank	0.000	11:32																																																																																																																																					
Control Test	0.049	11:33																																																																																																																																					
Air Blank	0.000	11:33																																																																																																																																					
Control Test Status																																																																																																																																							
Average	0.0490																																																																																																																																						
Std Dev	0.0000																																																																																																																																						
Rel Std Dev(%)	0.0000																																																																																																																																						
Air Blank	0.000	11:35																																																																																																																																					
Control Test	0.079	11:36																																																																																																																																					
Air Blank	0.000	11:37																																																																																																																																					
Control Test	0.079	11:38																																																																																																																																					
Air Blank	0.000	11:39																																																																																																																																					
Control Test	0.079	11:39																																																																																																																																					
Air Blank	0.000	11:39																																																																																																																																					
Control Test Status																																																																																																																																							
Average	0.0790																																																																																																																																						
Std Dev	0.0010																																																																																																																																						
Rel Std Dev(%)	0.0010																																																																																																																																						
Air Blank	0.000	11:44																																																																																																																																					
Control Test	0.199	11:45																																																																																																																																					
Air Blank	0.000	11:45																																																																																																																																					
Control Test	0.198	11:46																																																																																																																																					
Air Blank	0.000	11:47																																																																																																																																					
Control Test	0.199	11:47																																																																																																																																					
Air Blank	0.000	11:48																																																																																																																																					
Control Test Status																																																																																																																																							
Average	0.1987																																																																																																																																						
Std Dev	0.0006																																																																																																																																						
Rel Std Dev(%)	0.2906																																																																																																																																						
Air Blank	0.000	12:01																																																																																																																																					
Control Test	0.078	12:01																																																																																																																																					
Air Blank	0.000	12:02																																																																																																																																					
Control Test	0.079	12:02																																																																																																																																					
Air Blank	0.000	12:03																																																																																																																																					
Control Test	0.080	12:03																																																																																																																																					
Air Blank	0.000	12:03																																																																																																																																					
Control Test Status																																																																																																																																							
Average	0.0790																																																																																																																																						
Std Dev	0.0010																																																																																																																																						
Rel Std Dev(%)	1.2658																																																																																																																																						
<i>M. Magid</i> Operator's Signature	<i>M. Magid</i> Operator's Signature	<i>M. Magid</i> Operator's Signature	<i>M. Magid</i> Operator's Signature																																																																																																																																				

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: TALLAHASSEE PD
Time of Inspection: 17:44

Date of Inspection: 09/16/2025

Serial Number: 80-002989
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202406K Exp: 06/19/2026	0.08g/210L Test (g/210L) Lot#:202406L Exp: 06/19/2026	0.20g/210L Test (g/210L) Lot#:202406N Exp: 06/20/2026	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:AG429602 Exp: 10/22/2026
0.000	0.048	0.078	0.196	0.080
0.000	0.048	0.078	0.196	0.080
0.000	0.048	0.079	0.197	0.080
0.000	0.048	0.078	0.196	0.080
0.000	0.048	0.078	0.196	0.080
0.000	0.048	0.077	0.196	0.080
0.000	0.048	0.077	0.196	0.080
0.000	0.048	0.078	0.197	0.080
0.000	0.048	0.078	0.197	0.080
0.000	0.048	0.078	0.196	0.080


Standard Deviations	0.0000	0.0005	0.0004	0.0000
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0002 Number of Simulators Used: 5



Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.


LEANDRA HIGGINBOTHAM
Signature and Printed Name

09/16/2025
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-002989, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-002989</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>TALLAHASSEE PD</u>	0.050 g/ 210 L	0.004
Calibration Date:	<u>09/16/2025</u>	0.080 g/ 210 L	0.004
Calibration Time:	<u>17:44</u>	0.200 g/ 210 L	0.007
		0.080 g/ 210 L Dry Gas Control	0.005

All results are reported in g/ 210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$)

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards. This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

FDLE/ATP Form 69 October 2024

Issuing Authority: Alcohol Testing Program

09/16/2025

Date

LEANDRA HIGGINBOTHAM,

Department Inspector

Service • Integrity • Respect • Quality