

Alcohol Testing Program - Instrument Processing Sheet

Agency: PUTNAM COUNTY SO

Instrument Serial Number: 80-001129

Date In: 6/19/2026

DI Completion Date: 06/25/2026

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: <u>WKP</u> Date: <u>6/19/2026</u>		Quality Checks By: <u>SLH</u> Date: <u>6/23/2026</u>		Flow Adjustment By: _____																																						
☒ Annual ☐ Dropped Off ☐ Registration ☒ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>188</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP 105</u> 32 mm <u>0.156</u> (.139-.169) 36 mm <u>0.171</u> (.156-.190) 53 mm <u>0.238</u> (.228-.278) 103 mm <u>0.500</u> (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>34419</u> Gauge: <u>1016</u> Instrument: <u>1018</u> ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm (.139-.169) 36 mm (.156-.190) 53 mm (.228-.278) 103 mm (.447-.547)																																						
						<table border="1" style="width:100%; border-collapse: collapse;"> <thead> <tr> <th>Simulator</th> <th>Serial #</th> <th>Lot#/Exp</th> <th>Maintenance By:</th> <th>Date:</th> </tr> </thead> <tbody> <tr> <td rowspan="2">0.050</td> <td rowspan="2">MP5088</td> <td>202603A</td> <td rowspan="2">☐ Battery Replacement</td> <td rowspan="2"></td> </tr> <tr> <td>3/3/2028</td> </tr> <tr> <td rowspan="2">0.080</td> <td rowspan="2">MP5089</td> <td>202603B</td> <td rowspan="2">☐ Dry Gas Regulator Replacement</td> <td rowspan="2"></td> </tr> <tr> <td>3/3/2028</td> </tr> <tr> <td rowspan="2">0.200</td> <td rowspan="2">MP5090</td> <td>202603D</td> <td rowspan="2">☐ Tank Sensor Tare</td> <td rowspan="2"></td> </tr> <tr> <td>3/4/2028</td> </tr> <tr> <td rowspan="2">0.080 DGS</td> <td rowspan="2">N/A</td> <td>23825080A1</td> <td rowspan="2">☐ Breath Tube Replacement</td> <td rowspan="2"></td> </tr> <tr> <td>10/5/2027</td> </tr> <tr> <td colspan="2"></td> <td colspan="2"></td> <td colspan="2">☐ Other:</td> </tr> </tbody> </table>		Simulator	Serial #	Lot#/Exp	Maintenance By:	Date:	0.050	MP5088	202603A	☐ Battery Replacement		3/3/2028	0.080	MP5089	202603B	☐ Dry Gas Regulator Replacement		3/3/2028	0.200	MP5090	202603D	☐ Tank Sensor Tare		3/4/2028	0.080 DGS	N/A	23825080A1	☐ Breath Tube Replacement		10/5/2027					☐ Other:	
						Simulator	Serial #	Lot#/Exp	Maintenance By:	Date:																																
						0.050	MP5088	202603A	☐ Battery Replacement																																	
								3/3/2028																																		
0.080	MP5089	202603B	☐ Dry Gas Regulator Replacement																																							
		3/3/2028																																								
0.200	MP5090	202603D	☐ Tank Sensor Tare																																							
		3/4/2028																																								
0.080 DGS	N/A	23825080A1	☐ Breath Tube Replacement																																							
		10/5/2027																																								
				☐ Other:																																						
Optical Bench Adjustment By: <u>SLH</u>		Department Inspection By: <u>SLH</u>																																								
Barometric Pressure Gauge: <u>1014/1013*</u> ID#: <u>28662</u>		Barometric Pressure ID#: <u>34419</u>																																								
Simulator	Serial #	Lot #	Expiration	Gauge: <u>1014</u>	Instrument: <u>1013</u>																																					
0.000	MP6294	N/A	N/A	Mouth Alcohol Solution Lot #: <u>2025-D</u>	Exp: <u>09/25/2027</u>																																					
0.040	MP6295	25090	3/11/2027	Acetone Stock Solution Lot #: <u>2025-B</u>	Exp: <u>09/22/2027</u>																																					
0.100	MP6296	24390	10/29/2026	Simulator	Serial Number																																					
0.200	MP6297	25020	1/14/2027	0.000	MP6289																																					
0.300	MP6298	24430	12/10/2026	Interferent	MP5087																																					
0.080 DGS	N/A	AG526603	9/23/2027	0.050	MP6291																																					
☒ Post Optical Bench Adjustment Stability Checks				0.080	MP5089																																					
Simulator	Serial #	Lot #	Expiration	0.200	MP5090																																					
0.050	MP5088	202603A	3/3/2028	Attachments																																						
0.080	MP6289	202603B	3/3/2028	☒ Form 41	☒ Post-Stability Checks																																					
0.200	MP5090	202603D	3/4/2028	☒ Stability Checks	☐ Flow Adjustment																																					
0.080 DGS	N/A	23825080A1	10/5/2027	☒ Calibration Certificate	☐ Form 40																																					
Gauge ID #: <u>34419</u> Gauge: <u>1014/1013*</u> Instrument: <u>1012</u>				☒ Optical Bench Adjustment	☐ Other:																																					
Notes/Suggested Service: *A second OBA (optical bench adjustment) was performed due to post-stability results. All the same equipment and solutions used were the same. The pre and post barometric pressure gauge readings were the only difference. SLH 6/23/2026				☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC																																						
				☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use																																						
				☒ Conduct an Agency Inspection Before Evidentiary Use																																						
				Wen-Chi Pierson Digitally signed by Wen-Chi Pierson Date: 2026.06.30 13:39:55 -04'00'																																						
				Kaitlyn Spearin Digitally signed by Kaitlyn Spearin Date: 2026.07.15 13:48:57 -04'00'																																						
				Tech Review Admin Review																																						

Stability Checks

80-001129

0.050 g/210L 0.047 to 0.053 g/210L	0.080 g/210L 0.077 to 0.083 g/210L	0.200 g/210L 0.194 to 0.206 g/210L	DGS 0.080 g/210L 0.077 to 0.083 g/210L 50.003 g/210L of Wet
PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 06/23/2026 SN 80-001129 Software: 8100.27 Test g/210L Time Air Blank 0.000 10:19 Control Test 0.048 10:19 Air Blank 0.000 10:20 Control Test 0.048 10:20 Air Blank 0.000 10:21 Control Test 0.048 10:22 Air Blank 0.000 10:22 Control Test Stats Average 0.0480 Std Dev 0.0000 Rel Std Dev(%) 0.0000 <i>M. J. Jester</i> Operator's Signature	PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 06/23/2026 SN 80-001129 Software: 8100.27 Test g/210L Time Air Blank 0.000 10:24 Control Test 0.075 10:24 Air Blank 0.000 10:25 Control Test 0.076 10:25 Air Blank 0.000 10:26 Control Test 0.076 10:26 Air Blank 0.000 10:27 Control Test Stats Average 0.0757 Std Dev 0.0006 Rel Std Dev(%) 0.7630 <i>M. J. Jester</i> Operator's Signature	PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 06/23/2026 SN 81-001129 Software: 8100.27 Test g/210L Time Air Blank 0.000 10:13 Control Test 0.196 10:13 Air Blank 0.000 10:14 Control Test 0.196 10:15 Air Blank 0.000 10:15 Control Test 0.197 10:16 Air Blank 0.000 10:16 Control Test Stats Average 0.1963 Std Dev 0.0016 Rel Std Dev(%) 0.2941 <i>M. J. Jester</i> Operator's Signature	PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 06/23/2026 SN 80-001129 Software: 8100.27 Test g/210L Time Air Blank 0.000 10:44 Control Test 0.081 10:44 Air Blank 0.000 10:44 Control Test 0.081 10:45 Air Blank 0.000 10:45 Control Test 0.080 10:46 Air Blank 0.000 10:46 Control Test Stats Average 0.0807 Std Dev 0.0006 Rel Std Dev(%) 0.7157 <i>M. J. Jester</i> Operator's Signature

Verified instrument was able to read gas pressure.



PUTNAM COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001129
06/23/2026 14:02:09

#1

Auto Calibration

SW

pg 1 of 2

<<<< 3um >>>> <<<< 9um >>>>

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.1150	(-0.0220)	0.1130	(-0.0150)
Sample #2	0.1390	(-0.0230)	0.1270	(-0.0240)
Sample #3	0.1010	(-0.0160)	0.1460	(-0.0580)
Sample #4	0.0970	(-0.0040)	0.1300	(-0.0520)
Avg % Abs	0.1123	(-0.0143)	0.1343	(-0.0447)
STD DEV	0.0232	(0.0096)	0.0102	(0.0181)
REL STD DEV	20.635	(67.040)	7.604	(40.629)

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.7990	(-0.0150)	1.5160	(-0.0120)
Sample #2	0.7950	(0.0120)	1.5180	(-0.0110)
Sample #3	0.7820	(0.0200)	1.5160	(-0.0120)
Sample #4	0.7810	(-0.0050)	1.5050	(-0.0120)
Avg % Abs	0.7860	(0.0090)	1.5130	(-0.0117)
STD DEV	0.0078	(0.0128)	0.0070	(0.0006)
REL STD DEV	0.994	(141.857)	0.463	(4.949)

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	1.8330	(-0.0080)	3.5610	(0.0000)
Sample #2	1.8310	(-0.0140)	3.5770	(-0.0010)
Sample #3	1.8150	(0.0100)	3.5530	(0.0030)
Sample #4	1.8150	(0.0100)	3.5710	(0.0000)
Avg % Abs	1.8203	(0.0020)	3.5670	(0.0007)
STD DEV	0.0092	(0.0139)	0.0125	(0.0021)
REL STD DEV	0.507	(692.820)	0.350	(312.250)

Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	3.4700	(-0.0120)	6.8040	(-0.0110)
Sample #2	3.5320	(-0.0290)	6.8440	(-0.0110)
Sample #3	3.4990	(-0.0270)	6.7980	(-0.0210)
Sample #4	3.4960	(-0.0300)	6.8030	(-0.0260)
Avg % Abs	3.5090	(-0.0287)	6.8150	(-0.0193)
STD DEV	0.0200	(0.0015)	0.0252	(0.0076)
REL STD DEV	0.569	(5.329)	0.370	(39.505)

Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	5.0700	(-0.0180)	9.7940	(0.0000)
Sample #2	5.0800	(-0.0020)	9.8160	(-0.0100)
Sample #3	5.0910	(-0.0290)	9.8020	(-0.0310)
Sample #4	5.0610	(-0.0070)	9.7810	(-0.0100)
Avg % Abs	5.0773	(-0.0127)	9.7997	(-0.0170)
STD DEV	0.0152	(0.0144)	0.0176	(0.0121)
REL STD DEV	0.299	(113.402)	0.180	(71.320)

PUTNAM COUNTY SO
 Intoxilyzer - Alcohol Analyzer
 Model 8000 SN 80-001129
 06/23/2026 14:02:09

Auto Calibration

pg 2 of 2

```

<<<<<      3um      >>>>>
-----
Zero Order Coef   -262.92
First Order Coef  2677.57
Second Order Coef  36.17
  
```

```

<<<<<      9um      >>>>>
-----
Zero Order Coef   -152.73
First Order Coef  1316.96
Second Order Coef  15.82
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L) (g/210L)
0.000    0.001    -0.0008
0.040    0.039    0.0009
0.100    0.099    0.0006
0.200    0.201    -0.0011
0.300    0.300    0.0004
  
```

```

-----
Act      Fit      Residual
(g/210L) (g/210L) (g/210L)
0.000    0.001    -0.0005
0.040    0.039    0.0006
0.100    0.100    0.0003
0.200    0.201    -0.0007
0.300    0.300    0.0003
  
```

```

<<<<<      3um      >>>>>
-----
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
Sample
Sample #1      3144.00      3308.00
Sample #2      3148.00      3335.00
Sample #3      3154.00      3284.00
Sample #4      3143.00      3278.00
Avg            3148.3333    3299.0000
STD DEV        5.5076      31.3209
REL STD DEV    0.175      0.949
H2O adjust (mg/l*10k) 661      510
  
```

Barometric Pressure = 1012

*****CALIBRATION SUCCESSFUL*****

SW

Post Optical Bench Adjustment Stability Checks 80-001129 Sub

0.050 g/210L 0.047 to 0.053 g/210L	0.080 g/210L 0.077 to 0.083 g/210L	0.200 g/210L 0.194 to 0.206 g/210L	DGS 0.080 g/210L 0.077 to 0.083 g/210L 50.003 g/210L of Wet																																																																																																																																																												
☒ PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8100 SN 80-001129 06/23/2026 Software: 8100.27 <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>15:17</td></tr> <tr><td>Control Test</td><td>0.049</td><td>15:17</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:18</td></tr> <tr><td>Control Test</td><td>0.049</td><td>15:18</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:19</td></tr> <tr><td>Control Test</td><td>0.049</td><td>15:20</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:20</td></tr> <tr><td>Control Test</td><td>0.000</td><td>15:20</td></tr> <tr><td colspan="3">Control Test Stats</td></tr> <tr><td>Average</td><td>0.0490</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> Operator Signature: <i>M. J. Miller</i>	Test	g/210L	Time	Air Blank	0.000	15:17	Control Test	0.049	15:17	Air Blank	0.000	15:18	Control Test	0.049	15:18	Air Blank	0.000	15:19	Control Test	0.049	15:20	Air Blank	0.000	15:20	Control Test	0.000	15:20	Control Test Stats			Average	0.0490		Std Dev	0.0000		Rel Std Dev(%)	0.0000		☒ PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8100 SN 80-001129 06/23/2026 Software: 8100.27 <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>15:06</td></tr> <tr><td>Control Test</td><td>0.078</td><td>15:06</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:06</td></tr> <tr><td>Control Test</td><td>0.079</td><td>15:06</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:10</td></tr> <tr><td>Control Test</td><td>0.077</td><td>15:10</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:11</td></tr> <tr><td>Control Test</td><td>0.000</td><td>15:12</td></tr> <tr><td colspan="3">Control Test Stats</td></tr> <tr><td>Average</td><td>0.0780</td><td></td></tr> <tr><td>Std Dev</td><td>0.0010</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.2821</td><td></td></tr> </tbody> </table> Operator Signature: <i>M. J. Miller</i>	Test	g/210L	Time	Air Blank	0.000	15:06	Control Test	0.078	15:06	Air Blank	0.000	15:06	Control Test	0.079	15:06	Air Blank	0.000	15:10	Control Test	0.077	15:10	Air Blank	0.000	15:11	Control Test	0.000	15:12	Control Test Stats			Average	0.0780		Std Dev	0.0010		Rel Std Dev(%)	1.2821		☒ PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8100 SN 80-001129 06/23/2026 Software: 8100.27 <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>15:27</td></tr> <tr><td>Control Test</td><td>0.197</td><td>15:28</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:28</td></tr> <tr><td>Control Test</td><td>0.197</td><td>15:29</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:29</td></tr> <tr><td>Control Test</td><td>0.198</td><td>15:30</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:31</td></tr> <tr><td>Control Test</td><td>0.000</td><td>15:31</td></tr> <tr><td colspan="3">Control Test Stats</td></tr> <tr><td>Average</td><td>0.1973</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2926</td><td></td></tr> </tbody> </table> Operator Signature: <i>M. J. Miller</i>	Test	g/210L	Time	Air Blank	0.000	15:27	Control Test	0.197	15:28	Air Blank	0.000	15:28	Control Test	0.197	15:29	Air Blank	0.000	15:29	Control Test	0.198	15:30	Air Blank	0.000	15:31	Control Test	0.000	15:31	Control Test Stats			Average	0.1973		Std Dev	0.0006		Rel Std Dev(%)	0.2926		☒ PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8100 SN 80-001129 06/23/2026 Software: 8100.27 <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>15:22</td></tr> <tr><td>Control Test</td><td>0.081</td><td>15:23</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:23</td></tr> <tr><td>Control Test</td><td>0.082</td><td>15:24</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:24</td></tr> <tr><td>Control Test</td><td>0.081</td><td>15:24</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>15:25</td></tr> <tr><td>Control Test</td><td>0.000</td><td>15:25</td></tr> <tr><td colspan="3">Control Test Stats</td></tr> <tr><td>Average</td><td>0.0813</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.7099</td><td></td></tr> </tbody> </table> Operator Signature: <i>M. J. Miller</i>	Test	g/210L	Time	Air Blank	0.000	15:22	Control Test	0.081	15:23	Air Blank	0.000	15:23	Control Test	0.082	15:24	Air Blank	0.000	15:24	Control Test	0.081	15:24	Air Blank	0.000	15:25	Control Test	0.000	15:25	Control Test Stats			Average	0.0813		Std Dev	0.0006		Rel Std Dev(%)	0.7099	
Test	g/210L	Time																																																																																																																																																													
Air Blank	0.000	15:17																																																																																																																																																													
Control Test	0.049	15:17																																																																																																																																																													
Air Blank	0.000	15:18																																																																																																																																																													
Control Test	0.049	15:18																																																																																																																																																													
Air Blank	0.000	15:19																																																																																																																																																													
Control Test	0.049	15:20																																																																																																																																																													
Air Blank	0.000	15:20																																																																																																																																																													
Control Test	0.000	15:20																																																																																																																																																													
Control Test Stats																																																																																																																																																															
Average	0.0490																																																																																																																																																														
Std Dev	0.0000																																																																																																																																																														
Rel Std Dev(%)	0.0000																																																																																																																																																														
Test	g/210L	Time																																																																																																																																																													
Air Blank	0.000	15:06																																																																																																																																																													
Control Test	0.078	15:06																																																																																																																																																													
Air Blank	0.000	15:06																																																																																																																																																													
Control Test	0.079	15:06																																																																																																																																																													
Air Blank	0.000	15:10																																																																																																																																																													
Control Test	0.077	15:10																																																																																																																																																													
Air Blank	0.000	15:11																																																																																																																																																													
Control Test	0.000	15:12																																																																																																																																																													
Control Test Stats																																																																																																																																																															
Average	0.0780																																																																																																																																																														
Std Dev	0.0010																																																																																																																																																														
Rel Std Dev(%)	1.2821																																																																																																																																																														
Test	g/210L	Time																																																																																																																																																													
Air Blank	0.000	15:27																																																																																																																																																													
Control Test	0.197	15:28																																																																																																																																																													
Air Blank	0.000	15:28																																																																																																																																																													
Control Test	0.197	15:29																																																																																																																																																													
Air Blank	0.000	15:29																																																																																																																																																													
Control Test	0.198	15:30																																																																																																																																																													
Air Blank	0.000	15:31																																																																																																																																																													
Control Test	0.000	15:31																																																																																																																																																													
Control Test Stats																																																																																																																																																															
Average	0.1973																																																																																																																																																														
Std Dev	0.0006																																																																																																																																																														
Rel Std Dev(%)	0.2926																																																																																																																																																														
Test	g/210L	Time																																																																																																																																																													
Air Blank	0.000	15:22																																																																																																																																																													
Control Test	0.081	15:23																																																																																																																																																													
Air Blank	0.000	15:23																																																																																																																																																													
Control Test	0.082	15:24																																																																																																																																																													
Air Blank	0.000	15:24																																																																																																																																																													
Control Test	0.081	15:24																																																																																																																																																													
Air Blank	0.000	15:25																																																																																																																																																													
Control Test	0.000	15:25																																																																																																																																																													
Control Test Stats																																																																																																																																																															
Average	0.0813																																																																																																																																																														
Std Dev	0.0006																																																																																																																																																														
Rel Std Dev(%)	0.7099																																																																																																																																																														

* Repeat OBA Sub

PUTNAM COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001129
06/23/2026 15:43:05

#2

SW

Auto Calibration

pg 1 of 2

<<<< 3um >>>> <<<< 9um >>>>

Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.0990	(-0.0220)	0.1460	(-0.0060)
Sample #2	0.0930	(0.0070)	0.1470	(-0.0030)
Sample #3	0.0990	(0.0090)	0.1290	(-0.0040)
Sample #4	0.1050	(-0.0120)	0.1490	(-0.0120)
Avg % Abs	0.0990	(0.0013)	0.1417	(-0.0063)
STD DEV	0.0060	(0.0116)	0.0110	(0.0049)
REL STD DEV	6.061	(869.267)	7.775	(77.888)

Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	0.8000	(-0.0240)	1.5270	(-0.0140)
Sample #2	0.7590	(0.0010)	1.5220	(-0.0070)
Sample #3	0.7910	(-0.0230)	1.5340	(-0.0280)
Sample #4	0.7850	(-0.0020)	1.5290	(-0.0110)
Avg % Abs	0.7783	(-0.0080)	1.5283	(-0.0153)
STD DEV	0.0170	(0.0131)	0.0060	(0.0112)
REL STD DEV	2.185	(163.459)	0.394	(72.721)

Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	1.7930	(-0.0040)	3.5360	(-0.0190)
Sample #2	1.8150	(0.0010)	3.5430	(-0.0130)
Sample #3	1.8190	(-0.0190)	3.5470	(-0.0250)
Sample #4	1.8130	(-0.0130)	3.5460	(-0.0190)
Avg % Abs	1.8157	(-0.0103)	3.5453	(-0.0190)
STD DEV	0.0031	(0.0103)	0.0021	(0.0060)
REL STD DEV	0.168	(99.321)	0.059	(31.579)

Solution = 0.200 g/210L or 0.9524 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	3.4830	(-0.0140)	6.7790	(-0.0200)
Sample #2	3.4880	(-0.0040)	6.7960	(0.0050)
Sample #3	3.5010	(-0.0090)	6.8030	(0.0130)
Sample #4	3.4960	(0.0080)	6.7880	(0.0200)
Avg % Abs	3.4950	(-0.0017)	6.7957	(0.0127)
STD DEV	0.0066	(0.0087)	0.0075	(0.0075)
REL STD DEV	0.188	(524.214)	0.110	(59.254)

Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1

Sample	% Abs	(% Abs Ref)	% Abs	(% Abs Ref)
Sample #1	5.1260	(-0.0200)	9.8350	(-0.0070)
Sample #2	5.0920	(0.0140)	9.8410	(0.0040)
Sample #3	5.1070	(0.0220)	9.8350	(0.0170)
Sample #4	5.0660	(0.0320)	9.8270	(0.0210)
Avg % Abs	5.0883	(0.0227)	9.8343	(0.0140)
STD DEV	0.0207	(0.0090)	0.0070	(0.0089)
REL STD DEV	0.408	(39.787)	0.071	(63.487)

PUTNAM COUNTY SO
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-001129
06/23/2026 15:43:05

Auto Calibration

Sut

pg 2 of 2

<<<< 3um >>>>

Zero Order Coef -241.18
First Order Coef 2691.78
Second Order Coef 31.61

<<<< 9um >>>>

-182.76
1344.09
12.84

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.001	-0.0005
0.040	0.039	0.0007
0.100	0.100	0.0002
0.200	0.201	-0.0006
0.300	0.300	0.0002

Act (g/210L)	Fit (g/210L)	Residual (g/210L)
0.000	0.000	-0.0002
0.040	0.040	0.0001
0.100	0.100	0.0004
0.200	0.200	-0.0004
0.300	0.300	0.0002

<<<< 3um >>>> <<<< 9um >>>>

Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
Sample
Sample #1 3182.00 3352.00
Sample #2 3236.00 3355.00
Sample #3 3180.00 3350.00
Sample #4 3252.00 3339.00
Avg 3222.6667 3348.0000
STD DEV 37.8065 8.1854
REL STD DEV 1.173 0.244
H2O adjust (mg/l*10k) 587 461

Barometric Pressure = 1012

*****CALIBRATION SUCCESSFUL*****

#1

PULHAM COUNTY SO
Intoxilyzer - A Alcohol Analyzer
Model 8000
06/23/2026
SN 80-001129
Software: 81C0.27

Test	g/210L	Time
Air Blank	0.000	17:26
Control Test	0.078	17:27
Air Blank	0.000	17:27
Control Test	0.078	17:28
Air Blank	0.000	17:29
Control Test	0.078	17:29
Air Blank	0.000	17:30

Allyla

Internal
print paper
depleted. Re-run
SST

Post Optical Bench Adjustment Stability Checks #2

0.050 g/210L	0.080 g/210L	0.200 g/210L	DGS 0.080 g/210L
0.047 to 0.053 g/210L	0.077 to 0.083 g/210L	0.194 to 0.206 g/210L	0.077 to 0.083 g/210L 50.003 g/210L of Wet
PUNNAM COUNTY SO Intoxilyzer - Alcotest Analyzer Model 8000 06/23/2026 SN 80-001129 Software: 8100.27 Test g/210L Time Air Blank 0.000 17:19 Control Test 0.050 17:20 Air Blank 0.000 17:20 Control Test 0.049 17:21 Air Blank 0.000 17:21 Control Test 0.049 17:22 Air Blank 0.000 17:22 Control Test Stats Average 0.0493 Std Dev 0.0006 Rel Std Dev(%) 1.1703 Operator's Signature <i>M. Sigal</i>	PUNNAM COUNTY SO Intoxilyzer - Alcotest Analyzer Model 8000 06/23/2026 SN 80-001129 Software: 8100.27 Test g/210L Time Air Blank 0.000 17:33 Control Test 0.078 17:34 Air Blank 0.000 17:35 Control Test 0.078 17:35 Air Blank 0.000 17:35 Control Test 0.078 17:35 Air Blank 0.000 17:35 Control Test 0.078 17:35 Air Blank 0.000 17:37 Control Test Stats Average 0.0780 Std Dev 0.0008 Rel Std Dev(%) 0.0000 Operator's Signature <i>M. Sigal</i>	PUNNAM COUNTY SO Intoxilyzer - Alcotest Analyzer Model 8000 06/23/2026 SN 80-001129 Software: 8100.27 Test g/210L Time Air Blank 0.000 17:33 Control Test 0.158 17:34 Air Blank 0.000 17:35 Control Test 0.157 17:35 Air Blank 0.000 17:35 Control Test 0.157 17:35 Air Blank 0.000 17:35 Control Test 0.157 17:35 Air Blank 0.000 17:35 Control Test Stats Average 0.1577 Std Dev 0.0006 Rel Std Dev(%) 0.292 Operator's Signature <i>M. Sigal</i>	PUNNAM COUNTY SO Intoxilyzer - Alcotest Analyzer Model 8000 06/23/2026 SN 80-001129 Software: 8100.27 Test g/210L Time Air Blank 0.000 17:43 Control Test 0.080 17:43 Air Blank 0.000 17:44 Control Test 0.080 17:44 Air Blank 0.000 17:44 Control Test 0.079 17:44 Air Blank 0.000 17:45 Control Test Stats Average 0.0797 Std Dev 0.0006 Rel Std Dev(%) 0.7247 Operator's Signature <i>M. Sigal</i>

80-001129

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: PUTNAM COUNTY SO
Time of Inspection: 12:05

Date of Inspection: 06/25/2026

Serial Number: 80-001129
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23825080A1 Exp: 10/05/2027
0.000	0.050	0.078	0.195	0.080
0.000	0.050	0.078	0.195	0.080
0.000	0.050	0.079	0.194	0.080
0.000	0.050	0.078	0.194	0.080
0.000	0.049	0.078	0.193	0.081
0.000	0.049	0.079	0.193	0.081
0.000	0.050	0.079	0.196	0.081
0.000	0.050	0.079	0.198	0.081
0.000	0.050	0.078	0.198	0.080
0.000	0.050	0.078	0.198	0.081


Standard Deviations	0.0004	0.0005	0.0020	0.0005
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0008 Number of Simulators Used: 5

Remarks:

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.


I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.



Signature and Printed Name

LEANDRA HIGGINBOTHAM

06/25/2026
Date



ANSI National Accreditation Board
ACCREDITED
ISO/IEC 17025
FORENSIC CALIBRATION
LABORATORY

Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-001129, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-001129</u>	UNCERTAINTY* $\pm$	
Owning Agency:	<u>PUTNAM COUNTY SO</u>	0.050 g/210 L	0.004
Calibration Date:	<u>06/25/2026</u>	0.080 g/210 L	0.004
Calibration Time:	<u>12:05</u>	0.200 g/210 L	0.008
		0.080 g/210 L Dry Gas Control	0.005

All results are reported in g/210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration. *Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence ($k=3$).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/IEC 17025 standards. This document shall not be reproduced except in full without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

06/25/2026
Date

LEANDRA HIGGINBOTHAM,
Department Inspector

Alcohol Testing Program - Instrument Processing Sheet

Agency: PUTNAM COUNTY SO

Instrument Serial Number: 80-001129

Date In: 1/6/2026

DI Completion Date: N/A

☐ Ship

☐ P/U

☐ H/D

☒ CMI

☐ EE

Intake By: <u>KTS</u> Date: <u>1/6/2026</u>		Quality Checks By: <u>SLH</u> Date: <u>1/13/2026</u>		Flow Adjustment By: _____																																																																																								
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☐ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>93 *</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP 103 *</u> 32 mm <u>0.144</u> (<u>.139-.169</u>) 36 mm <u>0.160</u> (<u>.156-.190</u>) 53 mm <u>0.234</u> (<u>.228-.278</u>) 103 mm <u>0.488</u> (<u>.447-.547</u>) ☒ Barometric Pressure Check Gauge ID #: <u>28427</u> Gauge: <u>1018</u> Instrument: <u>1019</u> ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm (<u>.139-.169</u>) 36 mm (<u>.156-.190</u>) 53 mm (<u>.228-.278</u>) 103 mm (<u>.447-.547</u>)																																																																																								
						<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 20%;">Serial #</th> <th style="width: 25%;">Lot#/Exp</th> <th style="width: 40%;">Maintenance By: _____ Date: _____</th> </tr> </thead> <tbody> <tr> <td>0.050</td> <td rowspan="2">MP5088</td> <td>202406K</td> <td>☐ Battery Replacement</td> </tr> <tr> <td></td> <td>06/19/2026</td> <td>☐ Dry Gas Regulator Replacement</td> </tr> <tr> <td>0.080</td> <td rowspan="2">MP5089</td> <td>202406L</td> <td>☐ Tank Sensor Tare</td> </tr> <tr> <td></td> <td>06/19/2026</td> <td>☐ Breath Tube Replacement</td> </tr> <tr> <td>0.200</td> <td rowspan="2">MP5090</td> <td>202406N</td> <td>☐ Other:</td> </tr> <tr> <td></td> <td>06/20/2026</td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td rowspan="2">N/A</td> <td>AG510701</td> <td></td> </tr> <tr> <td></td> <td>04/17/2027</td> <td></td> </tr> </tbody> </table>		Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____	0.050	MP5088	202406K	☐ Battery Replacement		06/19/2026	☐ Dry Gas Regulator Replacement	0.080	MP5089	202406L	☐ Tank Sensor Tare		06/19/2026	☐ Breath Tube Replacement	0.200	MP5090	202406N	☐ Other:		06/20/2026		0.080 DGS	N/A	AG510701			04/17/2027																																																						
						Simulator	Serial #	Lot#/Exp	Maintenance By: _____ Date: _____																																																																																			
						0.050	MP5088	202406K	☐ Battery Replacement																																																																																			
								06/19/2026	☐ Dry Gas Regulator Replacement																																																																																			
0.080	MP5089	202406L	☐ Tank Sensor Tare																																																																																									
		06/19/2026	☐ Breath Tube Replacement																																																																																									
0.200	MP5090	202406N	☐ Other:																																																																																									
		06/20/2026																																																																																										
0.080 DGS	N/A	AG510701																																																																																										
		04/17/2027																																																																																										
Optical Bench Adjustment By: _____ Department Inspection By: _____ Barometric Pressure Gauge: _____ ID#: _____ Barometric Pressure ID#: _____																																																																																												
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 15%;">Simulator</th> <th style="width: 15%;">Serial #</th> <th style="width: 15%;">Lot #</th> <th style="width: 15%;">Expiration</th> <th style="width: 40%;">Gauge: _____ Instrument: _____</th> </tr> </thead> <tbody> <tr> <td>0.000</td> <td></td> <td>N/A</td> <td>N/A</td> <td>Mouth Alcohol Solution Lot #: _____ Exp: _____</td> </tr> <tr> <td>0.040</td> <td></td> <td></td> <td></td> <td>Acetone Stock Solution Lot #: _____ Exp: _____</td> </tr> <tr> <td>0.100</td> <td></td> <td></td> <td></td> <td>Simulator Serial Number</td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> <td>0.000</td> </tr> <tr> <td>0.300</td> <td></td> <td></td> <td></td> <td>Interferent</td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> <td>0.050</td> </tr> <tr> <td colspan="4">☐ Post Optical Bench Adjustment Stability Checks</td> <td>0.080</td> </tr> <tr> <td>Simulator</td> <td>Serial #</td> <td>Lot #</td> <td>Expiration</td> <td>0.200</td> </tr> <tr> <td>0.050</td> <td></td> <td></td> <td></td> <td rowspan="4" style="background-color: #d3d3d3;"> Attachments ☐ Form 41 ☐ Post-Stability Checks ☒ Stability Checks ☐ Flow Adjustment ☐ Calibration Certificate ☐ Form 40 ☐ Optical Bench Adjustment ☒ Other: RMA sheet </td> </tr> <tr> <td>0.080</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.200</td> <td></td> <td></td> <td></td> </tr> <tr> <td>0.080 DGS</td> <td>N/A</td> <td></td> <td></td> </tr> <tr> <td colspan="4">Gauge ID #: _____</td> <td rowspan="4" style="background-color: #d3d3d3;"> ☐ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use </td> </tr> <tr> <td colspan="4">Gauge: _____ Instrument: _____</td> </tr> <tr> <td colspan="4">Notes/Suggested Service: * Flow portion of quality checks performed 1/20/26 and is when noted R-value. **SLH RMA = return material authorization. SLH 1/28/2026 **Tech Review-added initials to 1st note on 1/29/26 SLH. </td> </tr> <tr> <td colspan="4"></td> </tr> <tr> <td colspan="4"></td> <td colspan="2" style="background-color: #d3d3d3;"> <table style="width: 100%;"> <tr> <td style="width: 50%; text-align: center;"> Digitally signed by Taylor Gutschow Date: 2026.01.29 14:13:13 -05'00' </td> <td style="width: 50%; text-align: center;"> Digitally signed by Wen-Chi Pierson Date: 2026.01.29 15:11:22 -05'00' </td> </tr> <tr> <td style="text-align: center;">Taylor Gutschow</td> <td style="text-align: center;">Wen-Chi Pierson</td> </tr> <tr> <td style="text-align: center;">Tech Review</td> <td style="text-align: center;">Admin Review</td> </tr> </table> </td> </tr> </tbody></table>		Simulator	Serial #	Lot #	Expiration	Gauge: _____ Instrument: _____	0.000		N/A	N/A	Mouth Alcohol Solution Lot #: _____ Exp: _____	0.040				Acetone Stock Solution Lot #: _____ Exp: _____	0.100				Simulator Serial Number	0.200				0.000	0.300				Interferent	0.080 DGS	N/A			0.050	☐ Post Optical Bench Adjustment Stability Checks				0.080	Simulator	Serial #	Lot #	Expiration	0.200	0.050				Attachments ☐ Form 41 ☐ Post-Stability Checks ☒ Stability Checks ☐ Flow Adjustment ☐ Calibration Certificate ☐ Form 40 ☐ Optical Bench Adjustment ☒ Other: RMA sheet	0.080				0.200				0.080 DGS	N/A			Gauge ID #: _____				☐ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use	Gauge: _____ Instrument: _____				Notes/Suggested Service: * Flow portion of quality checks performed 1/20/26 and is when noted R-value. **SLH RMA = return material authorization. SLH 1/28/2026 **Tech Review-added initials to 1st note on 1/29/26 SLH.												<table style="width: 100%;"> <tr> <td style="width: 50%; text-align: center;"> Digitally signed by Taylor Gutschow Date: 2026.01.29 14:13:13 -05'00' </td> <td style="width: 50%; text-align: center;"> Digitally signed by Wen-Chi Pierson Date: 2026.01.29 15:11:22 -05'00' </td> </tr> <tr> <td style="text-align: center;">Taylor Gutschow</td> <td style="text-align: center;">Wen-Chi Pierson</td> </tr> <tr> <td style="text-align: center;">Tech Review</td> <td style="text-align: center;">Admin Review</td> </tr> </table>		Digitally signed by Taylor Gutschow Date: 2026.01.29 14:13:13 -05'00'	Digitally signed by Wen-Chi Pierson Date: 2026.01.29 15:11:22 -05'00'	Taylor Gutschow	Wen-Chi Pierson	Tech Review	Admin Review
Simulator	Serial #	Lot #	Expiration	Gauge: _____ Instrument: _____																																																																																								
0.000		N/A	N/A	Mouth Alcohol Solution Lot #: _____ Exp: _____																																																																																								
0.040				Acetone Stock Solution Lot #: _____ Exp: _____																																																																																								
0.100				Simulator Serial Number																																																																																								
0.200				0.000																																																																																								
0.300				Interferent																																																																																								
0.080 DGS	N/A			0.050																																																																																								
☐ Post Optical Bench Adjustment Stability Checks				0.080																																																																																								
Simulator	Serial #	Lot #	Expiration	0.200																																																																																								
0.050				Attachments ☐ Form 41 ☐ Post-Stability Checks ☒ Stability Checks ☐ Flow Adjustment ☐ Calibration Certificate ☐ Form 40 ☐ Optical Bench Adjustment ☒ Other: RMA sheet																																																																																								
0.080																																																																																												
0.200																																																																																												
0.080 DGS	N/A																																																																																											
Gauge ID #: _____				☐ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☐ Return to/Place into Evidentiary Use ☒ Remain Out of Evidentiary Use ☐ Conduct an Agency Inspection Before Evidentiary Use																																																																																								
Gauge: _____ Instrument: _____																																																																																												
Notes/Suggested Service: * Flow portion of quality checks performed 1/20/26 and is when noted R-value. **SLH RMA = return material authorization. SLH 1/28/2026 **Tech Review-added initials to 1st note on 1/29/26 SLH.																																																																																												
				<table style="width: 100%;"> <tr> <td style="width: 50%; text-align: center;"> Digitally signed by Taylor Gutschow Date: 2026.01.29 14:13:13 -05'00' </td> <td style="width: 50%; text-align: center;"> Digitally signed by Wen-Chi Pierson Date: 2026.01.29 15:11:22 -05'00' </td> </tr> <tr> <td style="text-align: center;">Taylor Gutschow</td> <td style="text-align: center;">Wen-Chi Pierson</td> </tr> <tr> <td style="text-align: center;">Tech Review</td> <td style="text-align: center;">Admin Review</td> </tr> </table>		Digitally signed by Taylor Gutschow Date: 2026.01.29 14:13:13 -05'00'	Digitally signed by Wen-Chi Pierson Date: 2026.01.29 15:11:22 -05'00'	Taylor Gutschow	Wen-Chi Pierson	Tech Review	Admin Review																																																																																	
Digitally signed by Taylor Gutschow Date: 2026.01.29 14:13:13 -05'00'	Digitally signed by Wen-Chi Pierson Date: 2026.01.29 15:11:22 -05'00'																																																																																											
Taylor Gutschow	Wen-Chi Pierson																																																																																											
Tech Review	Admin Review																																																																																											

Stability Checks

0.050 g/210L 0.047 to 0.053 g/210L	0.080 g/210L 0.077 to 0.083 g/210L	0.200 g/210L 0.194 to 0.206 g/210L	DGS 0.080 g/210L 0.077 to 0.083 g/210L ≤0.003 g/210L of Wet																																																																																																																																																
✓ PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001129 01/13/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:50</td></tr> <tr><td>Control Test</td><td>0.050</td><td>12:51</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:51</td></tr> <tr><td>Control Test</td><td>0.049</td><td>12:52</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:53</td></tr> <tr><td>Control Test</td><td>0.049</td><td>12:53</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:54</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0493</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>1.1713</td><td></td></tr> </tbody> </table> <i>Shiglot</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:50	Control Test	0.050	12:51	Air Blank	0.000	12:51	Control Test	0.049	12:52	Air Blank	0.000	12:53	Control Test	0.049	12:53	Air Blank	0.000	12:54	Control Test Stats			Average	0.0493		Std Dev	0.0006		Rel Std Dev(%)	1.1713		✓ PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001129 01/13/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:31</td></tr> <tr><td>Control Test</td><td>0.080</td><td>12:31</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:32</td></tr> <tr><td>Control Test</td><td>0.080</td><td>12:33</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:33</td></tr> <tr><td>Control Test</td><td>0.080</td><td>12:34</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:34</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0800</td><td></td></tr> <tr><td>Std Dev</td><td>0.0000</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.0000</td><td></td></tr> </tbody> </table> <i>Shiglot</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:31	Control Test	0.080	12:31	Air Blank	0.000	12:32	Control Test	0.080	12:33	Air Blank	0.000	12:33	Control Test	0.080	12:34	Air Blank	0.000	12:34	Control Test Stats			Average	0.0800		Std Dev	0.0000		Rel Std Dev(%)	0.0000		✓ PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001129 01/13/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:43</td></tr> <tr><td>Control Test</td><td>0.195</td><td>12:43</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:44</td></tr> <tr><td>Control Test</td><td>0.195</td><td>12:44</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:45</td></tr> <tr><td>Control Test</td><td>0.196</td><td>12:46</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:46</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.1953</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.2956</td><td></td></tr> </tbody> </table> <i>Shiglot</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:43	Control Test	0.195	12:43	Air Blank	0.000	12:44	Control Test	0.195	12:44	Air Blank	0.000	12:45	Control Test	0.196	12:46	Air Blank	0.000	12:46	Control Test Stats			Average	0.1953		Std Dev	0.0006		Rel Std Dev(%)	0.2956		DGS PUTNAM COUNTY SO Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-001129 01/13/2026 Software: 8100.27 <table border="1"> <thead> <tr> <th>Test</th> <th>g/210L</th> <th>Time</th> </tr> </thead> <tbody> <tr><td>Air Blank</td><td>0.000</td><td>12:38</td></tr> <tr><td>Control Test</td><td>0.083</td><td>12:38</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:39</td></tr> <tr><td>Control Test</td><td>0.082</td><td>12:39</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:40</td></tr> <tr><td>Control Test</td><td>0.083</td><td>12:40</td></tr> <tr><td>Air Blank</td><td>0.000</td><td>12:40</td></tr> <tr><td>Control Test Stats</td><td></td><td></td></tr> <tr><td>Average</td><td>0.0827</td><td></td></tr> <tr><td>Std Dev</td><td>0.0006</td><td></td></tr> <tr><td>Rel Std Dev(%)</td><td>0.6984</td><td></td></tr> </tbody> </table> <i>Shiglot</i> Operator's Signature	Test	g/210L	Time	Air Blank	0.000	12:38	Control Test	0.083	12:38	Air Blank	0.000	12:39	Control Test	0.082	12:39	Air Blank	0.000	12:40	Control Test	0.083	12:40	Air Blank	0.000	12:40	Control Test Stats			Average	0.0827		Std Dev	0.0006		Rel Std Dev(%)	0.6984	
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:50																																																																																																																																																	
Control Test	0.050	12:51																																																																																																																																																	
Air Blank	0.000	12:51																																																																																																																																																	
Control Test	0.049	12:52																																																																																																																																																	
Air Blank	0.000	12:53																																																																																																																																																	
Control Test	0.049	12:53																																																																																																																																																	
Air Blank	0.000	12:54																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0493																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	1.1713																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:31																																																																																																																																																	
Control Test	0.080	12:31																																																																																																																																																	
Air Blank	0.000	12:32																																																																																																																																																	
Control Test	0.080	12:33																																																																																																																																																	
Air Blank	0.000	12:33																																																																																																																																																	
Control Test	0.080	12:34																																																																																																																																																	
Air Blank	0.000	12:34																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0800																																																																																																																																																		
Std Dev	0.0000																																																																																																																																																		
Rel Std Dev(%)	0.0000																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:43																																																																																																																																																	
Control Test	0.195	12:43																																																																																																																																																	
Air Blank	0.000	12:44																																																																																																																																																	
Control Test	0.195	12:44																																																																																																																																																	
Air Blank	0.000	12:45																																																																																																																																																	
Control Test	0.196	12:46																																																																																																																																																	
Air Blank	0.000	12:46																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.1953																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.2956																																																																																																																																																		
Test	g/210L	Time																																																																																																																																																	
Air Blank	0.000	12:38																																																																																																																																																	
Control Test	0.083	12:38																																																																																																																																																	
Air Blank	0.000	12:39																																																																																																																																																	
Control Test	0.082	12:39																																																																																																																																																	
Air Blank	0.000	12:40																																																																																																																																																	
Control Test	0.083	12:40																																																																																																																																																	
Air Blank	0.000	12:40																																																																																																																																																	
Control Test Stats																																																																																																																																																			
Average	0.0827																																																																																																																																																		
Std Dev	0.0006																																																																																																																																																		
Rel Std Dev(%)	0.6984																																																																																																																																																		

Return Material Authorization

Ship to:☒ CMI, Inc.☐ Enforcement ElectronicsShipment to repair facility authorized by: Trevor Bishop on 1/20/2026Items Returned: ☒ Instrument ☐ Supplies ☐ Other Describe: _____Instrument Model: Intoxilyzer 8000 Serial Number: 80-001129

Bill To Address:

Putnam County SO
130 Orié Griffin Blvd
Palatka, FL

Ship to Address:

FDLE Offsite Mail Facility
c/o FDLE HQ
813 B Lake Bradford Rd
Tallahassee, FL 32304

Reason for Return:

Replace flow sensor

I require an estimate **BEFORE** any repairs will be authorized and/or conductedPlease contact: Trevor Bishop Phone #: 386-240-8591Email: tbishop@putnamsheriff.orgATP Contact Name: LeAndra Higginbotham ATP Email: leandrahigginbotham@fdle.state.fl.us