

Alcohol Testing Program - Instrument Processing Sheet

Agency: SUMTER COUNTY SO

Instrument Serial Number: 80-000815

Date In: 8/12/2026

DI Completion Date: 8/21/26

☒ Ship ☐ P/U ☐ H/D ☐ CMI ☐ EE

Intake By: <u>KTS</u> Date: <u>8/12/2026</u>		Quality Checks By: <u>IAS</u> Date: <u>8/19/2026</u>		Flow Adjustment By: _____																																						
☒ Annual ☐ Dropped Off ☐ Registration ☐ Return from CMI / EE ☐ Training Instrument Visual Inspection ☒ Case ☒ Handle ☒ Keyboard ☒ Dry Gas Shelf ☒ Feet ☒ Breath Tube ☒ Ports ☒ Screws Tight Other Equipment/Accessories ☐ Power Cord ☐ Printer Cable ☒ Static Bag ☐ 12V DC Cable Notes:		☒ Breath Tube Screen ☒ Replace External O-Rings ☒ Instrument Set Up Verified ☒ R-Value: <u>190</u> ☒ Flow Verification (L/s) Flow Column #: <u>ATP103</u> 32 mm <u>.156</u> (.139-.169) 36 mm <u>.171</u> (.156-.190) 53 mm <u>.234</u> (.228-.278) 103 mm <u>.503</u> (.447-.547) ☒ Barometric Pressure Check Gauge ID #: <u>28427</u> Gauge: <u>1011</u> Instrument: <u>1012</u> ☒ Stability Checks		Flow Column #: _____ ☐ 5L/min – 17mm ☐ 15L/min – 53mm ☐ 30L/min – 103mm ☐ R-Value: _____ ☐ Post Adjustment Verification (L/S) Flow Column #: _____ 32 mm _____ (.139-.169) 36 mm _____ (.156-.190) 53 mm _____ (.228-.278) 103 mm _____ (.447-.547)																																						
		<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot#/Exp</th></tr></thead><tbody><tr><td rowspan="2">0.050</td><td rowspan="2">MP6291</td><td>202603A</td></tr><tr><td>3/3/2028</td></tr><tr><td rowspan="2">0.080</td><td rowspan="2">MP6292</td><td>202603B</td></tr><tr><td>3/3/2028</td></tr><tr><td rowspan="2">0.200</td><td rowspan="2">MP6293</td><td>202603D</td></tr><tr><td>3/4/2028</td></tr><tr><td rowspan="2">0.080 DGS</td><td rowspan="2">N/A</td><td>23825080A1</td></tr><tr><td>10/5/2027</td></tr></tbody></table>		Simulator	Serial #	Lot#/Exp	0.050	MP6291	202603A	3/3/2028	0.080	MP6292	202603B	3/3/2028	0.200	MP6293	202603D	3/4/2028	0.080 DGS	N/A	23825080A1	10/5/2027	<table border="1"><thead><tr><th>Maintenance</th><th>By:</th><th>Date:</th></tr></thead><tbody><tr><td>☐ Battery Replacement</td><td></td><td></td></tr><tr><td>☐ Dry Gas Regulator Replacement</td><td></td><td></td></tr><tr><td>☐ Tank Sensor Tare</td><td></td><td></td></tr><tr><td>☐ Breath Tube Replacement</td><td></td><td></td></tr><tr><td>☐ Other:</td><td></td><td></td></tr></tbody></table>		Maintenance	By:	Date:	☐ Battery Replacement			☐ Dry Gas Regulator Replacement			☐ Tank Sensor Tare			☐ Breath Tube Replacement			☐ Other:		
Simulator	Serial #	Lot#/Exp																																								
0.050	MP6291	202603A																																								
		3/3/2028																																								
0.080	MP6292	202603B																																								
		3/3/2028																																								
0.200	MP6293	202603D																																								
		3/4/2028																																								
0.080 DGS	N/A	23825080A1																																								
		10/5/2027																																								
Maintenance	By:	Date:																																								
☐ Battery Replacement																																										
☐ Dry Gas Regulator Replacement																																										
☐ Tank Sensor Tare																																										
☐ Breath Tube Replacement																																										
☐ Other:																																										
Optical Bench Adjustment By: _____		Department Inspection By: <u>KTS</u>																																								
Barometric Pressure Gauge: _____ ID#: _____		Barometric Pressure ID#: <u>28427</u>																																								
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.000</td><td></td><td>N/A</td><td>N/A</td></tr><tr><td>0.040</td><td></td><td></td><td></td></tr><tr><td>0.100</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.300</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.000		N/A	N/A	0.040				0.100				0.200				0.300				0.080 DGS	N/A			Gauge: <u>1012</u> Instrument: <u>1012</u> Mouth Alcohol Solution Lot #: <u>2025-D</u> Exp: <u>9/25/27</u> Acetone Stock Solution Lot #: <u>2025-B</u> Exp: <u>9/22/27</u>												
Simulator	Serial #	Lot #	Expiration																																							
0.000		N/A	N/A																																							
0.040																																										
0.100																																										
0.200																																										
0.300																																										
0.080 DGS	N/A																																									
☐ Post Optical Bench Adjustment Stability Checks		<table border="1"><thead><tr><th>Simulator</th><th>Serial Number</th></tr></thead><tbody><tr><td>0.000</td><td>MP5086</td></tr><tr><td>Interferent</td><td>MP5087</td></tr><tr><td>0.050</td><td>MP5088</td></tr><tr><td>0.080</td><td>MP5089</td></tr><tr><td>0.200</td><td>MP5090</td></tr></tbody></table>		Simulator	Serial Number	0.000	MP5086	Interferent	MP5087	0.050	MP5088	0.080	MP5089	0.200	MP5090																											
Simulator	Serial Number																																									
0.000	MP5086																																									
Interferent	MP5087																																									
0.050	MP5088																																									
0.080	MP5089																																									
0.200	MP5090																																									
<table border="1"><thead><tr><th>Simulator</th><th>Serial #</th><th>Lot #</th><th>Expiration</th></tr></thead><tbody><tr><td>0.050</td><td></td><td></td><td></td></tr><tr><td>0.080</td><td></td><td></td><td></td></tr><tr><td>0.200</td><td></td><td></td><td></td></tr><tr><td>0.080 DGS</td><td>N/A</td><td></td><td></td></tr></tbody></table>		Simulator	Serial #	Lot #	Expiration	0.050				0.080				0.200				0.080 DGS	N/A			Attachments ☒ Form 41 ☒ Stability Checks ☒ Calibration Certificate ☐ Optical Bench Adjustment																				
Simulator	Serial #	Lot #	Expiration																																							
0.050																																										
0.080																																										
0.200																																										
0.080 DGS	N/A																																									
Gauge ID #: _____ Gauge: _____ Instrument: _____		☐ Post-Stability Checks ☐ Flow Adjustment ☐ Form 40 ☒ Other: Calibration Certificate Calculations																																								
Notes/Suggested Service:		☒ Instrument Complies with Chapter 11D-8, FAC ☐ Instrument Does Not Comply with Chapter 11D-8, FAC ☒ Return to/Place into Evidentiary Use ☐ Remain Out of Evidentiary Use ☒ Conduct an Agency Inspection Before Evidentiary Use																																								
		Wen-Chi Pierson Digitally signed by Wen-Chi Pierson Date: 2026.08.31 09:16:21 -04'00' Taylor Gutschow Digitally signed by Taylor Gutschow Date: 2026.09.10 15:52:33 -04'00'																																								
		Tech Review Admin Review																																								

Stability Checks

0.050 g/210L	0.080 g/210L	0.200 g/210L	DGS 0.080 g/210L
0.047 to 0.053 g/210L	0.077 to 0.083 g/210L	0.194 to 0.206 g/210L	0.077 to 0.083 g/210L ≤0.003 g/210L of Wet
SUMTER COUNTY SO Intoxilyzer - Alcotest Analyzer Model 8000 SN 80-000815 08/19/2026 Software: 8100.27 Test g/210L Time Air Blank 0.000 10:24 Control Test 0.050 10:25 Air Blank 0.000 10:25 Control Test 0.049 10:26 Air Blank 0.000 10:26 Control Test 0.049 10:27 Air Blank 0.000 10:27 Control Test Stats Average 0.0493 Std Dev 0.0006 Rel Std Dev(%) 1.1703	SUMTER COUNTY SO Intoxilyzer - Alcotest Analyzer Model 8000 SN 80-000815 08/19/2026 Software: 8100.27 Test g/210L Time Air Blank 0.000 10:39 Control Test 0.078 10:40 Air Blank 0.000 10:40 Control Test 0.078 10:41 Air Blank 0.000 10:42 Control Test 0.077 10:42 Air Blank 0.000 10:43 Control Test Stats Average 0.0777 Std Dev 0.0016 Rel Std Dev(%) 0.7434	SUMTER COUNTY SO Intoxilyzer - Alcotest Analyzer Model 8000 SN 80-000815 08/19/2026 Software: 8100.27 Test g/210L Time Air Blank 0.000 10:05 Control Test 0.197 10:05 Air Blank 0.000 10:06 Control Test 0.195 10:06 Air Blank 0.000 10:07 Control Test 0.197 10:08 Air Blank 0.000 10:08 Control Test Stats Average 0.1963 Std Dev 0.0012 Rel Std Dev(%) 0.5891	SUMTER COUNTY SO Intoxilyzer - Alcotest Analyzer Model 8000 SN 80-000815 08/19/2026 Software: 8100.27 Test g/210L Time Air Blank 0.000 10:47 Control Test 0.081 10:48 Air Blank 0.000 10:48 Control Test 0.080 10:48 Air Blank 0.000 10:49 Control Test 0.079 10:49 Air Blank 0.000 10:50 Control Test Stats Average 0.0800 Std Dev 0.0010 Rel Std Dev(%) 1.2500
<i>Isabella Szygowska</i> Operator's Signature	<i>Isabella Szygowska</i> Operator's Signature	<i>Isabella Szygowska</i> Operator's Signature	<i>Isabella Szygowska</i> Operator's Signature

Florida Department of Law Enforcement Alcohol Testing Program

DEPARTMENT INSPECTION REPORT - INTOXILYZER 8000

Agency: SUMTER COUNTY SO
Time of Inspection: 11:30

Date of Inspection: 08/21/2026

Serial Number: 80-000815
Software: 8100.27

Check or Test	YES	NO	Check or Test	YES	NO
Diagnostic Check (Pre-Inspection): OK	Yes		Date and/or Time Adjusted		No
Minimum Sample Volume Check: OK	Yes		Barometric Pressure Sensor Check: OK	Yes	
Alcohol Free Subject Test: 0.000	Yes		Mouth Alcohol Test: Slope Not Met	Yes	
Interferent Detect Test: Interferent Detect	Yes		Diagnostic Check (Post-Inspection): OK	Yes	

Alcohol Free Test (g/210L)	0.05g/210L Test (g/210L) Lot#:202603A Exp: 03/03/2028	0.08g/210L Test (g/210L) Lot#:202603B Exp: 03/03/2028	0.20g/210L Test (g/210L) Lot#:202603D Exp: 03/04/2028	0.08 g/210L Dry Gas Std Test (g/210L) Lot#:23825080A1 Exp: 10/05/2027
0.000 / 0.000	0.050	0.078	0.198	0.080
0.000 / 0.000	0.050	0.078	0.198	0.080
0.000 / 0.000	0.051	0.079	0.198	0.080
INT / 0.000	0.050	0.078	0.198	0.080
/ 0.000	0.050	0.079	0.198	0.079
/ 0.000	0.051	0.079	0.198	0.080
/ 0.000	0.050	0.079	0.199	0.080
/ 0.000	0.050	0.079	0.198	0.080
/ 0.000	0.051	0.079	0.198	0.080
/ 0.000	0.050	0.079	0.199	0.080

Standard Deviations	0.0004	0.0004	0.0004	0.0003
---------------------	--------	--------	--------	--------

Average Standard Deviation of 0.05, 0.08 and 0.20 g/210L Tests: 0.0003 Number of Simulators Used: 5

Remarks:

00: Interferent Detect.

Re-ran Alcohol free test after turning on fans &
clearing area near instrument. *Katie* 8/21/26

The above instrument complies (☒) does not comply (☐) with Chapter 11D-8, FAC.

I certify that I performed this inspection in accordance with the provisions of Chapter 11D-8, FAC.

Katie

KATIE T SPEARIN

Signature and Printed Name

08/21/2026
Date



Calibration Certificate

Florida Department of Law Enforcement
Alcohol Testing Program
2331 Phillips Road
Tallahassee, FL 32308

This is to certify the calibration of Intoxilyzer 8000 serial number 80-000815, manufactured by CMI, Inc. was calibrated in accordance with FDLE/ATP Form 36 - Department Inspection Procedures - Intoxilyzer 8000.

Serial Number:	<u>80-000815</u>	Average of y ± Uncertainty*
Owning Agency:	<u>SUMTER COUNTY SO</u>	0.050 g/ 210 L 0.004
Calibration Date:	<u>08/21/2026</u>	0.079 g/ 210 L 0.004
Calibration Time:	<u>11:30</u>	0.198 g/ 210 L 0.008
		0.080 g/ 210 L Dry Gas Control 0.005

All results are reported in g/210 L.

Bias is limited by calibration acceptance criteria. All calibration results must be within ± 0.005 or 5%, whichever is greater, of the target alcohol concentration.

*Uncertainty is based on fleet-wide data and is expressed to a 99.73% level of confidence (k=3).

The instrument results before and after any adjustment are found in the associated pre and post stability checks.

TRACEABILITY INFORMATION

This instrument was calibrated using solutions prepared by Alcohol Countermeasure Systems, Inc. (ACS). ACS prepared and certified these CRMs in accordance with ISO 17034 and ISO/ IEC 17025 Standards.

Simulator temperatures are traceable to NIST. Simulator temperatures are checked with NIST traceable digital thermometers calibrated by Precision Metrology in accordance with ISO/ IEC 17025 standards.

Dry gas control measurements are traceable to NIST through the use of CRMs supplied by an accredited CRM supplier. The supplier of dry gas standard controls prepared and certified the CRMs in accordance with ISO Guide 34 and ISO/ IEC 17025 standards.

This document shall not be reproduced except in full, without written approval of the Florida Department of Law Enforcement Alcohol Testing Program.

Kaitlyn Spearin
Spearin
Digitally signed by Kaitlyn
Date: 2026.08.28 08:03:27 -04'00'

08/21/2026

Date

KATIE T SPEARIN,

Department Inspector

FDLE/ATP Form 69 August 2026
Issuing Authority: Alcohol Testing Program

Service • Integrity • Respect • Quality